



**PULSE PROGRAM
CATALOGUE:
I. 1D & 2D NMR
EXPERIMENTS**

Teodor Parella

Servei RMN, Universitat Autònoma de Barcelona
E-mail: teodor.parella@uab.cat



TOPSPIN v3.0
NMRGuide



UAB

Copyright © 2010 by Bruker BioSpin GmbH
All rights reserved. No part of this publication may be reproduced, stored in a retrieval system,
or transmitted, in any form, or by any means without the prior consent of the publisher.

15-04-2010
Part Number -----

Product names used are trademarks or registered trademarks of their respective holders.

VOLUME I: 1D & 2D NMR EXPERIMENTS

• Table of Contents	1
• List of NMR Building Blocks	5
• Introduction	7
• Basic 1D Experiments	
1. Pulse Calibration and Tests.....	25
2. Basic 1D pulse sequences.....	63
Single Pulse-Acquisition	
Double-Resonance Experiments	
Homodecoupling	
Heterodecoupling	
Aring, Udeflt	
3. T1 & T2 Relaxation.....	75
Inversion-Recovery (T1)	
CPMG (T2)	
4. Selective 1D Experiments.....	81
Selective Excitation Methods: SPFGGE and CSSF	
Selective ¹ H: Selective COSY, Selective RELAY, Selective TOCSY	
Selective NOESY, Selective ROESY	
Selective ¹³ C: Selective HSQC, Selective Long-range HSQC,	
Selective INADEQUATE, Selective INEPT (INAPT)	
5. Solvent suppression Methods	97
Presaturation	
Jump and Return	
Watergate	
Excitation Sculpting	
WET	
6. Single/Multiple Presaturation. LC-NMR Applications.....	107
7. ¹⁹ F-specific NMR experiments.....	113
8. ² H-specific NMR experiments.....	117
9. Basic 1D Gradients.....	119
10. ERETIC: Quantitative NMR.....	123
• Homonuclear 2D Experiments	
11. Introduction.....	125
12. 2D COSY.....	133
13. 2D COSY-DQF.....	141
14. 2D SECSY.....	147
15. 2D RELAY.....	149
16. 2D TOCSY.....	153
17. 2D NOESY.....	167
EXSY	
18. 2D ROESY.....	179
T-ROESY	
19. 1D & 2D Multiple-Quantum.....	187
Double-Quantum (DQ)	
Triple-Quantum (TQ)	
20. 2D J-Resolved	193
• Heteronuclear X-detected 2D experiments	
21. 1D INEPT.....	197

22. 1D DEPT	203
DEPT45, DEPT90, DEPT135	
DEPTQ	
23. Other ¹³ C editing methods.....	211
SEFT	
APT	
QUAT	
24. 2D Heteronuclear Correlations.....	213
HETCOR	
COLOC	
25. 2D Heteronuclear J-resolved.....	223
26. 2D HOESY.....	225
27. 1D & 2D INADEQUATE.....	229
• Inverse 2D Experiments	
28. Decoupler Pulse Calibration.....	235
Direct 2D Correlations	
29. HMQC.....	239
30. HSQC.....	253
31. Constant-time correlations.....	277
CT-HSQC	
CT-HMQC	
32. Inverse-HETCOR	285
Inverse INEPT	
Inverse DEPT	
Multiplicity-edited 2D Correlations	
33. Multiplicity-edited HSQC.....	289
34. DEPT-HMQC.....	299
Spin-State Edited 2D Correlations	
35. Spin-state edited HSQC.....	303
2D HSQC- α,β	
2D IPAP-HSQC	
36. TROSY	309
37. CRINEPT.....	319
38. IDIS-HSQC.....	321
2D HMQC hybrids	
39. HMQC-COSY.....	325
H2BC	
HAT-HMBC	
40. HMQC-TOCSY.....	331
41. HMQC-ROESY.....	337
42. HMQC-NOESY.....	341
2D HSQC hybrids	
43. HSQC-TOCSY	345
44. HSQC-ROESY	355
45. HSQC-NOESY	359
2D Long-Range Correlations	
46. HMBC	363
Multiplicity-Edited HMBC	
Constant-Time HMBC	
CIGAR-HMBC	
2J,3J-HMBC	
Simultaneous-CN-HMBC	
47. Measurement of long-range proton-carbon coupling constants	377
J-HMBC	

Long-range HSQC (HSQMBC)	
EXSIDE	
HETLOC	
HSQC-HECADE	
HSQC-TOCSY-IPAP	
Selective J-Resolved	
48. ADEQUATE	387
1,1-ADEQUATE	
1,n-ADEQUATE	
n,1-ADEQUATE	
n,n-ADEQUATE	
• Miscellaneous Experiments	
49. 1D, 2D & 3D Diffusion/DOSY.....	395
STE	
STEBP	
DSTE	
DSTEBP	
LED	
LEDBP	
DOSY-COSY	
DOSY-TOCSY	
DOSY-NOESY	
50. 1D & 2D Saturation Transfer Difference (STD).....	401
STD-TOCSY	
STD-NOESY	
STD-HSQC	
51. 1D & 2D Experiments using CLEANEX.....	413
CLEANEX-HSQC	
CLEANEX-TROSY	
52. Solid-State NMR Experiments.....	417
Basic 1D sequences	
Heteronuclear Cross-polarization (CP)	
Sidebands suppression: TOSS, SELTICS	
Multiplicity-editing: NQS	
Relaxation	
Double-CP	
Homonuclear/Exchange correlation	
CP-NOESY, CP-INADEQUATE, RFDR	
Heteronuclear correlation	
HETCOR, MAS-J-HMQC	
CS-CSA correlation	
PASS	
Multiple-quantum spectroscopy	
BABA, POST-C7, R14_2^6, SC14, SPC5	
Recoupling (REDOR) experiments	
MQMAS for quadrupolar nuclei	
Heteronuclear experiments for quadrupolar nuclei	
PISEMA	
Satellite-transition MAS (STMAS)	
CRAMPS	
BR24, MREV8	
• Appendix 1. Pulse Program Info.....	441
• Appendix 2. Pulse Program Parameters	446

-
- **Appendix 3. Relations with edprosol/getprosol** 454
 - **Appendix 4. Full List of Pulse Programs** 459

VOLUME II: BIOMOLECULAR NMR EXPERIMENTS.....467

NMR Building Blocks

Basic Elements

1. Pre-Scan Delay, the Read Pulse and the Acquisition Period.....	66
2. Broadband vs Selective Decoupling.....	69
3. Purge Element before d1 period.....	74

Filters

4. T2 vs T1rho filters.....	80
5. A Double-Quantum Filter	
6. A z-filter	
7. A Gradient-Based MQF	190
8. The BIRD(x) Operator.....	218
9. BIRD(-x)-Recovery Delay element.....	243
10. G-BIRD Element	
11. Isotope X-Filter using HMQC	
12. Low-pass Filters.....	367
13. Diffusion Filter.....	400
14. Saturation Loop for STD experiments.....	403

Selective Excitation Elements

15. Rectangular (Hard) vs Shaped (Selective) Pulses.....	85
16. Selective Pulsed-Field Gradient Echo (SPFGE)	85

Solvent Suppression Elements

17. Solvent Presaturation	99
18. WATERGATE 3-9-19.....	100
19. WATERGATE (90°selective)	101
20. Excitation Sculpting W5.....	102
21. Excitation Sculpting (180°selective)	103
22. WET.....	104

Evolution Periods

23. A Variable Evolution Period.....	136
24. 1H evolution and X-decoupling.....	216
25. Indirect Chemical Shift Evolution	
26. Constant-Time periods	280
27. Multiplicity-Editing.....	292

Homonuclear Mixing Blocks

28. COSY vs RELAY mixing times	
29. MLEV-17 as a TOCSY mixing time	158
30. Z-filtered DIPSI-2 as a TOCSY mixing time.....	159
31. NOESY mixing time.....	171
32. ROESY mixing time.....	182
33. T-ROESY mixing time.....	183

1H-to-X Transfer

34. Polarization Transfer via INEPT.....	200
35. The Refocused INEPT	
36. INEPT Train with a trim pulse.....	260
37. Purge Gradient in INEPT for zz-selection.....	260
38. Water Flip-back in INEPT. 90°-Gradient purge element.....	261
39. A CPMG-INEPT Block.....	274

X-to-1H Transfer

40. Reverse-INEPT: Normal vs WATERGATE	
41. Sensitivity-Improved (PEP) in HSQC.....	264
42. Spin-State Selection in Reverse INEPT.....	306
43. IPAP Editing in the indirect dimension.....	307
44. TROSY and Half-TROSY Blocks.....	312
45. TROSY and WATERGATE	
46. Clean-TROSY using Echo-Antiecho Gradient Selection	

Miscellaneous Elements

1. Generation of Homonuclear Multiple-Quantum Coherences.....	189
2. Generation of Heteronuclear Multiple-Quantum Coherences.....	242
3. Bipolar Gradients	
4. Bloch-Sieger Phase Shifts	

Triple-Resonance Elements

5. The Double ¹³ C, ¹⁵ N-Filter for Isotope-Edited Experiments.....	486
6. Concatenated Building Blocks in Triple-Resonance Experiments.....	542
7. Semi-Constant Time Periods in Triple-Resonance Experiments.....	543
8. Chemical Shift vs J Evolution	
9. The MUSIC Element	
10. T1 Filter in Heteronuclear Relaxation Experiments	
11. T2 Filter in Heteronuclear Relaxation Experiments	
12. T1rho Filter in Heteronuclear Relaxation Experiments	
13. 2H-decoupling in Triple-Resonance Experiments	
14. Band-Selective ¹³ C Pulses in Triple-Resonance Experiments	
15. Isotope-Filtered vs Isotope-Edited	

BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

INTRODUCTION

BRUKER Pulse Program Catalogue

written by Teodor Parella

This catalogue presents the pulse sequence diagram for all standard pulse programs included by default in TOPSPIN v3.0. This information is part of NMRGuide, which is also available for BRUKER AVANCE spectrometers.


These pulse programs are located in the

`/TOPSPINHOME/exp/stan/nmr/lists/pp/`

directory after conventional installation using **expinstall** and they can be visualized directly into the TOPSPIN program from the PulseProg section. Otherwise, alternative pulse program sequence graphical display is also available using the **showpp** program.

For more details on pulse programs, parameter sets, tutorials, experiment descriptions, bibliographic references and other related information, please refer to the electronic version of NMRGuide.


Starting NMRGuide ...




More Info using NMRGuide: The ghelp search tool

The **ghelp** program allows the direct user interaction with NMRGuide from the TOPSPIN command line. Otherwise, select NMRGuide from the Help menu to open the Explorer browser. The ghelp options are:

- tut**: Search into AVANCE Tutorials
- pp**: Search into pulse programs
- au**: Search into AU programs
- par**: Search into parameter sets
- exp**: Search into Experiment Wizards
- spec**: Search into Spectra Database
- doc**: Search into Bruker Documentation
- bib**: Search into Bibliography
- cmd**: Search into TOPSPIN commands



For instance, if you want information about pulse programs related to "HNCO" you can type "ghelp -pp hnco" into the command line as shown:



The output looks like this:

The screenshot shows the NMRGuide search tool interface. At the top, there is a search bar with 'hncq' entered and a 'GO' button. The page title is 'Search Tool'. Below the search bar, the results for '3D HNCQ experiment' are displayed. It lists available files, sensitivity improvement using PEP, and experiments using TROSY and WATERGATE. A chemical structure diagram is shown below the text, illustrating the HNCQ experiment setup with a backbone atom (N-C-H) and a side chain (C-N-H) connected by a vertical line. The diagram is labeled with $^1J_{NH} + ^1J_{NCO}$ and HNCQ. At the bottom of the interface, there are navigation buttons: a zoom icon, a copy icon, a refresh icon, and left/right arrow icons.

and a browser window will be automatically opened with the keywords found in the NMRGuide database. Then, further specific search can be performed from the buttons:

The screenshot shows the NMRGuide 4.0 search tool interface. At the top, there is a search bar with 'NMRGuide' entered and a 'GO' button. The page title is 'NMRGuide 4.0 - Search Tool'. Below the search bar, there is a table of navigation buttons. A red circle highlights the buttons: 'NMRGuide', 'AVANCE Experiment Tutorials', 'Wizards', 'TOPSPIN Commands', and 'NMR Spectra'. Below the table, there is a section titled 'NMR Experiments' with a list of links: '2Q-HMBC', 'ACCORD-HMBC', 'CIGAR-HMBC', 'HAT-HMBC', 'HMBC: Heteronuclear Multiple-Bond Connectivities', 'J-HMBC: J-resolved HMBC', 'CT-HMBC: Constant-Time HMBC', 'Selective HMBC', and 'Semi-selective HMBC'.

AUTOMATION INTO TOPSPIN




Several Options for automated Acquisition/Processing into TOPSPIN:

1. Using Standard Parameter Sets (rpar and getprosol)
2. Using macros (edmac)
3. Using Guides
4. Using ButtonNMR (buttonnmr & butselnmr)
5. Using ICONNMR (iconnmr)

Other Automation acquisition and processing protocols using AU Programs. (edau & xau)




Multiple sample/Multiple Experiment definition



Real-time sample status also available via www

Automation using NMR BIOTOOL: Automated Biomolecular NMR using a guided, full-automated process



Automation using AU Programs:

Au programs is the best option to perform full automated processes.
See NMRGuide and "AU Reference manual" for more information.

Some Features:

1. More than 400 programs by default
2. Written in C and python.
3. TOPSPIN commands are fully compatible
4. Edited with **edau** and executed automatically with the corresponding name.
6. Included into automated Acquisition and Processing modules under AUNM, and AUNMP.

Examples:

topshim: Automatic Shimming process
pulsecal: Performs automatic 90 proton pulse calibration
multizg: Performs multiple acquisition
multi_zgvd: Performs multiple acquisition at predefined times
multi_zgvt: Performs multiple acquisition at predefined temperatures
calcphomo: Performs automatic phase calibration in 2D homonuclear experiments


Automation using MACROS:

To avoid much and repetitive typing, the use of macros can be highly useful. Macros are invoked with the **edmac** command and are very simple to create. As an example, we create a macro called **proton** to acquire and process automatically a conventional 1D proton spectrum:


1. Type **edmac proton**
2. Write the following sequential commands
3. Execute the macro simply typing the name of the macro.

In this case, type **proton** and the macro will ask the necessary inputs and then will execute acquisition and processing commands.

```
edc
rpar teo_proto all
ii
ns
rga
zg
ft
apk
abs
```


aqguide




prguide




managuide



solaguide




t1guide



daisyguide



buttonnmr



butselnmr

1. Sample Preparation



- Filter the sample.
- Use clean NMR tubs of optimum quality.
- Use the required volume (0.6 ml for 5mm NMR tubes).
- Sample concentration as a function of NMR experiment sensitivity.




2. NMR Preparation

- 2.1 Put the sample into the magnet.
- 2.2 Select solvent.
- 2.3 Spin the sample.
- 2.3 Shimming.
- 2.4 Tuning and Matching .





3. Data Acquisition.

edc	: Define the file
rpar "filename" all	: Read starting parameter set
getprosol	: Get pulses and power levels
ased (or eda)	: Check original parameters and change if required
rga	: Adjust detector
zg	: Start Acquisition

4. Data processing

	Subroutine	1D	2D
Fourier Transformation		ft	xfb (xf2+xf1)
Phase Correction		apk	xfbp (xfb+pk)
Baseline Correction		abs	abs2 + abs1
Automatic Referencing		sref	

5. Data Analysis

-  Enter peak picking mode
-  Enter integration mode
-  Enter multiple display mode
-  Enter distance measurement mode

1. NMR Strategy

Molecular Size
Natural Abundance vs Isotopic Labeling
Structural vs Dynamic studies

2. Which NMR Experiment?

Sample Concentration
Available Probehead
Magnetic Field
NMR Parameters →
Which Nucleus?

Chemical Shift
Coupling Constants: J/D
NOE/ROE
Diffusion
Relaxation
Chemical Exchange

3. Which Version?

1D/2D/3D ???
Pulsed Field-gradients or phase-cycled
Qualitative or quantitative analysis?
Magnitude-mode or phase-sensitive
Sensitivity vs Resolution
Direct detection or Inverse Spectroscopy
Solvent Suppression?

4. Others

Spectrometer Set-up
Experiment Optimization: Pulse Calibration
Variable Temperature
Automation

Looking for a pulse program or NMR experiment using the **NMR Experiment Selector** tool

Which NMR experiment/version to use??:
Select the NMR Experiment

Molecule under study:

- Small-to-medium sized natural-abundance
- Unlabeled Protein | Unlabeled Nucleic Acid
- ¹⁵N-Labeled Protein
- ¹³C, ¹⁵N-Labeled Protein
- ²H, ¹³C, ¹⁵N-Labeled Protein
- Labeled Protein/Unlabeled Ligand Complex
- Labeled Nucleic Acid

Dimensionality: 1D 2D 3D

Use of Gradients?: No Yes

Data Representation: Magnitude-mode Phase-sensitive

Solvent Supression?:

- No
- Yes : Presaturation WATERGATE Excitation Sculpting WET

Inverse Spectroscopy Options:

- Using f2 channel Using f3 channel
- PEP TROSY ²H-decoupling
- Adiabatic Pulses

NOESY Options:

Isotope Editing/Filtering?

- No
- Yes: F1 F2 F1/F2

With simultaneous X,Y evolution


AVANCE Tutorials

NMR Assistance

NMR Experiment	Parameter Set	Pulse Program	Modify ...
Conventional 1D Experiments			
1D ¹H	PROTON	zg30	d1, ns, o1p, sw
1D ¹³C (¹H)	C13CPD	zpgg30	d1, ns, o1p, sw
1D DEPT-135	C13DEPT135	dept135	d1, ns, o1p, sw
1D DEPT-90	C13DEPT90	dept90	d1, ns, o1p, sw
Conventional 2D 1H-1H Experiments			
ge-2D ¹H-¹H COSY (magnitude)	COSYGPSW	cosygpqf	d1, ns, o1p, 1 sw, 2 sw
2D ¹H-¹H NOESY (phase sensitive)	NOESYPHSW	noesygpqh	ns, o1p, 1 sw, 2 sw, d8
2D ¹H-¹H ROESY (phase sensitive)	ROESYPHSW	roesyph	ns, o1p, 1 sw, 2 sw, p15
2D ¹H-¹H TOCSY (phase sensitive)	MLEVPHSW	mlevph	ns, o1p, 1 sw, 2 sw, d9
2D Inverse 1H-¹³C Experiments			
ge-2D ¹H-¹³C HMQC (magnitude)	HMQCGP	hmqcgp	ns, o1p, o2p, 1 sw, 2 sw
ge-2D ¹H-¹³C HSQC (phase sensitive)	HSQCGP	hsqcetgpsi2	ns, o1p, o2p, 1 sw, 2 sw
ge-2D ¹H-¹³C HSQC edited (phase sensitive)	HSQCEDETGP	hsqcedetgp	ns, o1p, o2p, 1 sw, 2 sw
ge-2D ¹H-¹³C HMBC (magnitude)	HMBCGP	hmbcgpplndqf	d6, ns, o1p, o2p, 1 sw, 2 sw
Selective 1D Experiments			
Selective 1D NOESY	SELNOGP	selnogp	ns, spoffs2, p12, sp2, d8
Selective 1D ROESY	SELROGP	selrogp	ns, spoffs2, p12, sp2, p15
Selective 1D TOCSY	SELMLGP	selmlgp	ns, spoffs2, p12, sp2, d9
31P Experiments			
1H-decoupled 31P	P31CPD	zpgg30	ns, sw, o1p
1H-coupled 31P	P31	zg30	ns, sw, o1p
31P-decoupled 1H spectrum	PROP31DEC	zggd30	ns, o2p
2D 1H-31P HMBC	HMBCGP	hmbcgpndqf	ns, o2p, 1 sw, 1 td, d6

Pulse Sequence Diagram: A brief Outlook.

This is a classical NMR pulse sequence representation:



We can observe the following Pulse Sequence Elements:

1. Different Rf Channels: ^1H , ^{13}C , ^{15}N , G_z ...
2. Radiofrequency (Rf) Pulses (90° , 180° ...) applied at specific power levels
3. Inter-pulse Fixed Delays (d_1 , d_4 , d_{24} ...) conveniently optimized to J , NOEs, relax ...
4. Variable delays that define dimensionality: d_0 and d_{10}
5. Multiple-pulse sequences: trim pulses, GARP decoupling, TOCSY transfers ...
6. Phase Cycling: Only the most important are usually described (see pulse program for a detailed description)
7. Gradients: Coherence Selection & Purge Elements
8. Starting Point: Important for sensitivity and repetition rates aspects
9. Final Point (see the FID): Important for sensitivity

Can I run a given NMR experiment on my spectrometer? This is a typical question for non-experienced users. It is needed to consider three independent aspects: Software, Hardware and Sample Requirements:


In principle, there is no limitations due to SOFTWARE requirements using TOPSPIN.

The main limitations came from HARDWARE requirements:

1. The number of available Rf channels
2. Inverse Spectroscopy capabilities
2. Availability of Pulsed-Field Gradients
3. Check for the available Probeheads!!!!

Another important question: SAMPLE requirements???



Many experiments are usually driven in natural abundance samples but in biomolecules are necessary some isotopic labeling strategies: 15N-labeled, doubly-15N,13C-labeled, partial or fully 2H combined with 15N,13C-labeled, or other selective labeling approaches




The NMR Pulse Sequence: Diagram vs Microprogram

Any NMR Experiment is closely related to a pulse sequence that can be analyzed in detail using TOPSPIN.

The pulse sequence diagram



The pulse program



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

STANDARD TESTS & PULSE CALIBRATIONS

BASIC TEST ON AVANCE SPECTROMETERS:

- Lineshape test for ^1H without rotation (lsnh)
- Lineshape test for ^1H with rotation (lsrh)
- Resolution test for ^1H (rsh)
- Determination of 90 degree ^1H high power transmitter (tph)
- Sensitivity test for ^1H (snh)
- Water suppression test (ws)
- Lineshape test for ^{13}C (lsc)
- Resolution test for ^{13}C (rsc)
- Determination of 90 degree ^{13}C high power transmitter pulse with popt (tpc)
- Determination 90 degree 1H pulse for high-power decoupling (CPD pulses) (dph)
- Determination 90 degree 1H pulse for low-power decoupling (CPD pulses) (cph)
- Sensitivity test for 13C ASTM (without 1H decoupling) (sna)
- Sensitivity test for 13C EB (with 1H CPD decoupling) (snc)
- Determination 90 degree 13C high power decoupling (puc)
- Determination 90 degree 13C low power decoupling (for GARP) (gac)
- Inverse spin-echo difference test (inv)
- Determination 90 degree 31P high power
- Sensitivity test for 31P (snp)
- Sensitivity test for 15N (snn)
- Determination 90 degree 15N high power decoupling (pun)
- Determination 90 degree 15N low power decoupling (GARP) (gan)
- Gradient Recovery Test
- Determination of selective shaped 90 degree 1H pulse
- Selective Excitation test
- Selective Experiment test
- Homodecoupling test (hde)
- Sensitivity test for 19F (snf)
- GARP decoupling 13C (garc)
- Decoupling test 19F (with 1H CPD decoupling) (decf)
- Triple Resonance experiment (tri)
- Dataset for quad image adjustment (qad)
- Decoupling test 19F (with 19F GARP decoupling) (garf)
- B1 homogeneity test 1H from HWT (b1h)

NMR Experiment: Lineshape test for ^1H without rotation (lsnh)
Basic Parameter Set: PROHUMP
Pulse Program: zg30

Sample: 1% CHCl_3 in Acetone- d_6

Spin: off


Basic acquisition parameters:

d1 2s


sw=1000Hz

o1p on the CHCl_3 resonance

ns=1 ds=0



Analysis: The linewidth is automatically determined at 0.55% of the chloroform line and at one fifth there of (0.11%) by the **humpcal** program



NMR Experiment: Lineshape test for ^1H with rotation (lsrh)
Basic Parameter Set: PROHUMP
Pulse Program: zg30

Sample: 1% CHCl_3 in Acetone- d_6

Spin: on


Basic acquisition parameters:

d1 2s


sw=1000Hz

o1p on the CHCl_3 resonance

Ns=4 ds=0



Analysis: The linewidth is automatically determined at 0.55% of the chloroform line and at one fifth there of (0.11%) by the **humpcal** program



NMR Experiment: Resolution test for ^1H (rsh)
Basic Parameter Set: PRORESOL
Pulse Program: zg30

Sample: 1% CHCl_3 in Acetone- d_6

Spin: on


Basic acquisition parameters:

d1 1s


sw=1000Hz

o1p on the CHCl_3 resonance

ns=1 ds=0



Analysis: The resolution is measured at half height of the chloroform line and calculated with **humpcal**



NMR Experiment: Determination of 90 degree ^1H high power transmitter pulse with popt (tph)
Basic Parameter Set: PROTON
Pulse Program: zg

Sample: 0,1% Ethylbenzene in CDCl_3

Spin: on

Basic acquisition parameters:


d1 120s

sw 3000Hz

o1p on the methylene resonance


ns=1 ds=0

Process with ef (lb=1)



Analysis: The pulse length strongly depends on the transmitter power and probe design. Remember to tune and match before you start calibration with **popt**. The first maximum is the 90° , the first null point is 180° and the second 360°

PULSE CALIBRATION



NMR Experiment: Sensitivity test for ^1H (snh)
Basic Parameter Set: PROSENS
Pulse Program: zg

Sample: 0,1% Ethylbenzene in CDCl_3

Spin: on

Basic acquisition parameters:


d1 60s

sw 3000Hz

ns=1 ds=0

Set p1 to the calibrated 90 degree value at p1 power level

Process with ef (lb=1)




Analysis: the signal-to-noise ratio is determined on the intensity of the quartet multiplet by **sino** and strongly depends on good resolution and lineshape.

SINO = 650.9

noise from 5.99
to 3.99 ppm

SINO = 776.2:1

noise from 4.47
to 4.07



NMR Experiment: Water suppression test (ws)
Basic Parameter Set: WATERSUP
Pulse Program: zgpr

Sample: 2mM sucrose with 0.5mM DSS, 2mM NaN₃ in 10% D₂O and 90% H₂O

Spin: off

Basic acquisition parameters:

d1 5s

set aq=1s and sw=12ppm (td is automatically calculated)


Set exactly o1p on the water resonance (about 4.7ppm)

by minimizing the FID in gs mode


ns=8 ds=4

Set pl9 between 55-65dB (50Hz rf field)

Process without window function



Analysis: The linewidth is measured at 50% height of the DSS signal and a 10% thereof. The S/N ratio is measured on the doublet left to the water signal using sino. On this doublet the splitting (resolution) is also measured.



NMR Experiment: Lineshape test for ^{13}C (lsc)
Parameter Set: C13HUMP
Pulse Program: zgcw30

Sample: ASTM (60% C6D6 / 40% p-Dioxane)

Spin: on

Basic acquisition parameters:

d1 10s


ns 4 and ds 0

o1p 66.5ppm

Optimize o2 (start with 3.39ppm) in gs mode by maximizing the FID


Process without window function

Analysis: The linewidth is measured automatically with humpcal at 0.55% and 0.11% levels.



```
- FINAL TEST - System: AV300 OrderNo.: JH030202 Customer: Engineer: GRR / SKD  
P/N Console: H203128/1203 Shim system: BOSSI original dataset: 8352_01551sc 2 1  
Probe: 5 mm QNP 1H/13C/31P/19F Z-GRD Z8352/0155 Sample depth:20 Gas: air  
Lineshape test for 13C with rotation; Sample: ASTM in C6D6 (P/N: Z10163)
```

```
Hump= 1.4/3.1 Hz  
P1= 11usec PL1= 0dB PL26= 20dB
```



NMR Experiment: Resolution test for ^{13}C (rsc)
Parameter Set: C13RESOL
Pulse Program: zgcw30

Sample: ASTM (60% C₆D₆ / 40% p-Dioxane)

Spin: on

Basic acquisition parameters:

d1 10s


ns 1 and ds 0

o1p 66.5ppm

Optimize o2 (start with 3.39ppm) in gs mode by maximizing the FID


Process without window function

Analysis: The resolution is measured automatically with humpcal at half height of the signal.



```
- FINAL TEST - System: AV300   OrderNo.: JH030202   Customer:   Engineer: GRR / SKO  
P/N Console: HZ03128/1203   Shim system: BOSSI   original dataset: 8352 0155rsc 2 1  
Probe: 5 mm QNP 1H/13C/31P/19F Z-GRD Z8352/0155   Sample depth:20   Gas: air  
Resolution test for 13C; Sample: ASTM in C6D6 (P/N: Z10163)
```

```
Resolution= 0.18 Hz  
P1= 11usec   PL1= 0dB   PL26= 20dB
```



NMR Experiment: Determination of 90 degree ^{13}C high power transmitter pulse with popt (tpc)

Basic Parameter Set: C13CPD

Pulse Program: zg

Sample: ASTM (60% C₆D₆ /40% p-Dioxane)

Spin: on

Basic acquisition parameters:


d1 80s (very important to use a long delay!!)

sw 15000Hz


o1p on-resonance

ns=1 ds=0

Process with ef (lb=3.5)



Analysis: The pulse length strongly depends on the transmitter power and probe design. Remember to tune and match before you start calibration with **popt**. The first maximum is the 90°, the first null point is 180° and the second 360°



NMR Experiment: Determination 90 degree 1H pulse for high-power decoupling (CPD pulses) (dph)
Basic Parameter Set: C13CPD
Pulse Program: decp90

Sample: ASTM (60% C6D6 / 40% p-Dioxane)

Spin: on

Basic acquisition parameters:


d1 45s

ns 1 and ds 0


o1p 66.5ppm; o2p=3.7ppm

CNST2=142Hz; d203.52ms

Process with ef (lb=3.5Hz)



Analysis: Without decoupling we observe an antiphase triplet, with the two outer lines in phase and the center line in opposite phase. Set pl2 (for instance, to -3dB) and varies p3 using **popt** until the triplet becomes zero.



NMR Experiment: Determination 90 degree 1H pulse for low-power decoupling (CPD pulses) (cph)
Basic Parameter Set: C13CPD
Pulse Program: decp90

Sample: ASTM (60% C6D6 / 40% p-Dioxane)

Spin: on


Basic acquisition parameters:

d1 45s

ns 1 and ds 0


o1p 66.5ppm; o2p=3.7ppm

CNST2=142Hz; d2=3.52ms



Process with ef (lb=3.5Hz)

Analysis: Without decoupling we observe an antiphase triplet, with the two outer lines in phase and the center line in opposite phase. Set p3=80us and start with a low power pl2=30dB. Increase pl2 using **popt** until the triplet becomes zero. These values will be used for CPD 1H decoupling experiments.



NMR Experiment: Sensitivity test for ^{13}C ASTM (without ^1H decoupling) (sna)
Basic Parameter Set: C13SENS
Pulse Program: zg

Sample: ASTM (60% C_6D_6 / 40% p-Dioxane)


Spin: on

Basic acquisition parameters:

d1 300s


p1 must be the calibrated 90° pulse at p11

ns 1 and ds 0



Process with ef (lb=3.5Hz)

Analysis: The S/N is determined on the triplet of the deuterated benzene at 128ppm using **sino**. The splitting of the 1:1:1 triplet should go lower than 9%.



NMR Experiment: Sensitivity test for ^{13}C EB (with ^1H CPD decoupling) (snc)

Basic Parameter Set: C13SENS

Pulse Program: zgdc

Sample: 10% Ethylbenzene in CDCl_3

Spin: on

Basic acquisition parameters:


d1 120s

p1 must be the calibrated 90° pulse at p1


ns 1 and ds 0

o1p=80ppm, o2p=5ppm, sw=200ppm and pl12 for CPD

Process with ef (lb=0.3Hz)



Analysis: The S/N is determined on the highest peak of the aromatic part. It is calculated using *sino* over a range of 40ppm between 30-125ppm.



NMR Experiment: Determination 90 degree ^{13}C high power decoupling (puc)
Basic Parameter Set: HMQC1D
Pulse Program: decp90

Sample: 100mM Urea 15N, 100mM CH_3OH ^{13}C in DMSO-d_6

Spin: off

Basic acquisition parameters:


d1 30s

ns 1 and ds 0


o1p=4ppm o2p=49.5ppm

CNST2=139Hz (d2=3.59ms)

Process with ef (lb=0.3Hz)



Analysis: Without decoupling we observe antiphase signals at 3.2ppm. Increase p3 at power level pl2 until the signals become zero.



NMR Experiment: Determination 90 degree ^{13}C low power decoupling (for GARP) (gac)
Basic Parameter Set: HMQC1D
Pulse Program: decp90

Sample: 100mM Urea 15N, 100mM CH₃OH ^{13}C in DMSO-d₆

Spin: off

Basic acquisition parameters:


d1 2s

ns 1 and ds 0


o1p=4ppm o2p=49.5ppm

CNST2=139Hz (d2=3.59ms)

Process with ef (lb=0.3Hz)



Analysis: Without decoupling we observe antiphase signals at 3.2ppm. Set p3=80us and optimize the power level p12 until the signals become zero.



Sample: 1% or 3% CHCl₃ in acetone-d₆

Spin: off

Basic acquisition parameters:

d1 2s

ns 1 and ds 0


o1p=4ppm sw=8ppm

o2p=77ppm



CNST2=214Hz (d2=2.336ms)

Process with ef (lb=0.3Hz)

Analysis: Without decoupling we observe antiphase satellites. Increase p3 until the ¹³C pulse satellites become zero. Use a suitable power level pl2 for 90 degree high power level pulse F2. Use all necessary filters.



NMR Experiment: Inverse spin-echo difference test (inv)
Basic Parameter Set: HMQC1D
Pulse Program: hmqcndrd1d



NMR Experiment: Determination 90 degree 31P high power
Basic Parameter Set: P31
Pulse Program: zg

Sample: 0.0485M Triphenylphosphate (TPP) in acetone-d6


Spin: on

Basic acquisition parameters:

d1 60s


ns 1 and ds 0

o1p=-16ppm



Process with ef (lb=5Hz)

Analysis: The pulse length strongly depends on the transmitter power and probe design. Remember to tune and match before you start calibration with **popt**. The first maximum is the 90°, the first null point is 180° and the second 360°



NMR Experiment: Sensitivity test for 31P (snp)
Basic Parameter Set: P31
Pulse Program: zg

Sample: 0.0485M Triphenylphosphate (TPP) in acetone-d₆

Spin: on


Basic acquisition parameters:

d1 60s

ns 1 and ds 0


o1p=-16ppm

Set p1 at the 90° degree pulse at power level p1.



Process with ef (lb=5Hz)

Analysis: The 31P sensitivity is determined on the fully coupled resonance line of TPP. S/N is calculated using sino over a range of 5 ppm.



NMR Experiment: Sensitivity test for 15N (snn)
Basic Parameter Set: N15IG
Pulse Program: zgig

Sample: 90% Formamide in DMSO-d6

Spin: on


Basic acquisition parameters:

d1 120s

ns 1 and ds 0

o1p=112.5ppm o2p=7.3ppm sw=20ppm

Set pl12 as determined for CPD 1H decoupling.




Process with ef (lb=0.3Hz)

Analysis: The nitrogen sensitivity is measured on the formamide resonance using **sino** over a range of 2ppm. Before this, you must determine the 90 degree pulse 15N (p1) at power level pl1 using **popt**.

```
- FINAL TEST - System: AV300 OrderNo.: ZH038702 Customer: Engineer: VOM  
P/N Console: H03128/1190 Shim system: BOSSI original dataset: 3934_0196snn 2 1  
Probe: 5 mm EBO BB-1H Z3934/0196 Sample depth:20 Gas:  
Sensitivity test for 15N; Sample: 90% Formamide in DMSO-d6 (P/N: Z10187)
```

```
Sino= 20:1 ( signal= 113 - 111 ppm noise= 116.17 - 114.17 ppm [2 ppm] noise range= 9 ppm )  
Pl= 10usec PL1= -6dB
```



NMR Experiment: Determination 90 degree 15N high power decoupling (pun)

Basic Parameter Set:

Pulse Program: `dec90`

Sample: 100mM urea 15N, 100mM Ch3OH 13C in DMSO-d6

Spin: off

Basic acquisition parameters:


`d1 24s`

`ns 1 and ds 0`


`o1p=4ppm o2p=76ppm sw=8ppm`

`CNST2=88.5Hz (d2=5.649ms)`

Process with `ef` (`lb=0.3Hz`)



Analysis: Without decoupling we observe an antiphase doublet. Increase `p3` (at power level `pl2`) until the doublet becomes zero. For three-channel system use `pulprog=dec90f3`.



NMR Experiment: Determination 90 degree 15N low power decoupling (GARP) (pun)
Basic Parameter Set:
Pulse Program: decp90

Sample: 100mM urea 15N, 100mM Ch3OH 13C in DMSO-d6

Spin: off


Basic acquisition parameters:

d1 2s

ns 1 and ds 0

o1p=4ppm o2p=76ppm sw=8ppm


CNST2=88.5Hz (d2=5.649ms)




Process with ef (lb=0.3Hz)

Analysis: Without decoupling we observe an antiphase doublet. Set p3=180-200us and decrease power level pl2 starting with a pl2=30dB until the doublet becomes zero. For three-channel system use pulprog=decp90f3.


p3 set to 120u





NMR Experiment: Gradient Recovery Test
Basic Parameter Set: PROTON
Pulse Program: preempg2





10 50 100u 200u 300u





**NMR Experiment: Determination of selective shaped
90 degree 1H pulse
Basic Parameter Set: SELZG
Pulse Program: selzg**



NMR Experiment: Selective Excitation test
Basic Parameter Set: SELZG
Pulse Program: selzg



NMR Experiment: Selective Experiment test
Basic Parameter Set: SELCO
Pulse Program: selco



NMR Experiment: Homodecoupling test
Basic Parameter Set: PROHOMODEC
Pulse Program: zgpd

Sample: 0.1% Ethylbenzene in CDCl₃

Spin: off



Basic acquisition parameters:

d1 3

ns 8 and ds 2

Set pl24 for hd decoupling (45-50dB)

Process with ef (lb=1Hz)



NMR Experiment: Sensitivity test for 19F (snf)
Basic Parameter Set: F19
Pulse Program: zg

Sample: 0.05% Trifluorotoluene (TFT) in CDCl3

Spin: on

Basic acquisition parameters:


d1 20s (T1 is about 3.5s)

ns 1 and ds 0

o1p=-63ppm sw=10ppm


Process with ef (lb=2Hz)

Analysis: Determine the 90 degree pulse p1. The fluorine sensitivity is measured on the TFT resonance by using the **sino** program over a range of 1 ppm.



```
- FINAL TEST - System: AV300   OrderNo.: JH030202   Customer:   Engineer: GRR / SKO
P/N Console: HZ03128/1203   Shim system: BOSSI
Probe: 5 mm QNP 1H/13C/31P/19F Z-GRD Z8352/0155   Sample depth:20   Gas: air
Sensitivity test for 19F; Sample: 0,05% TFT in CDCl3 (P/N: Z10234)

Sino= 229:1 ( signal= -61 - -65 ppm   noise= -61.42 - -62.42 ppm [1 ppm]   noise range= 4.5 ppm )
Pl= 5.75usec   PL1= -5dB
```



NMR Experiment: GARP decoupling 13C
Basic Parameter Set:
Pulse Program: zgig

Sample: 100mM Urea 15N, 100mM CH3OH 13C in DMSO-d6

Spin: off

Basic acquisition parameters:

d1 20s

ns 4 and ds 2


o1p=4ppm sw=8ppm AQ=1s

o2p=49.5ppm (13C)

pl12 is the determined power level for garp 13C decoupling


cpdprg2=garp

pcpd2=65...80us



Process with ef (lb=0.3Hz)

Analysis: Use all necessary filters and also an extra 2H-stop filter in the 13C channel at the probehead input. Be aware that a GARP decoupling sequence produces always modulation spikes around large signals on a level of <0.5%.



NMR Experiment: Decoupling test 19F (with 1H CPD decoupling) (decf)
Basic Parameter Set: F19CPD
Pulse Program: zgfhigqn

Sample: 0.05% Trifluorotoluene (TFT) in CDCl₃

Spin: on

Basic acquisition parameters:


d1 1s

ns 1 and ds 0

o1p=-63ppm sw=10ppm

o2p=5ppm


pl12 is the determined power level for cpd decoupling




Process with ft

Analysis: Use all necessary filters

```
- FINAL TEST - System: AV300 OrderNo.: JHD30202 Customer: Engineer: GRR / SKO  
P/N Console: H203128/1203 Shim system: BOSSI  
Probe: 5 mm QNP 1H/13C/31P/19F Z-GPD 28352/0155 Sample depth: 20 Gas: air  
Decoupling experiment 19F observe, 1H CPD-decoupling; Sample: 0,05% TFT in CDCl3 (P/N: Z10234)
```



No decoupling



1H CPD
decoupled

NMR Experiment: Triple Resonance experiment
Basic Parameter Set:
Pulse Program: zgfbig

Sample: 100mM Urea 15N, 100mM CH3OH 13C in DMSO-d6

Spin: on

Basic acquisition parameters:

d1 5s

ns 4 and ds 2

o1p=4ppm sw=10ppm


o2p=49.5ppm (13C)

o3p=76ppm (15N)

cpdprg2=cpdprg3=garp

pcpd2=65....80us


pcpd3=180...200us



Process with ef (lb=0.3Hz)

Analysis: It is important to use enough delay in between two scans to keep the temperature stable during the experiment. Use all necessary filters (13C-pass 2H-stop and 15N-pass 2H stop)

```
- FINAL TEST - System: AV600 OrderNo.: ZH054502 Customer: Beijing Inst. Chem. CN Engineer: VOM
P/N Console: H9118/1342 Shim system: E033 II
Probe: 5 mm TXI 1H-13C/15N XYZ-GHD Z0562/0002 Sample depth: 20 Gas:
Triple resonance experiment 1H, 13C, 15N (additional test): Sample: 0,1M Urea 15N, 0,1M CH3OH 13C in DMSO-d6 (P/N: Z10262)
```



NMR Experiment: Dataset for quad image adjustment
Basic Parameter Set:
Pulse Program: zg30

Sample: 0.1mg GdCl₃/ml + 1% H₂O

Spin: off


Basic acquisition parameters:

d1 1s

ns 1 and ds 0

o1p=7ppm sw=10ppm

AQMOD=QSIM




Process with ef (lb=1Hz)

Analysis: Start the Rx22 service tool from Unitool and adjust first DC offsets before phase and amplitude calibration. The residual quad image must be smaller than 1%.

```
- FINAL TEST - System: M300 OrderNo.: ZH038702 Customer: i Engineer: VOM  
P/N Console: H03128/1190 Shim system: BOSSI original dataset: 3934_0196quad 2 1  
Probe: 5 mm BBO BB-1H Z3934/0196 Sample depth:20 Gas:  
quad image adjustment (optional experiment) (additional test); Sample: 0,1 mg GdCl3/ml D2O 1* H2O (F
```

```
RX22 Gain Adjust Value = -2.52 *Gain  
RX22 Phase Adjust Value = -2.21 degrees  
Current Channel A Pattern: 127 (0x07F)  
Current Channel B Pattern: 127 (0x07F)
```

```
Quadimage RG=512, NS=1: QI < 0.1*  
P1= 9usec PL1= -2dB
```



NMR Experiment: Decoupling test ^{19}F (with ^{19}F GARP decoupling)
Basic Parameter Set: PROF19DEC
Pulse Program: zghfigqn

Sample: 0.05% Trifluorotoluene (TFT) in CDCl_3

Spin: on

Basic acquisition parameters:

d1 5s


ns 8 and ds 0

o1p=5ppm sw=10ppm

o2p=-63ppm


pl12 is the determined power level for garp decoupling

pcpd2=95us



Process with ft

Analysis: Use all necessary filters. Pulse according to formula $\text{pcpd} = 110 - (\text{BF1}_{1\text{H}}/20)$.
For instance, in a 300MHz: $\text{pcpd} = 110 - (300/20) = 95\text{us}$



NMR Experiment: B1 homogeneity test 1H from HWT (b1h)

Basic Parameter Set:

Pulse Program:

Sample: 0.1mg GdCl₃/ml D₂O +1% H₂O, 0.05% Meth.


Spin: off

Basic acquisition parameters:

d1 10s

ns 1 and ds 0

vdlist=syst1list




Process with ef (lb=1Hz)

Analysis: HWT; select HWTgen; select B1 homogeneity test H1.

The amplitude is measured at 810 degree as a percent level of the 90 degree pulse width

```
HWT Suite with xwin-nmr3.1 patchlevel 11; Probehead 1 : 5 mm EBI 1H-BB Z-CPD Z8202/0093
sysblhom - B1 homogeneity test H-1
Order Number: ZH022702; Instrument: av400; Sample: D20
100 experiments in 0.70 usec steps
Max[90] = 100.0%, Max[450] = 78.7%, Max[810] = 71.1%
```



BRUKER PULSE PROGRAM CATALOGUE


NMRGuide

BASIC 1D PULSE SEQUENCES

BASIC ACQUISITION PARAMETERS IN 1D

Once a parset has been loaded and getprosol optionally applied, the user can check and modify some acquisition parameters before start acquisition:

d1 is pre-scan delay
ns is the number of scans
o1p is the center of the spectrum
sw is the spectral width
td is the number of points to be acquired
aq is acquisition time
dw is dwell time
o2p is the decoupler frequency (optional)



BASIC PROCESSING PARAMETERS IN 1D

Processing parameters can be checked and modified from the ProcPars section:

si is the number of points to be processed (zero-filling)
wdw is the window function to be applied on the fid
ft: Fourier transformation
ef (em+ft): transform with a exponential using line broadening (**lb**)
gf (gm+ft): Transform with a Gaussian using **lb** and **gb**
apk is automatic phase correction
abs is automatic baseline correction
sr is the spectrum reference frequency

Basic 1D pulse sequences

- **Standard Experiments:**

Conventional ^1H spectrum (**zg30** / **zg** / **zg0** | PROTON)

Acquired as 2D (**zg2d**)

1D ^1H Homodecoupling (**zg0hd** / **zghd** / **zghd.2** | PROHOMODEC)

1D ^1H Band-selective homodecoupling (**zghc** / **zghc.2** / **zghc.3**)

NOEDIFF experiment:

Single irradiation (**zgf2pr**)

Using frequency list (**noediff** / **noediff.2** / **noedif.2** | NOEDIFF)

Irradiation multiplet frequencies within one multiplet (**noemul**)

^{13}C spectrum with selective ^1H decoupling using CW (**zgcw30** / **zgcw** / **zg0cw**)

^1H -decoupled ^{13}C spectrum (**zgd30** / **zgd** / **zg0dc** | c13cpd)

^1H -coupled ^{13}C spectrum (**zggd30** / **zggd** / **zg0gd** | c136d)

^1H -decoupled ^{13}C spectrum without NOE (**zgif30** / **zgif** / **zg0ig** | c1316)

^1H , ^{31}P -decoupled ^{13}C spectrum without NOE (**zgfbig**)

^1H , ^{31}P -decoupled ^{13}C spectrum without NOE in the channel f3 (**zgdcf2igf3**)

^1H , ^{31}P -decoupled ^{13}C spectrum without NOE (**zgif2igf3**)

Antiring sequence (**aring**, **aring2**)

1D sequence for suppression of background signals using composite pulse (**zgsb**)

UDEFT sequence with flip back pulse at the end of acquisition to decrease relaxation delay (**udeft**)

- **Multinuclear Applications:**

^{31}P -decoupled 1D ^1H spectrum (**zgif30** / **zgif** | PROP31DEC)

^{11}B -decoupled 1D ^1H spectrum (**zgif30** / **zgif** | PROB11DEC)

^1H -decoupled ^{15}N spectrum without NOE (**zgif** / **zgf3ig** | N1516)

^1H -coupled ^{15}N spectrum without NOE (**zg** | N15)

^1H -decoupled ^{31}P spectrum (**zpg30** | P31CPD)

^1H -coupled ^{31}P spectrum (**zg30** | P31)

- **Standard BRUKER parameter sets available for other nuclei:**

1D ^{11}B spectrum (**zg** | B11Z6)

1D ^{17}O spectrum (**zg** | O17Z6)

1D ^{23}Na spectrum (**zg** | NA23Z6)

1D ^{27}Al spectrum (**zg** | AL27ND)

1D ^1H -decoupled ^{29}Si spectrum (**zgif** | Si2916)

1D ^{35}Cl spectrum (**zg** | CL35Z6)

1D ^{37}Cl spectrum (**zg** | CL37Z6)

1D ^{71}Ga spectrum (**zg** | GA71Z6)

1D ^{71}Ga spectrum (**zg** | SE77Z6)

1D ^{103}Rh spectrum (**zg** | RH103Z6)

1D ^{111}Cd spectrum (**zg** | CD111Z6)

1D ^{113}Cd spectrum (**zg** | CD113Z6)

1D ^1H -decoupled ^{119}Sn spectrum (**zgif** | SN11916)

1D ^{195}Pt spectrum (**zg** | PT195Z6)

1D ^{199}Hg spectrum (**zpg** | HG199CPD)

Basic NMR Elements:

Pre-Scan Delay

Read Pulse

Acquisition Period


Phase Cycle

```

1 ze
2 30m
  d1
  p1 ph1
  go=2 ph31
  30m mc #0 to 2 F0(zd)
exit

ph1=0 2 2 0 1 3 3 1
ph31=0 2 2 0 1 3 3 1





```



- A) Pre-Scan Delay: Represented by **d1** (in seconds), is a period to allow relaxation to the original z position.
- B) Read Pulse: Represented by **p1** (in microseconds and applied at p1 power) is the time needed to create transverse magnetization from the original Iz magnetization.
- C) Acquisition Period: represented by **go** (the duration is **aq** in seconds) is the time required for the receiver to monitor the free-induction decay in the transverse plane.
- D) A phase Cycle is always included in all pulse programs to specify from which axis pulses and receivers are applied:


- 0 is x axis
- 1 is y axis
- 2 is -x axis
- 3 is -y axis

A Cyclops or Exorcycle loop is usually applied for better signal selection and artefact suppression



¹H NMR spectrum Analysis from: Chemical shifts, Multiplicities, Coupling Constants, Linewidths and Integration

Double-Resonance Experiments



Basic NMR Elements: Broadband Decoupling

```

...
1 ze
  d11 p112:f2
2 30m do:f2
  d11 cpd2:f2
  d1
  p1 ph1
  go=2 ph31
  30m do:f2 mc #0 to 2 F0(zd)
exit
....

;p112: f2 channel - power level for CPD/BB decoupling
;cpd2: decoupling according to sequence defined by cpdprg2
;pcpd2: f2 channel - 90 degree pulse for decoupling sequence

```



Selective Decoupling

```










...
1 ze
  d11 p126:f2
2 30m do:f2
  d11 cw:f2
  d1
  p1 ph1
  go=2 ph31
  30m do:f2 mc #0 to 2 F0(zd)
exit
....

;p126: f2 channel - power level for cw decoupling


```





13C NMR spectrum Analysis from Chemical shifts.



Selective irradiation with homodecoupling (see DIGMOD in eda)





zgdchd



Expansion of the homodecoupling spectrum of strychnine after irradiate at 1.29ppm. Note that J-coupled protons appear simplified.


Selective irradiation with CW




Also see: solvent suppression
(zgf2pr)

Miscellaneous Experiments


zgbs



aring




aring2



UDEFT:

M. Piotto, M. Bourdonneau, K. Elbayed, J.-M. Wieruszkeski & G. Lippens, *Magn. Reson. Chem.* 44, 943-947 (2006)

udeft



NMR Element: Purge Element Before d1

....


```
d12 p110:f1
p17 ph3
p17*2 ph4
d1 p11:f1
```

....

```
ph3=0
ph4=1
```

```
;p110: f1 channel - power level for TOCSY-spinlock
;p0 : f1 channel - 20 to 90 degree high power pulse
;p1 : f1 channel - 90 degree high power pulse
;p17: f1 channel - trim pulse
```

[2.5 msec]



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

T1 & T2 RELAXATION


Experiment Description:




T1 and T2 relaxation measurements are usually monitored by means of a series of 1D proton spectra acquired using inversion-recovery or CPMG pulse schemes, respectively.

T1 & T2 Relaxation

- ^1H T₁ measurements:
 - As 1D acquisition (**t1ir1d**)
 - As 2D acquisition (**t1ir || PROTON T1**)
- T₁ ^{13}C measurements (**t1irpg**) ^1H T₂ measurements
 - As 1D acquisition (**cpmg1d**)
 - As 2D acquisition (**cpmg**)
 - As 1D acquisition with presaturation (**cpmgpr1d**)
 - As 1D acquisition with presaturation and purge element (**zgpurge**)


Also see:
2D HSQC and 3D HNCO experiments for Backbone Dynamics in Proteins





1D Inversion-Recovery spectrum (top) of strychnine acquired with pulprog=t1ir1d (d7=0.5ms) showing the faster relaxation of methylenic CH₂ protons (resonances up) with respect to CH ones (resonances down)


T1 measurements:
The inversion-recovery experiment




$$dM_z/dt = R_1(M_{eq} - M_z) = (M_{eq} - M_z)/T_1$$

$R_1 = 1/T_1$ Hz is the longitudinal relaxational rate constant
 T_1 (in seconds) is the longitudinal relaxation time

$$I_z \xrightarrow{180^\circ I_x} -I_z \xrightarrow{d7} -I_z \exp(-d7/T_1) \xrightarrow{90^\circ I_x} I_y \exp(-d7/T_1)$$






A series of 1D Inversion-Recovery spectra of brucine automatically acquired with pulprog=t1ir and with a predefined delay list. 1D spectra can be automatically processed and then analyzed with the relaxation T1/T2 tool (**t1guide**) included into TOPSPIN. Otherwise, a simplified calculation can be done by detecting the so-called null point.

$$I_z = I_0(1 - 2\exp(-d/T_1))$$

↓

$$\ln(I_0 - I_z) = \ln(2I_0) - d/T_1$$

$$t_{\text{null}} = T_1 * \ln 2$$




$$dM_{xy}/dt = -R_2 M_{xy} = -M_{xy}/T_2$$

$R_2 = 1/T_2$ Hz is the transverse relaxation rate constant

T_2 (in seconds) is the spin-spin transverse relaxation time

T2 filter using CPMG (Carr-Purcell-Meiboom-Gill) sequence




NMR Building Block: A T2 Filter.

```

.....
3 d20
  p2 ph1
  d20
  lo to 3 times l4
.....

;d20: fixed echo time to allow elimination of J-mod. effects
;      d20 should be << 1/J ,but > (50 * P2)      [1-2 msec]
;l4: loop for T2 filter                          [4 - 20]
    
```




NMR Building Block: A T1(rho) Filter.

```


"p29=d29"

.....
4u p129:f1
  (p29 ph1)
.....

;p29: f1 channel - trim pulse
;d29: spinlock time                          [10 - 50 msec]
    
```



zgpurge



BRUKER PULSE PROGRAM CATALOGUE


NMRGuide

SELECTIVE 1D EXPERIMENTS

Experiment Description

For small molecules, it is usually only necessary to obtain specific information. For this reason, it is not necessary to record a time-consuming 2D dataset. A general approach is to convert the 2D experiment to a 1D counterpart by removing the variable t_1 evolution period and introducing some selective excitation element.

Selective excitation can be achieved by a selective 90 pulse (selzg) or by a selective gradient-based echo block that uses a selective 180 pulse (selgpse).



Sample Requirements

Selective experiments can be recorded on any type of sample. The only requirement is a well isolated resonance.

Hardware Requirements

In principle, selective experiments can be recorded on any probehead and on any modern spectrometer. An optional pulsed-field gradient coil (highly recommended) is required for gradient-based versions.

Set-up

For classical selective phase-cycled experiments:

p11 is the duration of the selective 90 pulse (in microseconds)

sp1 is the power level

spnam1 is the pulse shape (for instance, Gauss1.001)

spoffs1 is the offset (in Hz) to irradiate with respect to the central ω_1 frequency.

For gradient-based selective experiments:

p12 is the duration of the selective 180 pulse (in microseconds)

sp2 is the power level

spnam2 is the pulse shape (for instance, Gauss1.001)

spoffs2 is the offset (in Hz) to irradiate with respect to the central ω_1 frequency.

The result is a 1D spectrum equivalent to the specific row of the 2D analog. The advantages are: i) Acquisition is faster; ii) The resulting 1D spectrum offer better high-resolution properties than 2D, and iii) The experiment can be individually optimized.

Selective Excitation & 1D Selective Homonuclear experiments

- **Phase-Cycled:**

Using a shaped 90° pulse (**selzg** | **SELZG1H**)

Selective 1D COSY experiment (**selco** | **SELCO1H**)

Selective 1D RELAY experiment (**selcorl**)

Selective 1D TOCSY experiment (**selmlzf** | **SELMLZF1H**)

Selective 1D NOESY experiment

Using d8 random variation (**selno** | **SELNO1H**)

Using d8 and z-filter (**selnoz**)

Selective 1D ROESY experiment (**selro** | **SELRO1H**)

DANTE Excitation

2-2-6-2-2 DANTE-z scheme (**dazzg**)

3-6-3 DANTE-z scheme (**daz363zg**)

1-1 DANTE-z scheme (**daz11zg**)

- **Gradient-based using SPFGGE:**

Using selective pulsed-field-gradient spin-echo or SPFGGE (**selgpse** | **SELGPSE**)

Selective ge-1D COSY experiment (**selcogp** | **SELCPG0**)

Selective ge-1D TOCSY experiment:

using MLEV (**selmlgp** | **SELMLGP**)

using MLEV with ZQ suppression (**selmlgp.2**)

using DIPSI-2 (**seldigp**)

using DIPSI-2 with ZQ suppression (**seldigp.2**)

Selective ge-1D NOESY experiment (**selnogp** | **SELNOGP**)

with ZQ suppression (**selnogpzs**)

Selective ge-1D ROESY experiment (**selrogp** | **SELROGP**)

Selective ge-1D T-ROESY experiment (**selrogp.2**)

- **Gradient-based using Chemical Shift Selective Filter (CSSF):**

Using selective echo (**selcssf**)

Using selective echo and ZQ suppression (**selcssfzs**)

Using double selective echo (**selcssf.2**)

Using double selective echo and ZQ suppression (**selcssfzs.2**)

Selective ge-1D TOCSY experiment:

using DIPSI-2, single echo CSSF and ZQ suppression (**selcssfdizs**)


using DIPSI-2, double echo CSSF and ZQ suppression (**selcssfdizs.2**)

Selective ge-1D NOESY experiment


using single echo CSSF and ZQ suppression (**selcssfnzs**)

using double echo CSSF and ZQ suppression (**selcssfnzs.2**)


General Scheme Phase-cycled Selective



selzg




selco




selcor1




selmlzf




selro




selno



selnoz1




dazzg
daz11zg
daz363zg



NMR Building Block: Rectangular (Hard) vs Shaped (selective) Pulse


```

.....
4u p11:f1
p1 ph1
.....
  
```



```

.....
4u p10:f1
p11:sp1:f1 ph1:r
.....
  
```




General Scheme Gradient-Based SPFG

NMR Building Block: SPFG=Gradient-180sel-Gradient


```

50u UNBLKGRAD
.....
p16:gp1
d16 p10:f1
p12:sp2:f1 ph2:r
p16:gp1
d16 BLKGRAD
...
  
```




Common Parameters of Mixing Processes:
 p15 is the mixing time in ROESY experiments
 d8 is the mixing time in NOESY experiments
 d9 is the mixing time in TOCSY experiments
 d4 is the evolution time in COSY experiments

NMR Experiment: Calibration selective 180° shaped pulse
Parameter Set: SELGPSE
Pulse Program: selgpse



1. Record a ¹H spectrum
2. Define "pulprog selgpse"
3. Set o1p on the selected resonance
4. Set p12=20m
5. Set spnam2=Gauss1.1000
6. Run **paropt** varying **sp2** starting from 70 dB and with an increment of -2dB
7. The first maximum is the selective 180° pulse of 20ms of duration





Remember Linear Amplifiers:
Selectivity: $p12*2$ and $sp2+6dB$



- If the selective 180 pulse of 20 ms of duration needs 60 dB:
- The selective 90 pulse of 10 ms is 60 dB
 - The selective 180 pulse of 10 ms is 54dB
 - The selective 180 pulse of 40 ms is 66 dB




Selectivity in selective experiments:

Depends how isolated is the selected resonance. As a good starting conditions. A selective 180 pulse (spnam2=Gauss1.1000) of 20000 microseconds (p12) affords an effective field around 60Hz. More selectivity is achieved by using longer pulses, at expense of the presence of unwanted relaxation effects.







Use the shape tool (**stdisp**) to generate new shapes and to calculate a different number of parameters and features related to them.







Suppression of ZQC in selective 1D TOCSY experiments:




Selective ROESY vs NOESY in medium-sized molecules in a 500MHz




selcssf




selcssf.2




selcssfzs




selcssfzs.2





selcssfdizs



selcssfdizs.2





Selective Excitation & 1D Selective Heteronuclear experiments

¹³C Selective excitation using a shaped 90° pulse (**selzpgg**)

Selective 1D X-X COSY experiment (**selcpgg**)

Selective 1D INADEQUATE experiment (**selina**)

Selective ge-1D HSQC experiment (**selhsqcgpsp**)

Selective ge-1D HSQC experiment with PEP (**selhsqcgpsisp**)


Selective ge-1D Long-range HSQC experiment (HSQMBC) (**selhsqcgplrndsp**)

Selective ge-1D HSQC-NOESY experiment (**selhsqcgpnosp**)


Selective ge-1D HSQC-NOESY experiment without decoupling (**selhsqcgpndnosp**)

Selective 1D INEPT (INAPT) experiment (**selineptlrrdsp**)


selzpgg




selcpgg




selina




selhsqcgppsp





selhsqcgppsisp





selhsqcgplmdsp





selhsqcgpnosp





A complete Description of creating, analyzing and Manipulating
 RF and Gradient Shapes in the Shape Tool Manual
 (see **Help** menu in Topspin)

BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

SOLVENT SUPPRESSION

1D Solvent suppression

- **Classical:**

1D solvent presaturation:

- Conventional (**zgpr** / **zg0pr** | **ZGPR**)
- Using composite pulses (**zgcppr** | **ZGCCPR**)
- Using spoil gradient (**zggppr**)
- Using composite pulse and spoil gradient (**zgcpgppr**)
- From f2 channel (**zgf2pr** / **zg0f2pr**)
- Using shaped pulse for off-resonance presaturation (**zgps**)

Jump and return:

- 1-1 scheme (**p11**)
- 1-3-3-1 scheme (**p1331**)

- **Gradient-based:**

1D WATERGATE:

- Using 3-9-19 scheme (**p3919gp** | **P3919GP**)
- Using 3-9-19 and flip back pulse (**p3919fpgp**)
- Using 90° water-selective pulses (**zggpwg** | **ZGPPWG**)

1D Excitation Sculpting:

- Using 180° water-selective pulses (**zgesgp** | **ZGESGP**)
- Using 180° water-selective and flip back pulse (**zgesfpgp**)
- Using W5 pulse train (**zggpw5**)

1D WET Scheme:

- Conventional (**wet**)
- With ¹³C decoupling on f2 during WET and AQ (**wetdc** | **LC1BWTDc**)
- With ¹³C decoupling on f2 during WET (**wetdw**)
- With shape pulse and C-13 decoupling on f2 during WET and acquisition for LC isocratic runs (**lc2wetdc**)
- With shape pulse and C-13 decoupling on f2 during WET and AQ with intermediate preparation scan into second dataset for LC gradient runs with updated shapes (**lc2wetdcus** | **LC2BWTUS**)
- With shape pulse and C-13 decoupling on f2 during WET and acquisition with intermediate preparation scan into second dataset for LC gradient runs with updated shapes (**lc2grdonflow**)

Related Experiments:



- All these 1D experiments can be incorporated in any multidimensional NMR experiment. Please refer to each chapter to check the different possibilities for 2D and 3D solvent-suppressed experiments
- LC-NMR Experiments

NMR Building Block: Solvent Presaturation



```

...
d12 p19:f1
d1 cw:f1 ph29
d13 do:f1
d12 p11:f1
...


ph29=0
.....
;p19 : f1 channel - power
; level for presaturation
  
```



Also see: LC-NMR experiments




Also see: noediff



WATERGATE:

1. M. Piotto, V. Saudek & V. Sklenar, J. Biomol. NMR 2, 661 - 666 (1992)
2. V. Sklenar, M. Piotto, R. Leppik & V. Saudek, J. Magn. Reson., Series A 102, 241 -245 (1993)



NMR Element: WATERGATE (3-9-19)

...

```

50u UNBLKGRAD
p16:gp1
d16 p118:f1
p27*0.231 ph4
d19*2
p27*0.692 ph4
d19*2
p27*1.462 ph4
d19*2
p27*1.462 ph5
d19*2
p27*0.692 ph5
d19*2
p0*0.231 ph5
50u
p16:gp1
d16
4u BLKGRAD

```

.....

```

ph4=0
ph5=2

```

.....

```

;p11 : f1 channel - power level for pulse (default)
;p118: f1 channel - power level for 3-9-19-pulse (watergate)
;p0 : f1 channel - 90 degree pulse at p118
;
; use for fine adjustment
;p27: f1 channel - 90 degree pulse at p118
;d16: delay for homospoil/gradient recovery
;d19: delay for binomial water suppression
; d19 = (1/(2*d)), d = distance of next null (in Hz)

```


NMR Building Block: WATERGATE (water-selective 90°)

```

...
50u UNBLKGRAD
p16:gp1
d16 p10:f1
(p11:sp1 ph2:r):f1
4u
d12 p11:f1
(p2 ph3)
4u
d12 p10:f1
(p11:sp1 ph2:r):f1
46u
p16:gp1
d16
4u BLKGRAD
...


ph2=0
ph3=2

;p10 : 120dB
;sp1 : f1 channel - shaped pulse 90 degree
;p2 : f1 channel - 180 degree high power pulse
;p11: f1 channel - 90 degree shaped pulse
;p16: homospoil/gradient pulse
;d16: delay for homospoil/gradient recovery
    
```



1D Solvent suppression methods in a doubly-labeled ubiquitin sample dissolved in 95% H_2O /5% D_2O

M. Liu, X. Mao, C. He, H. Huang, J.K. Nicholson & J.C. Lindon, J. Magn. Reson. 132, 125 - 129 (1998)



NMR Element: Excitation Sculpting (W5)

```

...
50u UNBLKGRAD
p16:gp1
d16 p118:f1
p27*0.087 ph3
d19*2
p27*0.206 ph3
d19*2
p27*0.413 ph3
d19*2
p27*0.778 ph3
d19*2
p27*1.491 ph3
d19*2
p27*1.491 ph4
d19*2
p27*0.778 ph4
d19*2
p27*0.413 ph4
d19*2
p27*0.206 ph4
d19*2
p27*0.087 ph4
50u
p16:gp1
d16
4u


p16:gp2
d16
p27*0.087 ph5
d19*2
p27*0.206 ph5
d19*2
p27*0.413 ph5
d19*2
p27*0.778 ph5
d19*2
p27*1.491 ph5
d19*2
p27*1.491 ph6
d19*2
p27*0.778 ph6
d19*2
p27*0.413 ph6
d19*2
p27*0.206 ph6
d19*2
p27*0.087 ph6
p16:gp2
d16
.....
ph3=0
ph4=2
ph5=0
ph6=2
  
```

```

;p118: f1 channel - power level for 3-9-19-pulse (watergate)
;p16: homospoil/gradient pulse
;p27: f1 channel - 90 degree pulse at p118
;d16: delay for homospoil/gradient recovery
;d19: delay for binomial water suppression
;      d19 = (1/(2*d)), d = distance of next null (in Hz)

;gpz1: 34%
;gpz2: 22%
  
```


T.-L. Hwang & A.J. Shaka, *J. Magn. Reson., Series A* 112 275-279 (1995)



NMR Element: Excitation Sculpting (water-selective 180°)

```

...
50u UNBLKGRAD
p16:gp1
d16 p10:f1
(p12:sp1 ph2:r):f1
4u
d12 p11:f1

p2 ph3

4u
p16:gp1
d16
TAU
p16:gp2
d16 p10:f1
(p12:sp1 ph4:r):f1
4u
d12 p11:f1

p2 ph5

4u
p16:gp2
d16

....
ph2=0
ph3=2
ph4=0
ph5=2
.....

;p10 : 120dB
;sp1 : f1 channel - shaped pulse 180 degree
;p2 : f1 channel - 180 degree high power pulse
;p12: f1 channel - 180 degree shaped pulse (Squa100.1000) [2 msec]
;p16: homospoil/gradient pulse
;d16: delay for homospoil/gradient recovery

;gpz1: 31%
;gpz2: 11%

```

```

...
50u UNBLKGRAD
d12 p10:f1
d13
(p11:sp7 ph5):f1
4u
p16:gp21
d16 p10:f1
(p11:sp8 ph6):f1
4u
p16:gp22
d16 p10:f1
(p11:sp9 ph6):f1
4u
p16:gp23
d16 p10:f1
(p11:sp10 ph6):f1
4u
p16:gp24
d16
  
```




```

.....
ph5=0
ph6=1
.....
  
```

```

;p10 : 120 dB
;sp1: f1 channel - shaped pulse
;      powerlevel for 90 degree pulse
;sp7: f1 channel - shaped pulse
;sp8: f1 channel - shaped pulse
;sp9: f1 channel - shaped pulse
;sp10: f1 channel - shaped pulse
;p11: f1 channel - 90 degree shaped pulse
;p16: homospoil/gradient pulse [2msec]
;d16: delay for homospoil/gradient recovery
;gpz21: 80%
;gpz22: 40%
;gpz23: 20%
;gpz24: 10%
;use power levels
;
           sp7      sp8      sp9      sp10
           +0.87    -1.04    +2.27    -5.05
  
```


NMR Building Block: WET




How to generate a multi-frequency shape.

Step 1: Open the stdisp tool and select the appropriate shape (e.g: Gauss1.1000)


Step 2: Select Phase Modulation according offset in the Manipulate menu:




Step 3: Define type of manipulation and number of frequencies to be excited:



Step 4: Define frequencies:



Step 5: Store the resulting shape. It can be used in any selective-excitation or solvent-suppression pulse program using shaped (sp) pulses for multiexcitation (Hadamard spectroscopy) or multiple-solvent suppression (HPLC-NMR) purposes



A complete Description of creating, analyzing and Manipulating RF and Gradient Shapes in the Shape Tool Manual (see **Help** menu in Topspin)

BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

SINGLE/MULTIPLE SUPPRESSION LC-NMR EXPERIMENTS

Single-Multiple Presaturation LC-NMR Experiments

1D ¹H spectrum

- 1D ¹H with double presaturation (**lc1prf2** | **LC1D12**)
- 1D ¹H with triple presaturation (**lc1prft**)
- 1D ¹H with presaturation using shaped pulses, composite pulses and CW decoupling on f2 during acquisition (**lclcpwps**)

1D NOESY

- 1D NOESY with presaturation (**noesypr1d**)
- 1D NOESY with presaturation and CW decoupling on f2 (**lc1pncw**)
- 1D NOESY with double presaturation and CW decoupling on f2 (**lc1pncwfd**)
- 1D NOESY with presaturation using shaped pulse and CW decoupling on f2 (**lc1pncwps**)
- 1D NOESY with double presaturation (**lc1pnf2**)
- 1D NOESY with multiple presaturation (**lc1pnfr**)
- 1D NOESY with triple presaturation (**lc1pnft**)
- 1D NOESY with presaturation using shaped pulse (**lc1pnps**)

Pseudo-2D-sequence

- Pseudo-2D-sequence for lc-nmr on flow detection (**lc2**)
- Pseudo-2D-sequence for lc-nmr on flow detection with power-gated decoupling (**lc2pg**)
- Pseudo-2D-sequence for lc-nmr on flow detection with presaturation (**lc2pn**)
- Pseudo-2D-sequence for lc-nmr on flow detection with double presaturation (**lc2pnf2**)
- Pseudo-2D-sequence for lc-nmr on flow detection with solvent gradients (**lc2pnf2ul**)
- Pseudo-2D-sequence for lc-nmr on flow detection with solvent gradients (**lc2pnpl**)
- Pseudo-2D-sequence for lc-nmr on flow detection (**lc2pnps**)
- Pseudo-2D-sequence for lc-nmr on flow detection with solvent gradients (**lc2pnul**)
- Pseudo-2D-sequence for lc-nmr on flow detection with presaturation (**lc2pr**)
- Pseudo-2D-sequence for lc-nmr on flow detection with double presaturation (**lc2prf2**)
- Pseudo-2D-sequence for lc-nmr on flow detection with presaturation using shape pulse (**lc2ps**)

Also see in 1D solvent suppression methods

1D ¹H spectrum using WET solvent suppression

- 1D ¹H with WET (**wet**)
- 1D ¹H with WET and CW decoupling on f2 during WET and acquisition (**wetdc** | **LC1DWTDC**)
- 1D ¹H with WET and CW decoupling on f2 during WET (**wetdw**)
- 1D ¹H WET solvent suppression with shape pulse and C-13 decoupling on f2 during WET and acquisition for LC isocratic runs (**lc2wetdc**)
- 1D ¹H WET solvent suppression with shape pulse and C-13 decoupling on f2 during WET and AQ with intermediate preparation scan into second dataset for LC gradient runs with updated shapes (**lc2wetdcus** | **LC2DWTUS**)
- 1D ¹H WET solvent suppression with shape pulse and C-13 decoupling on f2 during WET and acquisition with intermediate preparation scan into second dataset for LC gradient runs with updated shapes (**lc2grdonflow**)

Also see LC-NMR related pulse programs in:

2D homonuclear J-resolved

- 2D J-resolved with double presaturation and cw-decoupling on f2 (lcjrescwf dprqf)
- 2D J-resolved with presaturation and cw-decoupling on f2 (lcjrescwpqrqf)
- 2D J-resolved with presaturation using shape pulse and cw-decoupling on f2 (lcjrescwpqsqf)
- 2D J-resolved with double presaturation (lcjresf2pqrqf)
- 2D J-resolved with presaturation (lcjrespqrqf)
- 2D J-resolved with presaturation using shape pulse (lcjrespsqf)

2D TOCSY

- 2D TOCSY with double presaturation and cw-decoupling on f2 (lcmlevcwf dpcph)
- 2D TOCSY with presaturation and cw-decoupling on f2 (lcmlevcwp cphps)
- 2D TOCSY with double presaturation using composite pulse (lcmlevf2pcph)
- 2D TOCSY with double presaturation (lcmlevf2p hpr | LCM12)
- 2D TOCSY with presaturation using shape pulse and composite pulse (lcmlevpcphps)
- 2D TOCSY with presaturation using composite pulse (lcmlevpcph)

2D Experiments using WET

- 2D COSY using WET (cosydcphwt | COSYDCPHWT)
- 2D TOCSY using WET (mlevdcphwt | MLEVDCPHWT)
- 2D HSQC using WET (hsqcetgpsiw t | HSQCETGPSIWT)


2D Experiments using single/multiple presaturation using shape pulse

- Phase-sensitive 2D COSY using using single/multiple presaturation (cosycwphps | COSYCWP HPS)
- Phase-sensitive 2D HSQC using using single/multiple presaturation (hsqcphps)
- 2D HMBC using using single/multiple presaturation (hmbcndpsqf)


Related experiments:

- Also see 1D Solvent suppression

lc2




lc2pr




lc2ps




lc2pg




lc1prf2 lc2prf2










lc1prft



lc1cpcwps





A Description of setting-up and running LC-NMR Experiment can be found in the [LC-NMR Manual](#) (see Help menu in Topspin)

List of Automation Protocols to run in LC-NMR Applications:

1. expsetup: Automation Setup
2. getinfo: Get information into the title.
3. multicmd: Perform Multiple Commands
4. au_lc1d: 1D spectra with solvent suppression
5. au_lc2d: 2D spectra with solvent suppression
6. au_lconfow: Onflow spectra with solvent suppression
7. lcprep: Acquire preparation spectrum
8. lcsetup: Search for solvents
9. lcabsf: Baseline Correction
10. lcsino: Signal-to-noise Calculation
11. proc_onflow: 1D & 2D processing
12. lctshim: Shimming with TOPSHIM
13. lcgshim: Shimming with gradshimau
14. lczshim: Iterative Shimming of the lineshape

BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

¹⁹F SPECIFIC EXPERIMENTS

^{19}F NMR Experiments

- **1D ^{19}F experiments:**

^1H -decoupled ^{19}F spectrum (**zgfhighq**n / **zgfhighq**n.2 | F19CPD)

^1H -coupled ^{19}F spectrum (**zgf1q**n | F19)

Using 19F lockswitch unit (**zg19f**)

^{19}F -homodecoupled ^{19}F spectrum (**zhflhdq**n)

^{19}F -decoupled 1D ^1H spectrum (**zghfigq**n / **zghfigq**n30 | PROF19DEC)


^1H spectrum with ^{19}F -presaturation (**zgf2hfpr**)

- **2D ^{19}F experiments:**


2D ^{19}F - ^1H HETCOR experiment (**hfcoqf**q)

2D ^{19}F - ^1H HOESY experiment (**hoesyfhqf**qrv)


zg19f
zgf1qn





zgfhighqn.2
zgfhighqn





zhflhdqn








zgf2hfpr



hfcoqfqm



hoesyfhqfqmrv



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

^2H SPECIFIC EXPERIMENTS

²H Experiments

- **1D spectra:**

1D ²H spectrum (**zg2h**)

Using ²H lockswitch unit (**zg2h.2**)

1D X-decoupled ²H spectrum (**zgig2h, zgig2hf4**)

- **2D spectra:**







Magnitude-mode 2D HETCOR with ²H-decoupling (**hxcoqf2h**)

- **Miscellaneous:**

High-power 90° ²H decouple pulse calibration (**decp902h, decp902hf4**)

Also see:

2H-decoupled 3D triple-resonance experiments



A Description of Gradient-Shimming interface and the User Manual for the Automatic Shimming can be found in the **Gradient Shimming Manual** and **Topshim Manual**, respectively (see **Help** menu in Topspin)

BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

BASIC 1D GRADIENTS

Basic 1D Gradients

- **Standard:**

Gradient-enhanced 1D Echo experiment (*zggeg*)
Gradient-enhanced 1D Spin-Echo experiment (*zggpse*)

- **Gradient Calibration:**

Gradient Strength Calibration (*calibgp*)
Gradient Preemphasis Adjustment. Gradient Recovery Test (*preempgp2*)

- **Gradient shimming:**

1D Gradient Echo for gradshim-procedure (*imgegp1d*)
Using 2H lockswitch unit (*imgegp1d2h*)
Using 19F lockswitch unit (*imgegp1d19f*)
Using selective pulse (*imgegpsp1d*)

1D convection-compensated Gradient Echo for gradshim-procedure (*imgegpvcv1d*)
Using 2H lockswitch unit (*imgegpvcvsp1d2h*)

2D Gradient Echo for gradshim-procedure (*imgegp2d*)

3D Gradient Echo for gradshim-procedure (*imgegp3d*)
Using BSMS RCB board(*imrcbgegp3d*)

- **Zero-Quantum Gradient Calibration:**


ZQ setup calibration (*zs_setup*)

Basic Syntax for a Gradient:


```
#include <Grad.incl>

....
50u UNBLKGRAD
p16:gp1
d16
20u BLKGRAD
....


;p16: homospoil/gradient pulse
;d16: delay for homospoil/gradient recovery
;gpz1: gradient strength in %
;gpnam1: SMSQ10.100
```



imgegp2d imgegp1d
 imgegp3d imgegp1d19f
 imgegp1d2h




imgegpsp1d




imgegpspcv1d

if ZGOPTNS= -DLABEL_NOCOMP




else




imgegpspcv1d2h


if ZGOPTNS= -DLABEL_NOCOMP



else



zs_setup








BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

1D ERETIC SEQUENCES

ERETIC 1D pulse sequences

- 1D ERETIC for quantitative measurements and using f2 channel for ERETIC signal
 - Conventional 1D (**eretic**)
 - Conventional 1D using read pulse of 30° (**eretic30**)
 - Conventional using presaturation (**ereticpr**)
 - Conventional using presaturation and read pulse of 30° (**ereticpr30**)
 - Using 1D NOESY and preaturation (**ereticgppn**)




BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

GENERAL SCHEMES FOR
2D/3D/4D EXPERIMENTS

NMR Pulse Sequence: Definition of Time Periods

A NMR pulse sequence can be splitted in several different and independent parts, namely, preparation, evolution, mixing and acquisition periods:




The relative sensitivity of a given NMR pulse scheme depends basically of two factors: i) the starting nucleus and ii) the detected nucleus as a function of the following relationship:

$$S/N \propto \gamma_{\text{excited}} \gamma_{\text{detected}}^{3/2}$$


Thus, due to the highest γ value for ^1H , it is generally preferable to start from ^1H and to detect ^1H when possible.

GENERAL 2D HOMONUCLEAR SCHEMES


1. Phase-Cycled homonuclear



2. Gradient-Enhanced homonuclear




3. Solvent-Suppressed homonuclear




The preparation period can also include a solvent-suppression element

HOW 2D NMR WORKS?


Consider the most basic 1D experiment:




A target spin 1 evolves during t_1 :



Imagine an experiment in which the information of spin 1 is transferred to spin 2




The resulting 2Dmap looks like this




HOW 2D NMR WORKS?: COSY

The experiment is repeated td times using a variable period with an initial d_0 value and an increment of in_0 , giving the maximum value of $d_0 * td$




GENERAL 2D 1H-X HETERONUCLEAR


1. Old X-detected (HETCOR, COLOC)



2. Phase-Cycled Inverse (old HMQC, old HSQC...)




3. Gradient-Enhanced Inverse (Gradient-Enhanced HMQC, HSQC..)




Solvent-suppression elements can be included
in both preparation and mixing periods



GENERAL 2D HYBRID INVERSE





GENERAL 2D RELAXATION



GENERAL 3D DOUBLE-RESONANCE




Example: 3D HSQC-TOCSY experiment




Example: 3D TOCSY-HSQC experiment

GENERAL 3D TRIPLE-RESONANCE

Out-and-back
 (Example: 3D HNCO)




Out-and-stay




ANALYZING 3D SPECTRA

The same as 2D NMR but with two independent variable evolution periods. The first planes can be recorded from the same dataset by setting the number of TD in each indirect dimension to 1.



GENERAL 4D

4 different frequencies can be sampled during the same experiment



In principle, nD spectroscopy is possible by introducing (n-1) variable evolution periods.
Major advantage: Excellent dispersion and maximum information from a single spectrum.

Major inconvenients are:

- i) the sequence becomes very long and sensitivity losses due to fast transverse relaxation is usually a great problem;
- ii) The overall experimental acquisition time is increased.

Some solutions:

1. Combination of several 3D experiments
2. Use of time-sharing evolution approach: Reduced dimensionality ...

Classification of NMR Experiments. The 3 key points:



J_{HH} NOE

COSY NOESY
RELAY ROESY
TOCSY (EXSY)

Initial Nucleus: Determine d1
Indirect Evolving Nucleus: Determine F1 dimension
Detected Nucleus: Determine F2 dimension



$^1J_{CH}$ NOE

DEPT/INEPT
HETCOR HOESY



$^nJ_{CH}$

COLOC



$^1J_{CH}$ HMQC
HSQC
Edited HSQC
TROSY



$^nJ_{CH}$
HMBC



$^1J_{CH} + J_{HH}$

HSQC-COSY
HSQC-TOCSY

$^1J_{CH} + \text{NOE}$

HSQC-NOESY
HSQC-ROESY



$J_{HH} + ^1J_{CH}$

TOCSY-HSQC

$\text{NOE} + ^1J_{CH}$

NOESY-HSQC



$^1J_{CC}$

INADEQUATE



$^1J_{CH} + ^1J_{CC}$

ADEQUATE

BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D COSY EXPERIMENTS

Experiment Description

A COSY (CORrelation Spectroscopy) experiment is a high-sensitive tool that provides through-bond proton-proton connectivities by means of scalar J coupling constants.

Sample Requirements

COSY experiments can be recorded on any type of sample. Solvent-suppressed versions are required for samples dissolved in H₂O.

Hardware Requirements

In principle, COSY experiments can be recorded on any probehead. An optional pulsed-field gradient coil (highly recommended) is required for gradient-based versions.

NMR Spectrum


The experiment yields a 2D proton-proton correlation map with two different types of signals: i) autocorrelation diagonal peaks, and ii) Off-diagonal cross-peaks correlating J-coupled spins. COSY spectra are usually represented in magnitude-mode, and a COSY-DQF experiment is recommended when phase-sensitive data is required,

References:

W.P. Aue, E. Bartholdi, R.R. Ernst, J. Chem. Phys. 64, 2229 (1976)
K. Nagayama et al., J. Magn. Reson. 40, 321 (1980)

Related Experiments

Selective 1D COSY
2D COSY-DQF version
2D Relayed and TOCSY experiments



2D COSY Experiments

- **Phase-cycled:**

Magnitude-mode 2D COSY (**cosyqf** | *cosy45sw* / *cosy90sw*)
Magnitude-mode 2D COSY using a 45 pulse (**cosyqf45** | *cosy45sw*)
Magnitude-mode 2D COSY using a 90 pulse (**cosyqf90** | *cosy90sw*)
Magnitude-mode 2D COSY using purge pulses before d1 (**cosyppqf**)
Phase-sensitive 2D COSY (**cosyph**)

Magnitude-mode Long-Range optimized 2D COSY (**cosylrqf**)

Constant-Time 2D COSY (**cosyjdqf**)

- **Phase-cycled and solvent suppression:**

Magnitude-mode 2D COSY with presaturation (**cosyprqf**)
Phase-sensitive 2D COSY with presaturation (**cosyphpr** | *cosyphpr*)

- **Gradient-based:**

Magnitude-mode ge-2D COSY (**cosygpqf** | *cosygpsw*)
Magnitude-mode ge-2D COSY using presaturation (**cosygpprqf**)
Magnitude-mode ge-2D COSY using purge pulses before d1 (**cosygpppqf**)
Phase-sensitive ge-2D COSY using echo-antiecho (**cosyetgp**)
Phase-sensitive 2D z-COSY using ZQ suppression (**cosygpzhfzs**)

- **Gradient-based and solvent suppression:**

Magnitude-mode ge-2D COSY using presaturation (**cosygpprqf**)
2D COSY using WET (**cosydcphwt** | *cosydcphwt*)

- **Miscellaneous:**

Phase-sensitive w_1 -region-selective 2D COSY (**sosyph**)
Phase-sensitive w_1 -region-selective 2D COSY with refocusing (**sosyphrd**)


Phase-sensitive 2D COSY with off-resonance single or multiple presaturation (**cosycwphps** | *cosycwphps*)
Magnitude-mode ge-2D COSY with off-resonance single or multiple presaturation (**cosycwgppsqf**)

Magnitude-mode 2D ^{13}C - ^{13}C COSY (**cosydcqf**)
Magnitude-mode long-range optimized 2D ^{13}C - ^{13}C COSY (**cosydlrqf**)
Phase-sensitive 2D ^{13}C - ^{13}C COSY (**cosydcph**)

NMR Building Block: A Variable Evolution Period.



```

"i0=inf1"
"d0=3u"
...
d0
...
d1 mc #0 to 2 F1QF(caldel(d0, +i0))
...
;d0 : incremented delay (2D)           [3 usec]
;i0 : 1/SW = 2 * DW
;i1 : 1/(1 * SW) = 2 * DW
;nd0: 1
    
```





$$I_{1z} \xrightarrow{90^\circ I_x} -I_{1y} \xrightarrow{2\pi\omega_1 t_1 I_z} -I_{1y} \cos \omega_1 t_1 + I_{1x} \sin \omega_1 t_1$$

$$2\pi J_{12} t_1 I_{1z} I_{2z}$$



$$-I_{1y} \cos \omega_1 t_1 \cos \pi J_{12} t_1 + 2I_{1x} I_{2z} \cos \omega_1 t_1 \sin \pi J_{12} t_1 + I_{1x} \sin \omega_1 t_1 \cos \pi J_{12} t_1 + 2I_{1y} I_{2z} \sin \omega_1 t_1 \sin \pi J_{12} t_1$$




$$90^\circ I_x$$


$-I_{1z} \cos \omega_1 t_1 \cos \pi J_{12} t_1$	see NOESY
$-2I_{1x} I_{2y} \cos \omega_1 t_1 \sin \pi J_{12} t_1$	see Double-Quantum Filtered COSY
$+I_{1x} \sin \omega_1 t_1 \cos \pi J_{12} t_1$	Autocorrelated (Diagonal) COSY
$-2I_{1z} I_{2y} \sin \omega_1 t_1 \sin \pi J_{12} t_1$	COSY cross-peak

Phase-Cycled COSY


cosyph
 cosyqf




cosyqf45




cosyqf90




cosyppqf




cosyphpr
 cosyprqf








cosylrqf




cosyjdqf






COSY with Gradient Selection




Region-Selective COSY




X-Decoupled COSY



COSY with WET Solvent Suppression



Homonuclear Mixing Times



COSY		
COSY-MQF		Also with Gradients
RELAY		Also with Multiple Relay Steps
TOCSY		Also with DIPSI-2 ...
ROESY		Also with T-ROESY ..
NOESY		Also with Purge Gradients

BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D COSY-DQF EXPERIMENTS

Experiment Description

A COSY-MQF experiment is a modified COSY pulse scheme that incorporates a multiple-quantum filter just prior acquisition.

It provides through-bond proton-proton connectivities by means of scalar J coupling constants

Sample Requirements

COSY-DQF experiments can be recorded on any type of sample. Solvent-suppressed versions are required for samples dissolved in H₂O.

Hardware Requirements

In principle, COSY-DQF experiments can be recorded on any probehead. An optional pulsed-field gradient coil is required for gradient-based versions.

NMR Spectrum


The experiment yields a 2D proton-proton correlation map with two different types of signals: i) autocorrelation diagonal peaks, and ii) Off-diagonal cross-peaks correlating J-coupled spins. Although COSY-DQF can be represented in magnitude-mode, phase-sensitive data is most suitable for analyzing the anti-phase character of cross-peaks.

Related Experiments

Selective 1D COSY
2D COSY experiment
2D Relayed and TOCSY experiments

References:

1. A.A. Shaw, C. Salaun, J.-F. Dauphin & B. Ancian, J. Magn. Reson. A 120, 110-115 (1996)
2. B. Ancian, I. Bourgeois, J.-F. Dauphin & A.A. Shaw, J. Magn. Reson. A 125, 348- 354 (1997)



2D COSY-MQF Experiments

- **Phase-cycled:**

Magnitude-mode 2D COSY with DQF (**cosydfqf**)
Magnitude-mode 2D COSY with TQF (**cosyqftf**)
Phase-sensitive 2D COSY with DQF (**cosydfph** | **cosyDQFPHSW**)
Phase-sensitive 2D COSY with DQF and purge pulse (**cosydfphpp**)
Phase-sensitive 2D COSY with TQF (**cosyphftf**)

Phase sensitive 2D E.COSY -KcMAX=3 (**ecos3nph**)
Complementary Phase sensitive 2D E.COSY - KcMAX=3 (**ecos3cph**)

- **Phase-cycled and solvent suppression:**

Phase-sensitive 2D COSY with DQF & presaturation (**cosydfphpr**)

- **Gradient-based:**

Magnitude-mode ge-2D COSY with multiple-quantum filter (**cosygpmf qf** | **cosyGPMFSW**)
Magnitude-mode ge-2D COSY with multiple-quantum filter and purge pulse (**cosygpmf ppqf**)
Phase-sensitive ge-2D COSY with multiple-quantum filter (**cosygpmf ph** | **cosyGPDPHSW**) Phase-sensitive ge-2D COSY with multiple-quantum filter and purge pulse (**cosygpmf phpp**)
Phase-sensitive ge-2D COSY with DQF using echo-antiecho (**cosydfetgp.1**)
Phase-sensitive ge-2D COSY with gradient-based DQF using echo-antiecho (**cosydfetgp.2**)
Phase-sensitive ge-2D COSY with gradient-based DQF using echo-antiecho and purge pulse (**cosydfetgppp.2**)


Gradient E.COSY (**ecosygpph**)

- **Gradient-based and solvent suppression:**


Phase-sensitive 2D COSY-DQF with WATERGATE using 3-9-19 (**cosydfgpph19**)
Phase-sensitive 2D COSY-DQF with Excitation Sculpting using 180 water-selective pulse (ES element) (**cosydfesgpph**)
Using purge pulses before d1 (**cosydfesgpphpp**)

Also See:
2D COSY Experiments


Phase-Cycled COSY-QMF




Gradient-Enhanced COSY-QMF



cosygpmpfqf




cosygpmpfpqf




16:12:40 for DQ filter
 4:12:40 for TQ filter

cosygpmpfph




cosygpmpfphpp




10:20 for DQ filter
 10:30 for TQ filter


cosydfetgp.1



cosydfetgp.2




cosydfetgppp.2




COSY-QMF with Gradient Solvent Suppression


cosydfgpph19



cosydfesgpph



cosydfesgpphpp



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D SECSY EXPERIMENTS

Experiment Description

A SECSY (Spin-Echo Correlated Spectroscopy) experiment is actually an obsolete tool that provides through-bond proton-proton connectivities by means of scalar J coupling constants in the same way as described for COSY experiment.

Sample Requirements

SECSY experiments can be recorded on any type of sample.

Hardware Requirements

SECSY experiments can be recorded on any probehead.

NMR Spectrum

The experiment yields a 2D correlation map in which the conventional diagonal axis is rotated to $F1=0$. J Cross-peaks are correlated about an axis of 135° relative to $F1=0$

Related Experiments

See 2D COSY experiments

References:

K. Nagayama et al., J. Magn. Reson. 40, 321 (1980)



2D SECSY Experiments

Magnitude-mode 2D SECSY (**secsyqf**)

Magnitude-mode long-range optimized 2D SECSY (**secsylrqf**)

Also See:

2D COSY and 2D COSY-DQF Experiments



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D RELAY EXPERIMENTS

Experiment Description

RELAY experiments are extensions of the conventional COSY experiment and provide through-bond proton-proton connectivities by means of scalar J coupling constant in a step-by-step way. Thus, one-step, two-step, or three-step RELAY versions are available for the analysis of individual subspin systems into a molecular framework.

Sample Requirements

RELAY experiments can be recorded on any type of sample.

Hardware Requirements

RELAY experiments can be recorded on any probehead.

NMR Spectrum

The experiment yields a 2D proton-proton correlation map with different types of signals: i) autocorrelation diagonal peaks, ii) Off-diagonal cross-peaks correlating directly J-coupled spins (COSY peaks), and iii) relayed off-diagonal cross-peaks correlating spins belonging to the same spin system but not J-coupled.

Related Experiments

Selective 1D RELAY
2D COSY and 2D TOCSY experiment


References:

1. G. Wagner, J. Magn. Reson. 55, 151 (1983)
2. A. Bax & G. Drobny, J. Magn. Reson. 61, 306 (1985)


2D RELAY Experiments

Magnitude-mode one-step 2D RELAY (*cosyqfr1*)
Magnitude-mode one-step 2D RELAY with incremented mixing times (*cosyimqfr1*)
Magnitude-mode two-step 2D RELAY (*cosyqfr2*)
Magnitude-mode two-step 2D RELAY with incremented mixing times (*cosyimqfr2*)
Magnitude-mode three-step 2D RELAY (*cosyqfr3*)

Also See:
2D COSY and 2D COSY-DQF Experiments
2D TOCSY Experiment





cosyqfrl




cosyimqfrl




The delays d2, d3, and d4 define the one-step, two-step, and three-step RELAY periods. They are optimized to $1/4J(\text{HH})$







cosyqfr2



cosyimqfr2



cosyqfr3



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D TOCSY EXPERIMENTS

Experiment Description

A TOCSY (TOtal Crrelation Spectroscopy) experiment is a high-sensitive tool that allows the analysis of individual J-coupled spin systems.

Sample Requirements

TOCSY experiments can be recorded on any type of sample. Solvent-suppressed versions are required for samples dissolved in H₂O.

Hardware Requirements

TOCSY experiments can be recorded on any probehead. Optionally, a pulsed-field gradient coil (highly recommended) is required for gradient-based versions.

NMR Spectrum

The experiment yields a 2D proton-proton correlation map with two different types of signals: i) autocorrelation diagonal peaks, and ii) Off-diagonal cross-peaks correlating spins belonging to the same spin system. TOCSY spectra are usually represented in phase-sensitive mode.

Related Experiments

Selective 1D TOCSY
2D COSY experiment
2D Relayed experiments
2D HMQC-TOCSY and HSQC-TOCSY experiments

References:

1. A. Bax & D.G. Davis, *J. Magn. Reson.* 65, 355-360 (1985)
2. J. Cavanagh & M. Rance, *J. Magn. Reson.* 88, 72-85 (1990)

The delay d9 defines the mixing period in all TOCSY experiments, independently if DIPSI-2 or MLEV are used. Typical values range from 0 to 120 ms and transfer efficiency depends on spin system topologies.

2D TOCSY Experiments

- **Phase-cycled**

Phase-sensitive 2D TOCSY using MLEV (**mlevph** | **mlevphsw**)
Phase-sensitive 2D TOCSY using MLEV with purge pulses before d1 (**mlevphpp**)
Phase-sensitive 2D TOCSY using DIPSI-2 (**dipsi2ph**)

- **Phase-cycled and solvent suppression**

Phase-sensitive 2D TOCSY with presaturation using MLEV (**mlevphpr** | **mlevphpr**)
Phase-sensitive 2D TOCSY with presaturation using MLEV only using first trim pulse (**mlevphpr.2** | **h2osupmlev**)
Phase-sensitive 2D TOCSY with presaturation using MLEV and spoil gradient (**mlevgpphrzs**)
Phase-sensitive 2D TOCSY with presaturation using DIPSI-2 (**dipsi2phpr**)
Phase-sensitive 2D TOCSY with presaturation using DIPSI-2 and spoil gradient (**dipsi2gpphr**)
Phase-sensitive 2D Clean-TOCSY with presaturation using MLEV (**clmlevphpr**)

- **Gradient-based**

Phase-sensitive ge-2D TOCSY with MLEV using echo-antiecho (**mlevetgp**)
Phase-sensitive ge-2D TOCSY with DIPSI-2 using echo-antiecho (**dipsi2etgp**)
Phase-sensitive ge-2D TOCSY with DIPSI-2 using PEP (**dipsi2etgps**)
Phase-sensitive ge-2D TOCSY with DIPSI-2 and Zero-Quantum suppression (**dipsi2gpphzs**)



- **Gradient-based and solvent suppression**

Phase-sensitive 2D TOCSY with WATERGATE (3-9-19) using MLEV (**mlevgpph19** | **mlevgpph19sw**)
Phase-sensitive 2D TOCSY with WATERGATE (3-9-19) using DIPSI-2 (**dipsi2gpph19**)
Phase-sensitive sensitivity-improved 2D TOCSY with WATERGATE (3-9-19) and using DIPSI-2 (**dipsi2etgps19**)
Phase-sensitive 2D Adiabatic TOCSY with WATERGATE (3-9-19) using X_M16 sequence (**atocsygpph19**)

Phase-sensitive 2D TOCSY with excitation sculpting (W5) using MLEV (**mlevgpph5**)
Phase-sensitive 2D TOCSY with excitation sculpting (180 water-selective pulse-ES element) using MLEV (**mlevesgpph**)
Phase-sensitive 2D TOCSY with excitation sculpting (180 water-selective pulse-ES element) using DIPSI-2 (**dipsi2esgpph**)
Phase-sensitive 2D TOCSY with excitation sculpting (180 water-selective pulse-ES element) using DIPSI-2 and optional ¹³C,¹⁵N-decoupling (**dipsi2esfbgpph**)

- **Single/Multiple Solvent Suppression (LC-NMR):**


2D TOCSY with double presaturation and cw-decoupling on f2 (**lcmlevcwfdpcph**)
2D TOCSY with presaturation and cw-decoupling on f2 (**lcmlevcwpcphps**)
2D TOCSY with double presaturation using composite pulse (**lcmlevf2pcph**)
2D TOCSY with double presaturation (**lcmlevf2phpr** | **LCML12**)
2D TOCSY with presaturation using shape pulse and composite pulse (**lcmlevpcphps**)
2D TOCSY with presaturation using composite pulse (**lcmlevpcph**)
2D TOCSY using WET (**mlevdcphwt** | **mlevdcphwt**)




2D TOCSY spectrum of sucrose (mixing time d9 of 90ms). Note that the glucose (red) and fructose (blue) sub-residues are clearly identified.

Phase-Cycled TOCSY


mlevphpp




mlevphpr




mlevphpr.2




clmlevphpr



dipsi2ph



dipsi2phpr



NMR BUILDING BLOCK: MLEV-17

```

"p5=p6*.667"
"p7=p6*2"

"SCALEF=p7*2/p5"
"FACTOR1=((d9-p17*2)/(p6*64+p5))/SCALEF"
"l1=FACTOR1*SCALEF"
.....

4u p110:f1
  (p17 ph26)
;begin MLEV17
4 (p6 ph22 p7 ph23 p6 ph22)
  (p6 ph24 p7 ph25 p6 ph24)
  (p6 ph24 p7 ph25 p6 ph24)
  (p6 ph22 p7 ph23 p6 ph22)
  (p6 ph24 p7 ph25 p6 ph24)
  (p6 ph24 p7 ph25 p6 ph24)
  (p6 ph22 p7 ph23 p6 ph22)
  (p6 ph22 p7 ph23 p6 ph22)
  (p6 ph24 p7 ph25 p6 ph24)
  (p6 ph22 p7 ph23 p6 ph22)
  (p6 ph22 p7 ph23 p6 ph22)
  (p6 ph24 p7 ph25 p6 ph24)
  (p6 ph22 p7 ph23 p6 ph22)
  (p6 ph22 p7 ph23 p6 ph22)
  (p6 ph24 p7 ph25 p6 ph24)
  (p6 ph24 p7 ph25 p6 ph24)
  (p5 ph23)
  lo to 4 times l1
;end MLEV17
  (p17 ph26)

.....
ph22=3
ph23=0
ph24=1
ph25=2
ph26=0

;p110: f1 channel - power level for TOCSY-spinlock
;p5 : f1 channel - 60 degree low power pulse
;p6 : f1 channel - 90 degree low power pulse
;p7 : f1 channel - 180 degree low power pulse
;p17: f1 channel - trim pulse [2.5 msec]
;d1 : relaxation delay; 1-5 * T1
;d9 : TOCSY mixing time
;l1: loop for MLEV cycle: ((p6*64) + p5) * l1) + (p17*2) = mixing time
  
```

The Basic Element of the MLEV (acronym of **Malcom Levitt's** CPD sequence) is the composite R=90(x)-180(y)-90(x) cluster.

The 16 elements in the MLEV-16 supercycle are cycled as RRRR RRRR RRRR RRRR where **R** is the inverse of R.

Reference:

MH Levitt, R. Freeman, T. Frenkiel, J. Magn. Reson. , 47, 328-330 (1982).

NMR BUILDING BLOCK: z-filtered DIPSI-2

```

"FACTOR1=(d9/(p6*115.112))/2"
"l1=FACTOR1*2"

p1 ph2
  d20 p110:f1

;begin DIPSI2
4 p6*3.556 ph23
  p6*4.556 ph25
  p6*3.222 ph23
  p6*3.167 ph25
  p6*0.333 ph23
  p6*2.722 ph25
  p6*4.167 ph23
  p6*2.944 ph25
  p6*4.111 ph23

  p6*3.556 ph25
  p6*4.556 ph23
  p6*3.222 ph25
  p6*3.167 ph23
  p6*0.333 ph25
  p6*2.722 ph23
  p6*4.167 ph25
  p6*2.944 ph23
  p6*4.111 ph25

  p6*3.556 ph25
  p6*4.556 ph23
  p6*3.222 ph25
  p6*3.167 ph23
  p6*0.333 ph25
  p6*2.722 ph23
  p6*4.167 ph25
  p6*2.944 ph23
  p6*4.111 ph25

  p6*3.556 ph23      ;p6 : f1 channel - 90 degree low power pulse
  p6*4.556 ph25      ;d0 : incremented delay (2D)
  p6*3.222 ph23      ;d1 : relaxation delay; 1-5 * T1
  p6*3.167 ph25      ;d9 : TOCSY mixing time
  p6*0.333 ph23      ;d12: delay for power switching           [20 usec]
  p6*2.722 ph25      ;d20: first z-filter delay             [2 msec]
  p6*4.167 ph23      ;d21: second z-filter delay            [3 msec]
  p6*2.944 ph25      ;l1: loop for DIPSI cycle: ((p6*115.112) * l1) = mixing time
  p6*4.111 ph23
  lo to 4 times l1
;end DIPSI2

  d21 p11:f1
  p1 ph3


  ph23=3
  ph25=1
  
```

DIPSI = Decoupling In the Presence of Scalar Interactions



R Element= 320(x), 410(-x), 290(x), 285(-x), 30(x), 245(-x), 375(x), 265(-x), 370(x)

Ref: AJ Shaka, CJ Lee, A. Pines, J. Magn. Reson. 77, 274-293 (1988).

Gradient-based TOCSY using echo/antiecho




TOCSY using WATERGATE




Adiabatic TOCSY:

1. E. Kupce, P. Schmidt, M. Rance & G. Wagner, *J. Magn. Reson.* 135, 361-367 (1998)
2. W. Peti, C. Griesinger & W. Bermel, *J. Biomol. NMR* 18, 199 - 205 (2000)




TOCSY using WATERGATE and echo-antiecho




TOCSY using excitation sculpting solvent suppression


mlevgpwh5




mlevesgpwh




dipsi2esgpwh




dipsi2esfbgpph




dipsi2gpphpr



dipsi2gpphzs




mlevgpphprzs




TOCSY using singlemultiple solvent suppression (LC-NMR applications)


lcmlevpcph




lcmlevpcphps




lcmlevf2phpr




lcmlevf2pcph



lcmlevcwpcphps




lcmlevcwfdpoph




TOCSY using WET solvent suppression

mlevdcphwt



20 AMINOACIDS: 20 DIFFERENT TOCSY PATTERNS



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D NOESY EXPERIMENTS

Experiment Description

A NOESY (Nuclear Overhauser Spectroscopy) experiment provides through-space proton-proton connectivities.

EXSY (Exchange Spectroscopy) experiments use the same pulse sequences.

Sample Requirements

NOESY experiments can be recorded on any type of sample. Solvent-suppressed versions are required for samples dissolved in H₂O.

Hardware Requirements


NOESY experiments can be recorded on any probehead. A pulsed-field gradient coil (highly recommended) is required for gradient-based versions.

NMR Spectrum

The NOESY experiment yields a 2D proton-proton correlation map with two different types of signals: i) autocorrelation diagonal peaks, and ii) Off-diagonal cross-peaks correlating spins close each other. NOESY spectra are represented in phase-sensitive mode in order to distinguish exchange cross-peaks and unwanted anti-phase COSY contributions.

Related Experiments

Selective 1D NOESY
2D ROESY experiment
2D HMQC-NOESY and HSQC-NOESY experiments



2D NOESY spectrum of strychnine. In small molecules, NOE cross-peaks have opposite phase with respect to diagonal autocorrelation peaks.

2D NOESY Experiments

- **1D Version:**

1D NOESY (**noesy1d**)

1D NOESY with presaturation during d1 and mixing time (**noesypr1d**)

1D NOESY with presaturation during d1 and mixing time and spoil gradients (**noesygppr1d**)

1D NOESY with presaturation during d1 and spoil gradients (**noesygppr1d.2**)

1D NOESY with X-decoupling during acquisition (**noesyig1d**)

Also see in "LC-NMR Experiments":

1D NOESY experiments with single/multiple presaturation

- **Phase-cycled:**

Phase-sensitive 2D NOESY (**noesyph** | **NOESYPHSW**)

Phase-sensitive 2D NOESY using purgue pulses before d1 (**noesyphpp**)

Phase-sensitive 2D NOESY using random mixing time (**noesyphrv**)

- **Phase-cycled and solvent suppression:**

Phase-sensitive 2D NOESY with presaturation (**noesyphpr** | **H2OSUPNOESY**)

Phase-sensitive 2D NOESY with presaturation using random mixing time (**noesyphprrv**)

Phase-sensitive 2D NOESY with 1-1 solvent suppression (**noesyph11**)

Phase-sensitive ge-2D NOESY with presaturation using optional ^{13}C , ^{15}N -decoupling (**noesyfbphpr**)

Phase-sensitive 2D NOESY using jump-and-return and optional ^{13}C , ^{15}N decoupling (**noesygpphjrrs**)

- **Gradient-based:**

Phase-sensitive ge-2D NOESY (**noesygpph**)

Phase-sensitive ge-2D NOESY using purgue pulses before d1 (**noesygpphpp**)

Phase-sensitive ge-2D NOESY with z-spoil (**noesygpphzs**)

Phase-sensitive ge-2D NOESY using echo-antiecho (**noesyetgp**)

Phase-sensitive ge-2D NOESY using presaturation (**noesygpphpr**)

- **Gradient-based and solvent suppression:**

Phase-sensitive 2D NOESY with WATERGATE:

Using 3-9-19 (**noesygpph19** | **NOESYGPPH19SW**)

Using water flip-back and 3-9-19 (**noesyfpgpph19**)

Using water flip-back and water-selective 90 pulses (**noesyfpgpphgw**)

Using water flip-back, 3-9-19 and PFG in t_1 (**noesyfpgpphrs19**)

Using water flip-back, water-selective 90 pulses and PFG in t_1 (**noesyfpgpphrswg**)

Phase-sensitive 2D NOESY with excitation sculpting:

Using W5 (**noesygp^hw5**)

Using 180 water-selective pulse (ES element) (**noesyeg^hpph**)

Using 180 water-selective pulse (ES element) and ZQ suppression (**noesyeg^hpph^zs**)

Using 180 water-selective pulse (ES element) and optional ^{13}C , ^{15}N -decoupling

(**noesyef^hbg^hpph**)




Using 180 water-selective pulse (ES element), water flip-back and optional ^{13}C , ^{15}N -decoupling (**noesyef^hg^hpph^rs**)

- **Related Experiment:**

Phase-sensitive 2D NOESY with RELAY and DQF (NOESY-RELAY experiment)
(**noesydf^hphrl**)

Also see:

Selective 1D NOESY, 2D ROESY, 2D HSQC-NOESY & 2D HMQC-NOESY, 3D NOESY-HSQC & 3D HSQC-NOESY-HSQC, 2D & 3D X-filtered NOESY.



The delay d8 defines the mixing period in all NOESY experiments,


NMR Element: NOESY Block

```

...
p1 ph2
d8
p1 ph3
....

ph2=0
ph3=0

;d8 : mixing time
  
```




```

"TAU=d8*0.5-p16-d16-50u"


...
p1 ph2
TAU
50u UNBLKGRAD
p16:gp1
d16
3u
(p2 ph4) :f1
3u
p16:gp1*-1
d16
50u BLKGRAD
TAU
p1 ph3
...

ph2=0
ph3=0
  
```




1D NOESY with presaturation (also see related pulse programs in lc-nmr experiments)


noesypr1d



noesyppr1d




noesyppr1d.2




Classical 2D NOESY with presaturation


noesyph
 noesyphrv




noesyphpp




noesyphpr
 noesyphprv



noesyph11




noesydfphrl




2D NOESY using Gradients


noesygpph




noesygpphpp




noesygpphpr



noesyetgp




noesygpphzs




2D NOESY using WATERGATE


noesygpph19




noesyfpgpph19




noesyfpgpphrs19



noesyfpgpphwg




noesyfpgpphrswg




Classical 2D NOESY using excitation sculpting


noesygp-phw5



noesyegpph




noesyegpphzs




10:40:31:11

2D NOESY with Heteronuclear Decoupling (for labelled compounds)


noesyfbphpr




noesygpphjrs



noesyefbgpph



noesyefpgpphrs




NOE in small molecules is negative

NOE in large molecules is positive

There is a MW range near to NOE null point: **Need for ROESY Experiment!!!**

In NOESY possibility of exchange cross-peaks (positive) that can originate false NOEs (negative), three-spin effects (positive) and COSY cross-peaks (antiphase-J)

In ROESY experiments, the ROE cross-peaks are always negative: There is also possibility of exchange cross-peaks (positive) that can originates false ROEs (negative), three-spin effects (positive) and TOCSY cross-peaks (inphase-J)



NOE transfer strongly depends on Temperature and Solvent (viscosity),
Magnetic field, molecular weight and mixing time

BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D ROESY EXPERIMENTS

Experiment Description

A ROESY (Rotating-Frame Overhauser Spectroscopy) experiment provides through-space proton-proton connectivities in an alternative way to NOESY experiment in null-NOE experimental conditions.

Sample Requirements

ROESY experiments can be recorded on any type of sample.
Solvent-suppressed versions are required for samples dissolved in H₂O.

Hardware Requirements

ROESY experiments can be recorded on any probehead.
A pulsed-field gradient coil (highly recommended) is required for gradient-based versions.


NMR Spectrum

The ROESY experiment yields a 2D proton-proton correlation map with two different types of signals: i) autocorrelation diagonal peaks, and ii) Off-diagonal cross-peaks correlating spins close each other. ROESY spectra are represented in phase-sensitive mode in order to distinguish exchange cross-peaks and unwanted TOCSY contributions.

Related Experiments

Selective 1D ROESY
2D NOESY experiment
2D HMQC-ROESY and HSQC-ROESY experiments

The pulse p15 (in microseconds), applied at power level pl11, defines the mixing period in all ROESY experiments.



2D ROESY Experiments

- **Phase-cycled:**

Phase-sensitive 2D ROESY (**roesyph** | ROESYPHSW)
Phase-sensitive 2D ROESY using purgic pulses before d1 (**roesyphpp**)
Phase-sensitive 2D T-ROESY (**roesyph.2**)
Phase-sensitive 2D T-ROESY using purgic pulses before d1 (**roesyphpp.2**)
Phase-sensitive 2D ROESY with compensation (**croesyph**)
Phase-sensitive off-resonance 2D ROESY (**troesyph**)

- **Phase-cycled and solvent suppression:**

Phase-sensitive 2D ROESY with presaturation (**roesyphpr** | ROESYPHPR)
Phase-sensitive 2D T-ROESY with presaturation (**roesyphpr.2**)
Phase-sensitive 2D ROESY with compensation and presaturation (**croesyphpr**)
Phase-sensitive off-resonance 2D ROESY with presaturation (**troesyphpr**)

- **Gradient-based:**

Phase-sensitive ge-2D ROESY using echo-antiecho (**roesyetgp**)
Phase-sensitive ge-2D ROESY with T-ROESY using echo-antiecho (**roesyetgp.2**)
EASY-ROESY: Phase-sensitive ge-2D ROESY using adiabatic spin-locks and presaturation (**roesyadjsphpr**)

- **Gradient-based and solvent suppression:**

Phase-sensitive 2D ROESY with WATERGATE using 3-9-19 (**roesygp19**)
Phase-sensitive 2D T-ROESY with WATERGATE using 3-9-19 (**roesygp19.2**)
Phase-sensitive 2D ROESY with excitation sculpting using 180 water-selective pulse (ES element) (**roesyegpph**)

Also see:

Selective 1D ROESY Experiments

2D NOESY Experiments



2D HSQC-ROESY Experiments

References:

1. A. Bax & D.G. Davis, *J. Magn. Reson* 63, 207-213 (1985)
2. T.-L. Hwang & A.J. Shaka, *J. Am. Chem. Soc.* 114, 3157-3159 (1992)
3. J. Schleucher, J. Quant, S. Glaser & C. Griesinger, *J. Magn. Reson A* 112, 144-151 (1995)

NMR Element: ROESY Element

```
...  
 4u p111:f1  
  p15 ph2  
....  
  
ph2=0  
  
;p111: f1 channel - power level for ROESY-spinlock  
;p15: f1 channel - pulse for ROESY spinlock
```



NMR Element: T-ROESY Element





```

"l4=p15/(p25*2)"





....
4u p127:f1
4 p25 ph2
  p25 ph3
  lo to 4 times l4
....

ph2=0 2 0 2 1 3 1 3
ph3=2 0 2 0 3 1 3 1
ph31=0 2 2 0 1 3 3 1


;p127: f1 channel - power level for pulsed ROESY-spinlock
;p15: f1 channel - pulse for ROESY spinlock
;p25: f1 channel - 180 degree pulse at p127
;l4: loop for spinlock = p15 / p25*2
  
```




2D ROESY with presaturation



2D ROESY with echo/antiecho




2D ROESY using WATERGATE or excitation sculpting



2D ROESY using adiabatic ramps

C.M. Thiele, K. Petzold & J. Schleucher, *Chem. Eur. J.* 15, 585-588 (2009)



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

1D & 2D MULTIPLE-QUANTUM EXPERIMENTS

Experiment Description

Multiple-Quantum experiments permit to study J-coupled homonuclear spin systems by the generation of MQ coherences. The creation of particular double-quantum (DQ), triple-quantum (TQ)... coherences can be achieved by specific phase cycling or appropriate gradient ratios.

Sample Requirements

MQ experiments can be recorded on any type of sample.
Solvent-suppressed versions are required for samples dissolved in H₂O.

Hardware Requirements

MQ experiments can be recorded on any probehead.
A pulsed-field gradient coil is required for gradient-based versions.

NMR Spectrum

The MQ experiment yields a 2D correlation map correlating MQ frequencies with conventional chemical shifts. For instance, in a DQ experiment, non-coupled resonances (singlets) are filtered out and they do not appear in the spectrum.

Related Experiments

1D & 2D INADEQUATE

References:

1. U. Piantini et al., J. Am. Chem. Soc. 104, 6800 (1982)
 2. T.H. Mareci & R. Freeman, J. Magn. Reson. 51, 531 (1983)
- 2D DQ:
3. C. Dalvit and J-M. Boehlen, J. Magn. Reson. B111, 76 (1996)
 4. C. Dalvit and J-M. Boehlen, J. Magn. Reson. B113, 195 (1996)

The delay d_4 defines the initial echo period that allows to generate the corresponding anti-phase state that is suitable to afford MQ coherences after a second 90 pulse.

1D & 2D DQ Experiments

- **1D Double-Quantum (DQ)**

1D Double-Quantum experiments (*dqs1d*)
1D Multiple Quantum Filter (*mqsqpld* | *mqsqpld2*)

- **2D Phase-cycled Double-Quantum (DQ)**

Magnitude-mode 2D Double-Quantum (DQ) (*dqsqf*)
Phase-sensitive Double-Quantum (DQ) (*dqsph*)
Phase-sensitive 2D Double-Quantum (DQ) with presaturation (*dqsphpr*)

- **2D Gradient-based Double-Quantum (DQ)**

Phase-sensitive ge-2D Double-quantum using echo-antiecho, 45/135 degree conversion pulse for better sensitivity and remote peak minimisation (*dqseagp135*)
Phase-sensitive ge-2D Double-quantum using echo-antiecho (*dqseagp90*)

Also see:
1D & 2D INADEQUATE

NMR Building Block:

Generation of Homonuclear Multiple-Quantum Coherences

```

...
d1
  p1 ph1
  d4
  p2 ph2
  d4
  p1 ph1
  d13
  p0 ph3
  go=2 ph31
  30m mc #0 to 2 F0(zd)
exit

```

```


ph1=0 2
ph2=0 2 0 2 0 2 0 2 0 2 0 2 0 2 0
ph3=0 0 2 2 3 3 1 1
ph31=0 0 2 2 1 1 3 3

```

```

...
;d4 : 1/(4J)


```





NMR Building Block: Gradient-Based MQ Filter

```


  ...
  50u UNBLKGRAD
  p16:gp1
  d16
  p1 ph2
  d13
  p16:gp2
  d16
  p1 ph3
  50u
  p16:gp3
  d16
  4u BLKGRAD
  .....
```



mqsgp1d2




mqsgp1d




10:20 for double quantum filter
 10:30 for triple quantum filter


dqsph
dqsqf




dqsphpr



dqseagp90



dqseagp135



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D J-RESOLVED EXPERIMENTS

Experiment Description

A 2D J-Resolved experiment is a simple variable spin-echo period that allows to decode all J-coupling information from chemical shift information

Sample Requirements

J-Resolved experiments can be recorded on any type of sample. Solvent-suppressed versions are required for samples dissolved in H₂O. In dilute samples, X-decoupling is needed to efficiently remove satellites from the strong solvent resonances

Hardware Requirements


J-Resolved experiments can be recorded on any probehead.

NMR Spectrum

The J-Resolved experiment yields a 2D correlation map correlating J-coupling with conventional chemical shifts.

Related Experiments

2D Heteronuclear J-Resolved



2D J-Resolved Experiments

- **Classical:**


Magnitude-mode 2D J-Resolved (**jresqf**)
Magnitude-mode 2D J-Resolved with f2 decoupling (**jresdcqf**)

- **Gradients and presaturation:**


ge-2D J-Resolved with presaturation and gradients (**jresgpprqf**)

- **With single/multiple solvent suppression (LC-NMR):**


2D J-Resolved with presaturation (**lcjresprqf**)
2D J-Resolved with presaturation using shape pulse (**lcjrespsqf**)
2D J-resolved with double presaturation and cw-decoupling on f2 (**lcjrescwfprqf**)
2D J-resolved with presaturation and cw-decoupling on f2 (**lcjrescwprqf**)
2D J-resolved with presaturation using shape pulse and cw-decoupling on f2 (**lcjrescwpsqf**)
2D J-resolved with double presaturation (**lcjresf2prqf**)




lcjrescwprqf




lcjrescwpsqf



lcjresf2prqf



lcjrescwfdprqf



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

1D INEPT EXPERIMENTS

Experiment Description

The INEPT experiment is based on a heteronuclear pulse sequence in which signal intensities are strongly enhanced by polarization transfer.

Sample Requirements

INEPT experiments can be recorded on any type of sample.

Hardware Requirements

INEPT experiments can be recorded on any probehead.

NMR Spectrum

The non-refocused INEPT version affords anti-phase multiplets with enhanced sensitivity. On the other hand, refocused INEPT versions afford in-phase multiplets that can be decoupled.

Related Experiments

1D DEPT experiments and Other 1D heteronuclear editing methods
2D HETCOR and 2D HSQC correlation experiments


1D INEPT Experiments

INEPT without refocusing (**ineptnd**)
Refocused INEPT with decoupling (**ineptrd**)
Refocused INEPT with decoupling using adiabatic pulses (**ineptrdsp**)
INEPT+ without decoupling (**ineptrnd**)


Non-refocused ¹H-coupled ¹⁵N spectrum using INEPT (**ineptnd**)
¹H-decoupled ¹⁵N spectrum using INEPT (**ineptrd** | **N15INEPT**)
Refocused ¹H-coupled ¹⁵N spectrum using INEPT+ (**ineptrnd**)

1D X-relayed H,X-COSY (**ineptrl1** / **ineptrl2**)

Also see:
DEPT Experiments




Adiabatic Pulses in Refocused INEPT



NMR BUILDING BLOCK: POLARIZATION TRANSFER VIA INEPT

```
"d4=1s/(cnst2*4)"
....
d1
(p3 ph1):f2
d4
(center (p4 ph1):f2 (p2 ph1) )
d4
(p3 ph3):f2 (p1 ph1)
.....
ph1=0
ph2=1
```




The delay d4 is optimized to $1/4J(XH)$.

NOE $1 + \gamma_I/2 * \gamma_S$
 INEPT γ_I/γ_S

X	¹¹ B	¹³ C	¹⁵ N	²⁹ Si	⁵⁷ Fe	¹⁰³ Rh	¹⁰⁹ Ag	¹¹⁹ Sn	¹⁸³ W
NOE	2.56	2.99	-3.94	-1.52	16.48	-16.89	-9.75	-0.41	13.02
INEPT	3.12	3.98	9.87	5.03	30.95	31.77	21.50	2.81	24.04

REFOCUSED INEPT



Delay d3 :
 $1/(6J(XH))$ XH, XH2, XH3 positive
 $1/(4J(XH))$ XH only
 $1/(3J(XH))$ XH, XH3 positive, XH2 negative

$$-2I_z S_y \xrightarrow{d3-180^\circ I_x-180^\circ S_x-d3} S_x$$

Antiphase S Magnetization


In-Phase S Magnetization


CH spin systems $-2I_{1z}S_y \cos + S_x \sin$

CH₂ spin systems $-2I_{1z}S_y \cos \cos + S_x \sin \cos$
 $-2I_{1z}I_{2z}S_y \cos \sin + 2I_{2z}S_y \sin \sin$


CH₃ spin systems $-2I_{1z}S_y \cos \cos \cos + S_x \sin \cos \cos$
 $-4I_{1z}I_{2z}S_y \cos \sin \cos + 2I_{2z}S_y \sin \sin \cos$
 $+4I_{1z}I_{3z}S_x \cos \cos \sin + 2I_{3z}S_y \sin \cos \sin$
 $+8I_{1z}I_{2z}I_{3z}S_x \cos \sin \sin - 4I_{2z}I_{3z}S_x \sin \sin \sin$

Signal intensity of
a $I_n S$ peak depends of
 $\sin(\pi J_{IS}\tau) \cos^{n-1}(\pi J_{IS}\tau)$






ineptr1



Also see inepin

ineptr2



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

1D DEPT EXPERIMENTS

Experiment Description

The DEPT (Distortionless Enhancement Polarization Transfer) experiment is based on a heteronuclear pulse sequence in which signal intensities are enhanced by polarization transfer and their phases depend of a proton flip angle.

Sample Requirements

DEPT experiments can be recorded on any type of sample.

Hardware Requirements

DEPT experiments can be recorded on any probehead.


NMR Spectrum

A DEPT experiment affords a ^{13}C spectrum in which the phase of each signal is edited as a function of carbon multiplicity. For instance, a conventional DEPT135 spectrum gives CH and CH₃ with a positive phase and CH₂ with a negative phase.

Related Experiments

1D INEPT experiments and Other 1D heteronuclear editing methods
2D HETCOR and 2D DEPT-HMQC correlation experiments

Key: Evolution of Multiple-Quantum Coherences
Result: Multiplicity Information



IS Spin Systems:

$$\begin{array}{c}
 \xrightarrow{d2} -2I_x S_y \xrightarrow[180^\circ(x)-S]{\phi^\circ(y)-I} 2I_x S_y \cos(\phi) - 2I_z S_y \sin(\phi) \xrightarrow{d2} S_x \sin(\phi) + \dots
 \end{array}$$


I₁I₂S Spin Systems:

$$\begin{array}{c}
 \xrightarrow{d2} -4I_{1x} I_{2z} S_x \xrightarrow[180^\circ(x)-S]{\phi^\circ(y)-I} -4I_{1x} I_{2z} S_x \cos(\phi) \cos(\phi) + 4I_{1z} I_{2z} S_x \sin(\phi) \cos(\phi) \\
 -4I_{1x} I_{2x} S_x \cos(\phi) \sin(\phi) + 4I_{1z} I_{2x} S_x \sin(\phi) \sin(\phi) \\
 \xrightarrow{d2} S_x \sin(\phi) \cos(\phi) + \dots
 \end{array}$$

I₁I₂I₃S Spin Systems:

$$\begin{array}{c}
 \xrightarrow{d2} -8I_{1x} I_{2z} I_{3z} S_y \xrightarrow[180^\circ(x)-S]{\phi^\circ(y)-I} \\
 8I_{1x} I_{2z} I_{3z} S_y \cos(\phi) \cos(\phi) \cos(\phi) - 8I_{1z} I_{2z} I_{3z} S_y \sin(\phi) \cos(\phi) \cos(\phi) \\
 8I_{1x} I_{2x} I_{3z} S_y \cos(\phi) \sin(\phi) \cos(\phi) - 8I_{1z} I_{2x} I_{3z} S_y \sin(\phi) \sin(\phi) \cos(\phi) \\
 8I_{1x} I_{2z} I_{3y} S_y \cos(\phi) \cos(\phi) \sin(\phi) - 8I_{1z} I_{2z} I_{3y} S_y \sin(\phi) \cos(\phi) \sin(\phi) \\
 8I_{1x} I_{2x} I_{3y} S_y \cos(\phi) \sin(\phi) \sin(\phi) - 8I_{1z} I_{2x} I_{3y} S_y \sin(\phi) \sin(\phi) \sin(\phi) \\
 \xrightarrow{d2} S_x \sin(\phi) \cos(\phi) \cos(\phi) + \dots
 \end{array}$$

Signal intensity of a I_nS peak in the DEPT pulse scheme depends of $\sin(\phi)\cos^{n-1}(2\phi)$



The delay d2 is optimized to $1/2J(XH)$.

1D DEPT Experiments

Classical

- DEPT (**dept**)
- DEPT-45 (**dept45** | $c13DEPT45$)
- DEPT-90 (**dept90** | $c13DEPT90$)
- DEPT-135 (**dept135** | $c13DEPT135$)

Using Adiabatic pulses

- DEPT-45 with adiabatic pulses (**deptsp45** / **deptsp**)
- DEPT-90 with adiabatic pulses (**deptsp90**)
- DEPT-135 with adiabatic pulses (**deptsp135**)

Using Composite pulses

- DEPT with composite pulses (**deptcp**)
- DEPT-45 with composite pulses (**deptcp45**)
- DEPT-90 with composite pulses (**deptcp90**)
- DEPT-135 with composite pulses (**deptcp135**)

Without Decoupling


- DEPT without 1H -decoupling (**deptnd**)
- DEPT++ without 1H -decoupling (**deptppnd**)





DEPTQ for quaternary carbons

- Using adiabatic pulses (**deptqsp**)
- Using gradients and adiabatic pulses (**deptqgpsp**)
- Using gradients, adiabatic pulses and NOE (**deptqgpsp.2**)


Also see:







INEPT Experiments





Adiabatic pulse in DEPT:







DEPTQ

R. Burger & P. Bigler, J. Magn. Reson. 135, 529-534 (1998)



deptqsp



deptqgpsp



deptqgpsp.2



Gradients in DEPTQ

gp1:gp2:gp3
31:31:31 for all C
0:31:31 for CHn only
31:31:11 for Cquat only

BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

OTHER ¹³C EDITING EXPERIMENTS

Other 1D ¹³C-editing Experiments

Spin-Echo or SEFT (**jmod**)
Conventional APT (**apt** | **c13APT**)
APT with J-compensation (**aptjc**)


Quaternary-carbons with decoupling (**quatd**)
Quaternary-carbons without decoupling (**quat**)

Also see:
1D INEPT and 1D DEPT experiments





APT: S.L. Patt & J.N. Shoolery, *J. Magn. Reson.* 46, 535-539 (1982)

ATJ-J Compensated: A.M. Torres, T.T. Nakashima & R.E.D. McClung, *J. Magn. Reson.* A101, 285-294 (1993)

QUAT: M.R. Bendall & D.T. Pegg, *J. Magn. Reson.* 53, 272 (1983)



$d20: 1/(J(XH))$ X, XH2 positive, XH, XH3 negative
 $1/(2J(XH))$ X only



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D X-DETECTED HETCOR EXPERIMENTS

Experiment Description

The HETCOR (Heteronuclear CORrelation) experiment is a 2D version of the 1D INEPT or DEPT pulse schemes. Therefore, it is a X-detected experiment providing heteronuclear correlation between ^1H and X heteronuclei via the scalar coupling constant, $J(\text{XH})$.

Sample Requirements

HETCOR experiments can be recorded on any type of sample but sensitivity suffers of the low receptivity of the detected X nucleus.

Hardware Requirements

X- ^1H HETCOR experiments can be recorded on any probehead that can be tuned to X but a direct probe is recommended.

NMR Spectrum


A 2D HETCOR map correlates ^1H and X chemical shifts via $1J(\text{XH})$.

Related Experiments

1D INEPT and DEPT experiments
2D HMQC and HSQC correlation experiments
2D Inverse-INEPT experiment
2D COLOC experiment

References:

A. Bax & G.A. Morris, *J. Magn. Reson.* 42, 501 (1981)
A. Bax, *J. Magn. Reson.* 53, 517 (1983)
V. Rutar, *J. Magn. Reson.* 58, 306 (1984)
J.A. Wilde & P.H. Bolton, *J. Magn. Reson.* 59, 343 (1984)



2D HETCOR Experiments


- **INEPT-Based HETCOR**
 - Magnitude-mode 2D HETCOR (**hxcoqf** | **hccosw**)
 - Magnitude-mode 2D HETCOR with ²H-decoupling (**hxcoqf2h**)
 - Magnitude-mode 2D HETCOR using composite pulses (**hxcocpqf**)
 - Magnitude-mode 2D HETCOR with ¹H-¹H decoupling in F1 using BIRD (**hxcobiqf**)
 - Magnitude-mode 2D HETCOR with ¹H-¹H decoupling in F1 using BIRD and composite pulses (**hxcobicpqf**)
 - Magnitude-mode 2D HETCOR with refocusing of chemical shifts (**hxinepqf**)
 - Phase-sensitive 2D HETCOR with refocusing of chemical shifts (**hxinepph**)
- **DEPT-based HETCOR**
 - Magnitude-mode DEPT-based 2D HETCOR (**hxdeptqf**)
 - Phase-sensitive DEPT-based 2D HETCOR (**hxdeptph**)
 - Magnitude-mode DEPT-based 2D HETCOR with ¹H-¹H decoupling in F1 using BIRD (**hxdeptbiqf**)
 - Phase-sensitive DEPT-based 2D HETCOR with ¹H-¹H decoupling in F1 using BIRD (**hxdeptbiph**)
 - Phase-sensitive DEPT-based TOCSY-HETCOR experiment (**hxdeptmlph**)
- **2D H-relayed HETCOR experiment**
 - Magnitude-mode 2D H-relayed HETCOR (**hxxcoqf** / **hxxcoqf.2**)
- **2D X-relayed HETCOR experiment**
 - Magnitude-mode 2D X-relayed HETCOR (**hxxcoqf**)
 - Also see in 2H-specific Experiments:
Magnitude-mode 2D HETCOR with ²H-decoupling (**hxcoqf2h**)
- **2D COLOC experiment**
 - Magnitude-mode 2D COLOC (**colocqf** | **hccolocsw**)

Also see in 19F-specific Experiments:
2D ¹⁹F-¹H HETCOR experiment (**hfcoqfn**)

Also see in 2D Inverse-INEPT experiments:
Phase-sensitive ge-2D Inverse INEPT using echo-antiecho (**xhcoetgp**)


NMR BUILDING BLOCK: 1H EVOLUTION and X-DECOUPLING

All HETCOR-type experiments start from 1H and finish with X-detection. The first step is the 1H chemical evolution during the variable d0 period:




Variable period
 Overall Duration: 2d0

^1H Chemical shift: $d0 + d0 = 2d0$
 X Chemical shift: $d0 - d0 = 0$
 $J(\text{XH})$ Coupling Constant: $d0 - d0 = 0$
 $J(\text{HH})$ and $J(\text{XX})$ Coupling Constant: $d0 + d0 = 2d0$




```


...
(p3 ph1):f2
d0
p2 ph4
d0
d4
(center (p4 ph2):f2 (p2 ph4) )
d4
(p3 ph3):f2 (p1 ph5)
d3
(center (p4 ph2):f2 (p2 ph6) )
d3 p112:f2
go=2 ph31 cpd2:f2
d11 mc #0 to 2 F1QF(caldel(d0, +in0))
...
  
```




hxcoqf



hxcocpqf



hxinepqf
hxinepph




d3:
 $1/(6J(XH))$ XH,XH2,XH3 positive
 $1/(4J(XH))$ XH only
 $1/(3J(XH))$ XH,XH3 positive,XH2 negative


NMR Element: BIRD(x) Element

```


....
(p3 ph1):f2
d2
(center (p2 ph1) (p4 ph1):f2 )
d2
(p3 ph1):f2
d0
....
ph1=0
;d2 : 1/(2J(XH))
    
```



BIRD(x) inverts 1H-13C whereas 1H-12C is not affected:




or




BIRD:

A. Bax and S. Subramanian, *J. Magn. Reson.* 67, 565-569 (1986)

BIRD-HETCOR: To minimize J(HH) evolution during the d0 period






P.H. Bolton, *J. Magn. Reson.* 48, 336 (1982)



A. Bax, *J. Magn. Reson.* 53, 149 (1983)

H. Kessler et al., *J. Am. Chem. Soc.* 105, 6944 (1983)

hxxcoqf





hxxcoqf.2




H. Kessler, W. Bermel & C. Griesinger, *J. Magn. Reson.* 62, 573 (1985)

hxxcoqf







M.R. Bendall & D.T. Pegg, *J. Magn. Reson.* 53, 144 (1983)
 T.T. Nakashima et al., *J. Magn. Reson.* 59, 124 (1984)




p0
 45 degree - all positive
 90 degree - XH only
 135 degree - XH, XH3 positive, XH2 negative





H. Kessler, G. Gemmecker, M. Koeck, R. Osowski & P. Schmieder,
Magn. Reson. Chem. 28, 62 - 67 (1990)





Experiment Description

The COLOC experiment is a X-detected experiment providing long-range heteronuclear correlations between ^1H and X heteronuclei via the scalar coupling constant, $J(\text{XH})$. It can be understood as a long-range optimized constant-time version of the HETCOR experiment.

References

H. Kessler et al., J. Magn. Reson. 57, 331 (1984)



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D HETERONUCLEAR J-RESOLVED EXPERIMENT

Experiment Description

The HETJRES (Heteronuclear J-Resolved) experiment is the heteronuclear version of the J-Resolved experiment. Basically, it consists of a variable ^{13}C spin-echo period in which $J(\text{CH})$ evolves whereas ^{13}C chemical shift is refocused.

Sample Requirements

Heteronuclear J-Resolved experiments can be recorded on any type of sample but its sensitivity is related to ^{13}C detection.

Hardware Requirements

Heteronuclear J-Resolved experiments can be recorded on any probehead but better sensitivity will be achieved in direct X probes.

NMR Spectrum

A basic Heteronuclear J-Resolved experiment yields a 2D map in which ^{13}C chemical shift is displayed in the F2 dimension whereas $J(\text{CH})$ is displayed in the F1 dimension. Thus, methine carbons will display a doublet, a methylene carbon a triplet, and a methyl carbon a typical quartet.

The proton-selective version allows the measurement of small, long-range proton-carbon coupling constants. The only requisite is a well-isolated proton resonance.


Related Experiments

Homonuclear 2D J-Resolved experiment

2D Heteronuclear J-Resolved Experiments

Magnitude-mode 2D Heteronuclear J-Resolved (*hjresqf*)

Also see chapter: Measurement of long-range proton-carbon coupling constant.
Magnitude-mode proton-selective 2D Heteronuclear J-Resolved (*seljresqfsp*)



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D HOESY EXPERIMENTS

Experiment Description

The HOESY (Heteronuclear nOESY) experiment is the heteronuclear version of the conventional NOESY experiment.

Sample Requirements

HOESY experiments can be recorded on any type of sample but its sensitivity is very low.

Hardware Requirements

HOESY experiments can be recorded on any probehead but better sensitivity will be achieved in direct X probes. However, HOESY pulse scheme can be also applied with proton detection in inverse probes.

NMR Spectrum

A HOESY experiment provides a 2D heteronuclear correlation map. Directly-attached CH systems will give a strong peak and the important long-range NOE peaks are usually very weak.

Related Experiments

2D NOESY experiment

References

C. Yu & G. Levy, J. Am. Chem. Soc. 106, 6533 (1984)




The delay d8 defines the mixing time in HOESY experiments

2D HOESY Experiments

Magnitude-mode 2D ^1H - ^{13}C HOESY(hoesyqfrv)
Phase-sensitive 2D ^1H -X HOESY (hoesyph)

Also see in Chapter 19F experiments:
2D ^{19}F - ^1H HOESY experiment (hoesyfhqfrv)

Also see:
2D NOESY Experiment



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

1D & 2D INADEQUATE EXPERIMENTS

Experiment Description

The INADEQUATE experiment is the ^{13}C version of the Double-quantum experiment.

Sample Requirements

The INADEQUATE experiment is a very insensitive experiment because two adjacent ^{13}C nuclei are needed. Therefore, a highly concentrated sample is needed for molecules at natural abundance. Otherwise, ^{13}C -labeled material is highly recommended.

Hardware Requirements

The INADEQUATE experiment can be recorded on any probehead but better sensitivity will be achieved in direct X probes.

NMR Spectrum

As similar to the DQ experiment, the INADEQUATE affords a 2D map correlating ^{13}C chemical shift in F2 and double-quantum frequencies ($^{13}\text{C}_\text{a}+^{13}\text{C}_\text{b}$) in the indirect F1 dimension.


Related experiments

1D and 2D Multiple-quantum experiments

References

D.L. Turner, J. Magn. Reson. 49, 175 (1982)

M. Bourdonneau & B. Ancian, J. Magn. Reson. 132, 316-327 (1998)



The delay d_4 is optimized to $1/4J(\text{CC})$.

2D INADEQUATE Experiments

- 1D INADEQUATE


1D INADEQUATE without refocusing (*inad1d*)
1D INADEQUATE using composite pulses (*inadcp1d*)
1D INADEQUATE with refocusing (*inadrd1d*)
1D INADEQUATE using initial INEPT (*inopin*)

- 2D INADEQUATE


Magnitude-mode 2D INADEQUATE (*inadqf* / *inadqf.2* | INAD)
Phase sensitive 2D INADEQUATE (*inadph*)
Phase sensitive 2D INADEQUATE using adiabatic pulse (*inadphsp*)
Magnitude-mode symmetric 2D INADEQUATE (*inadqfsy*)
Magnitude-mode 2D INADEQUATE using adiabatic pulse and gradients (*inadgpqfsp*)

Also see: 1D & 2D Multiple-Quantum Experiments



inad1d




inadcp1d




inadrd1d






inepin




Also see inept11




inadqf.2
 inadph
 inadqf




inadphsp



inadphsp



inadqfsy




BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

DECOUPLER PULSE CALIBRATION

Experiment Description:

In heteronuclear experiments it is necessary to calibrate the duration of the hard 90/180 pulses applied from the decoupler (the non-observed nucleus) as well as the required power level for effective CPD decoupling. A series of simple pulse sequences are used for these purposes. Variants of these experiments also allow to calibrate duration and power levels for shaped pulses as, for instance, adiabatic pulses.



Example: Determination 90 decoupler pulse needed for 1H decoupling (using WALTZ16) whereas observing 13C:


Basic acquisition parameters:

Sample: ASTM (60% C6D6 / 40% p-Dioxane)

Spin: on, d1 45s, ns 1 and ds 0, o1p 66.5ppm; o2p=3.7ppm, CNST2=142Hz; d2=3.52ms

Process with ef (lb=3.5Hz)

Analysis: Without decoupling we observe an antiphase triplet, with the two outer lines in phase and the center line in opposite phase. Set p3=80us and start with a low power p12=30dB. Increase p12 using **popt** until the triplet becomes zero. These values will be used for CPD 1H decoupling experiments.



Decoupler Pulse Calibration

Calibration of the 90 decoupler pulse (**dec90**, **dec90f3**)

Calibration of the 90 decoupler shaped pulse (**dec90sp**)


Calibration of the 180 decoupler pulse (**dec180**)

Calibration of the shaped 180 decoupler pulse (**dec180sp**)

Calibration of the 180 decoupler pulse using presaturation (**dec180pr**, **dec180f3pr**)

Also See:

Pulse Calibration and Tests



Example: Determination 90 decoupler pulse needed for 15N whereas observing 1H:


Basic acquisition parameters:

Sample: 100mM urea 15N, 100mM Ch3OH 13C in DMSO-d6


Spin: off, d1 24s, ns 1 and ds 0, o1p=4ppm o2p=76ppm sw=8ppm, CNST2=88.5Hz (d2=5.649ms)

Process with ef (lb=0.3Hz)


Analysis: Without decoupling we observe an antiphase doublet. Increase p3 (at power level pl2) until the doublet becomes zero. For three-channel system use pulprog=dec90f3.




dec180




dec180sp



dec180pr



dec180f3pr



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D HMQC EXPERIMENT

Experiment Description

The HMQC (Heteronuclear Multiple-Quantum Correlation) experiment is a proton-detected experiment designed to obtain heteronuclear correlation between ^1H and X heteronuclei via the scalar coupling constant, $1J(\text{XH})$

Sample and Hardware Requirements

HMQC experiments can be recorded on any type of sample. It can be recorded on any probehead but an inverse probe equipped with gradients is strongly recommended.

NMR Spectrum

A 2D HMQC map correlates ^1H and X chemical shifts via $1J(\text{XH})$.


Related Experiments

1D Inverse, 2D HSQC and 2D HMBC experiments

Also see other 2D HMQC-related experiments: HMQC-TOCSY, HMQC-NOESY ...

References

1. A. Bax, R.H. Griffey & B.L. Hawkins, *J. Magn. Reson.* 55, 301 (1983)
2. A. Bax and S. Subramanian, *J. Magn. Reson.* 67, 565-569 (1986)
3. V. Sklenar & A. Bax, *J. Magn. Reson.* 74, 469 (1987)
4. A.G. Palmer III, J. Cavanagh, P.E. Wright & M. Rance, *J. Magn. Reson.* 93, 151-170 (1991)



2D HMQC Experiments

- **Phase-cycled:**

Magnitude-mode 2D HMQC (**hmqcqh** | **hmqc**)
Magnitude-mode 2D HMQC without decoupling (**hmqcndqh**)
Magnitude-mode 2D HMQC using BIRD (**hmqcbiqf** | **hmqcBI**)
Magnitude-mode 2D HMQC using BIRD without decoupling (**hmqcbindqh**)
Phase-sensitive 2D HMQC (**hmqcph** | **hmqcPH**)
Phase-sensitive 2D HMQC without decoupling (**hmqcndph**)
Phase-sensitive 2D HMQC using BIRD (**hmqcbiph** | **hmqcBIPH**)
Phase-sensitive 2D HMQC using BIRD without decoupling (**hmqcbindph**)

- **Phase-cycled and solvent suppression**

From f2 channel:

Phase-sensitive 2D HMQC with presaturation (**hmqcphpr** | **hmqcPHPR**)
Phase-sensitive 2D HMQC using BIRD and presaturation (**hmqcbiphpr**) / **hmqcbiphpr2**)
Phase-sensitive 2D HMQC with 1-1 water suppression (**hmqcph11**)

From f3 channel:

Phase-sensitive 2D ^1H - ^{15}N HMQC (**hmqcf3ph**)
Phase-sensitive 2D ^1H - ^{15}N HMQC using presaturation (**hmqcf3phpr**)
Phase-sensitive 2D ^1H - ^{15}N HMQC using BIRD (**hmqcbif3ph**)
Phase-sensitive 2D ^1H - ^{15}N HMQC using decoupling in a third f2 channel (**hmqcfbph**)

- **Gradient-based:**

From f2 channel:

Magnitude-mode ge-2D HMQC (**hmqcgpqh** | **hmqcGP**)
Phase-sensitive ge-2D HMQC using echo-antiecho (**hmqcetgp**)
Phase-sensitive ge-2D HMQC using echo-antiecho with adiabatic refocusing (**hmqcetgp.2**)
Phase-sensitive ge-2D HMQC using PEP (**hmqcetgps**)
Phase-sensitive ge-2D HMQC using PEP and shorter overall timing (**hmqcetgps.2**)

From f3 channel:

Phase-sensitive ge-2D ^1H - ^{15}N HMQC using echo-antiecho (**hmqcetf3gp**)
Phase-sensitive ge-2D ^1H - ^{15}N HMQC using PEP (**hmqcetf3gps**)
Phase-sensitive ge-2D ^1H - ^{15}N HMQC using PEP and shorter overall timing (**hmqcetf3gps.2**)


- **Gradient-based and solvent suppression**

Phase-sensitive ge-2D ^1H - ^{15}N HMQC using WATERGATE (3-9-19) (**hmqcf3gpqh19**)
Phase-sensitive ge-2D ^1H - ^{15}N HMQC using WATERGATE (3-9-19) in the middle of t1 (**hmqcf3gpqh19.2**)
Phase-sensitive ge-2D ^1H - ^{15}N HMQC using WATERGATE (selective 90 pulses) in the middle of t1 (**hmqcf3gpqh19wg**)

Also see in chapter "SOFAST/BEST experiments":

NH-selective 2D Sofast-HMQC experiment for rapid pulsing (**sfhmqcf3gpqh** | **SFHMQCF3GPPH**)
NH-selective 2D Sofast-HMQC with inversion of water/aliphatic protons (**hetsfhmqcf3gpqh**)
NH-selective 2D Sofast-HMQC with sensitivity improved (**sfhmqcf3gpqhiasi**)

Also see other related correlation experiments: 2D HMQC-DEPT, 2D HSQC, 2D CT-HMQC...




NMR Building Block:
Generation of Heteronuclear Multiple-Quantum Coherences




```

...
2 d1
  d12 p11:f1 p12:f2
3 p1 ph1
  d2 p12:f2
  p3:f2 ph3
....
ph1=0
ph3=0

;d2 : 1/(2J)XH
    
```



Evolution of **Multiple-Quantum** Coherences!!!




NMR Element: BIRD(-x)-Recovery Delay Element

```


.....
p1 ph1
d2
(center (p2 ph1) (p4 ph1):f2 )
d2
p1 ph2
d7
.....

ph1=0
ph2=2

;d2 : 1/(2J)XH
;d7 : delay for inversion recovery
;      optimize to give null for protons bound to C-12
    
```




BIRD(-x) inverts 1H-12C whereas 1H-13C is not affected:




BIRD:


A. Bax and S. Subramanian, J. Magn. Reson. 67, 565-569 (1986)




hmqcbiph
hmqcbiqf




hmqcbindph
hmqcbindqf




hmqcbiphpr




hmqcbiphpr2






hmqcf3ph




p116

hmqcf3phpr




p116

hmqcbif3ph



p116


hmqcfbph



p116


p112

hmqcgpqf




gp 1 : gp 2 : gp 3
 50 : 30 : 40.1 for C-13
 70 : 30 : 50.1 for N-15


hmqcetgp




hmqcetgp.2




hmqcetgpsi




hmqcetgpsi.2




hmqcetf3gp



hmqcetf3gpsi




hmqcetf3gpsi.2




HMQC with WATERGATE


hmqcf3gpph19



hmqcf3gpph19.2



hmqcf3gpphwg



1D HMQC

- Phase-Cycled:


- 1D HMQC with refocusing but not decoupling (`hmqcndrd1d` | `HMQC1D`)
- 1D HMQC without refocusing and without decoupling (`inv3nd1d`/ `hmqcnd1d`)
- 1D HMQC with refocusing and decoupling (`hmqcrd1d`)
- 1D HMQC using BIRD without refocusing and without decoupling (`hmqcbind1d`)
- 1D HMQC using BIRD with refocusing and without decoupling (`hmqcbindrd1d`)
- 1D HMQC using BIRD with refocusing and decoupling (`hmqcbird1d`)




- Gradient-based:

- ge-1D HMQC with refocusing but not decoupling (`hmqcgpnd1d`)

Inverse 1D HMQC

A. Bax, R.H. Griffey & B.L. Hawkins, J. Magn. Reson. 55, 301 (1983)





gp 1 : gp 2 : gp 3
 50 : 30 : 40.1 for C-13
 70 : 30 : 50.1 for N-15

BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D HSQC EXPERIMENT

Experiment Description

The HSQC (Heteronuclear Single-Quantum Correlation) experiment is a proton-detected experiment designed to obtain heteronuclear correlation between ^1H and X heteronuclei via the scalar coupling constant, $1J(\text{XH})$

Sample Requirements

HSQC experiments can be recorded on any type of sample.

Hardware Requirements

HSQC experiments can be recorded on any probehead but an inverse probe equipped with gradients is strongly recommended.

NMR Spectrum

A 2D HSQC map correlates ^1H and X chemical shifts via $1J(\text{XH})$.




Related Experiments

1D Inverse, 2D HMQC and 2D HSQC-edited experiments

Also see other 2D HSQC-related experiments: HSQC-TOCSY, HSQC-NOESY ...

References

G. Bodenhausen & D.J. Ruben, Chem. Phys. Lett. 69, 185 (1980)



2D HSQC Experiments (f2 channel)

- **Phase-cycled:**

Phase-sensitive 2D HSQC (**hsqcph**)

- **Phase-cycled and solvent suppression:**

Phase-sensitive 2D HSQC with presaturation (**hsqcphpr**)

Phase-sensitive 2D HSQC with off-resonance presaturation (**hsqcphps**)

- **Gradient-based:**

ge-1D HSQC with refocusing and no decoupling (**hsqcgpnd1d**)

Phase-sensitive ge-2D HSQC using z-filter and selection before t1 (**hsqcgpph** | **HSQC6P**)

Phase-sensitive ge-2D HSQC using z-filter and selection after t1 (**hsqcgpph2**)

Phase-sensitive ge-2D HSQC using echo-antiecho (**hsqcetgp** | **HSQCET6P**)

Phase-sensitive ge-2D HSQC using echo-antiecho and adiabatic pulses for inversion

(**hsqcetgpsp**)

Phase-sensitive ge-2D HSQC using echo-antiecho and adiabatic pulses for inversion and

refocusing(**hsqcetgpsp.2**)

Phase-sensitive ge-2D HSQC using echo-antiecho and adiabatic pulses for inversion and

refocusing and BS effects(**hsqcetgpsp.3**)

Phase-sensitive ge-2D HSQC using PEP (**hsqcetgpsi**)

Phase-sensitive ge-2D HSQC using PEP with gradients in back-inept (**hsqcetgpsi2**)

Phase-sensitive ge-2D HSQC using PEP with gradients in back-inept and water flip-back

(**hsqcetfpgpsi2**)

Phase-sensitive ge-2D HSQC using PEP and adiabatic pulses for inversion (**hsqcetgpsisp** |

HSQCET6PSISP)

Phase-sensitive ge-2D HSQC using PEP and adiabatic pulses for inversion with gradients in

back-inept (**hsqcetgpsisp2**)

Phase-sensitive ge-2D HSQC using PEP and adiabatic pulses for inversion and refocusing

(**hsqcetgpsisp.2** | **HSQCET6PSISP.2**)

Phase-sensitive ge-2D HSQC using PEP and adiabatic pulses for inversion and refocusing with

gradients in back-inept (**hsqcetgpsisp2.2**)

Phase-sensitive ge-2D HSQC using PEP and adiabatic pulses for inversion and refocusing with

gradients in back-inept and presaturation (**hsqcetgpprsisp2.2**)

Phase-sensitive ge-2D HSQC using PEP and adiabatic pulses for inversion and refocusing with

gradients in back-inept with suppression of HSQC-COSY peaks (**hsqcetgpsisp23**)

ge-2D ¹H-X HSQC experiment with X-Y-decoupling during acquisition and with selective C β /C=O

decoupling. (**hsqcdhetgpsp**)

Phase-sensitive ge-2D HSQC with WET solvent suppression (**hsqcetgpsiwT** | **HSQCET6PSIWt**)

2D HSQC Experiments (f3 channel)

- **Phase-cycled:**

Phase-sensitive 2D ^1H - ^{15}N HSQC (**hsqcf3ph**)

- **Phase-cycled and solvent suppression:**

Phase-sensitive 2D ^1H - ^{15}N HSQC using presaturation (**hsqcf3phpr**)

- **Gradient-based:**

Phase-sensitive ge-2D ^1H - ^{15}N HSQC using echo-antiecho (**hsqcetf3gp** | **HSQCETf3GP**)

Phase-sensitive ge-2D ^1H - ^{15}N HSQC using PEP (**hsqcetf3gpsi** | **HSQCETf3GPSt**)

Phase-sensitive ge-2D ^1H - ^{15}N HSQC using PEP with gradients in back-inept (**hsqcetf3gpsi2**)

Phase-sensitive ge-2D ^1H - ^{15}N HSQC using XY16-CPMG (**hsqcetf3gpxy**, **hsqcetf3gpxy.2**)

- **Gradient-based and solvent suppression**

Phase-sensitive ge-2D ^1H - ^{15}N HSQC using water flip-back and echo-antiecho (**hsqcetf3gp** | **HSQCETf3GP**)

Phase-sensitive ge-2D ^1H - ^{15}N HSQC using water flip-back and PEP (**hsqcetf3gpsi** | **HSQCETf3GPSt**)

Phase-sensitive ge-2D ^1H - ^{15}N HSQC using water flip-back and PEP with gradients in back-inept (**hsqcetf3gpsi2**)

Phase-sensitive ge-2D ^1H - ^{15}N HSQC using WATERGATE (3-9-19) (**hsqcf3gp**)

Fast-HSQC, Phase-sensitive ge-2D ^1H - ^{15}N HSQC using WATERGATE (3-9-19) (**fhqcf3gp** | **FHSQCF3GPPH**)

Phase-sensitive ge-2D ^1H - ^{15}N HSQC using water flip-back and WATERGATE (selective pulse) (**hsqcf3gp** | **HSQCF3GPPHWG**)

- **Gradient-based simultaneous CN-HSQC**


Phase-sensitive ge-2D simultaneous ^1H - ^{15}N / ^1H - ^{13}C HSQC using echo-antiecho (**hsqcetgpsismsp**)

Also see in chapter SOFAST/BEST experiments:

NH-selective SOFAST/BEST-HMQC experiment for rapid pulsing (**b_hsqcetf3gpsi**)

Also see other related correlation experiments (HSQC, TROSY ...), hybrid experiments (HSQC-TOCSY, HSQC-NOESY...) or related experiments for relaxation or coupling constants measurements

HSQC: Sequential Four-Step Experiment



- 1** 1H-to-X INEPT transfer optimized to $1/4J(XH)$
- 2** $\delta(X)$ evolution and $J(CH)$ refocusing during variable d_0 period
- 3** X-to-1H retro-INEPT transfer optimized to $1/4J(XH)$
- 4** 1H detection with optional X-decoupling


Options for HSQC included in this chapter:

- Use of Trim Pulses in step 1
- Use of Gradients for Coherence Selection (step 2 and 3)
- Use of Gradients for purge (in steps 1 and/or 3)
- Water-flip back pulses (in zz-filters and step 3)
- Sensitivity Improved (PEP) version (step 3)
- Use of WATERGATE solvent suppression (in step 3)
- Use of Adiabatic 180° Pulses instead of hard 180 pulses in steps 1,2 and 3.


Options for HSQC included in other chapters:

- Multiplicity-editing incorporated in step 2
- Constant-Time periods incorporated in step 2
- Spin-state selection incorporated in steps 2 and/or 3
- Relaxation blocks in steps 2 or 3
- Additional mixing processes between steps 3 and 4


hsqcph




hsqcf3ph




hsqcphpr




hsqcf3phpr



hsqcphps




hsqcgpnd1d




gp 1 : gp 2 : gp 3
 80 : 30 : 20.1 for C-13
 80 : 30 : 8.1 for N-15

hsqcgpqh




p112

hsqcgpqh2




p112

gp 1 : gp 2 : gp 3
 30:80:20.1 for C-13
 30:80:8.1 for N-15



NMR Building Block: Modified INEPT Transfers in HSQC-Type Experiments

A. INEPT with a trim pulse




```

.....
d12 p11:f1 p12:f2
3 (p1 ph1)
d4
(center (p2 ph1) (p4 ph1):f2 )
d4
p28 ph1
4u
(p1 ph2) (p3 ph1):f2
.....
  
```

➔ An optional short trim pulse of duration p28=1msec applied at high power from the x axis removes any residual Iy contribution.

B. INEPT with zz-purge Gradient:




```

.....
d12 p11:f1 p12:f2
3 (p1 ph1)
d4
(center (p2 ph1) (p4 ph1):f2 )
d4 UNBLKGRAD
(p1 ph2)
3u
p16:gp1
d16
(p3 ph1):f2
.....
  
```

A gradient G1 applied just after the 90°y(1H) removes any residual transverse magnetization. IzSz coherences are not dephased by this purge gradient.

The trim pulse and the purge gradient are independents and they can be applied simultaneously

C. INEPT with Water Flip-Back and zz-purge Gradient:




```


.....
d12 p11:f1 p12:f2
3 (p1 ph1)
d4
(center (p2 ph1) (p4 ph1):f2 )
d4 UNBLKGRAD
(p1 ph2)
4u p10:f1
(p11:sp1 ph1:r):f1
4u p11:f1
p16:gp1
d16
(p3 ph1):f2
.....
    
```

A water-selective 90° pulse creates transverse magnetization for the water whereas other signals are unaffected. The gradient G1 dephases water signal


D. INEPT with Adiabatic 180° pulse:




Hard 180° ^{13}C pulses suffer of unwanted off-resonance effects. In 500MHz spectrometers (and higher) it is recommended to use a 500ms long adiabatic pulse as a inversion ^{13}C element to improve sensitivity.




HSQC with Gradient Selection




gp 1 : gp 2
 80:20.1 for C-13
 80:8.1 for N-15




hsqcetgppsp.2




hsqcetgppsp.3



hsqcdhetgppsp



HSQC: Sensitivity Improved Version



- 1 $J(\text{CH})$ evolution using a 1H-to-X INEPT transfer
- 2 $\delta(\text{C})$ evolution and $J(\text{CH})/J(\text{NH})$ refocusing
Refocused-Echo to accomodate gradient G1 encoding
- 3 $J(\text{CH})$ evolution using a X-to-1H retro-INEPT transfer
- 4 $J(\text{CH})$ evolution using a X-to-1H retro-INEPT transfer (PEP methodology)
- 5 Echo to accomodate gradient G2 decoding
- 6 1H detection with optional X-decoupling




References SI:

1. A.G. Palmer III, J. Cavanagh, P.E. Wright & M. Rance, J. Magn. Reson. 93, 151-170 (1991)
2. L.E. Kay, P. Keifer & T. Saarinen, J. Am. Chem. Soc. 114, 10663-5 (1992)
3. J. Schleucher, M. Schwendinger, M. Sattler, P. Schmidt, O. Schedletzky, S.J. Glaser, O.W. Sorensen & C. Griesinger, J. Biomol. NMR 4, 301-306 (1994)


SI version affords a theoretical sensitivity gain of $S/N=100\%$ in IS spin systems

HSQC with PEP


gp 1 : gp 2
80:20.1 for C-13
80:8.1 for N-15




hsqcetgpsisp




hsqcetgpsisp2




hsqcetgpsisp2.3




hsqcetgpsisp.2




hsqcetgpsisp2.2




hsqcetgpprsisp2.2






HSQC with WET solvent Suppression




HSQC from f3 channel




hsqcetf3gp




hsqcetf3gpsi




hsqcetf3gpsi2




¹H- ¹⁵N HSQC Ubi quit i na 1mM 500 MHz
Optional ¹³C decoupling during t1
in doubly-labeled proteins
zgoptns -DLABEL_CN in eda



hsqcetfpf3gpsi





hsqcetfpf3gpsi2





HSQC using WATERGATE

FAST-HMQC:

S. Mori, C. Abeygunawardana, M. O'Neil-Johnson & P.C.M. van Zijl, J. Magn. Reson. B 108, 94-98 (1995)






CPMG-HSQC:


F.A.A. Mulder, C.A.E.M. Spronk, M. Slijper, R. Kaptein & R. Boelens, *J. Biomol. NMR* 8, 223-228 (1996)

**NMR Building Block:
 A CPMG-INEPT Block**




```
"DELTA1=(d26*2/16-p22)/2"
....
4 DELTA1
  (center (p2 ph6) (p22 ph16):f3 )
  DELTA1 ipp6 ipp16
  lo to 4 times 16
....
ph6=0
ph16=0
```

hsqcetf3gpxy





hsqcetf3gpxy.2



Simultaneous C,N-HSQC:

M. Sattler, M. Maurer, J. Schleucher & C. Griesinger, *J. Biomol. NMR* 5, 97-102 (1995)



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D CONSTANT-TIME HSQC AND HMQC EXPERIMENTS

Experiment Description

The Constant-time version of the HSQC experiment (CT-HSQC) is designed to remove the $J(CC)$ evolution during the variable d_0 period.

Sample Requirements

CT-HSQC experiments are recorded in ^{13}C -labeled compounds. Pulse programs also includes suitable ^{15}N decoupling and therefore they are ready to be applied on ^{15}N -labeled compounds.

Hardware Requirements

CT-HSQC experiments can be recorded on any probehead but triple-resonance probes equipped with gradients is strongly recommended.

NMR Spectrum

A 2D CT-HSQC map correlates 1H and X chemical shifts via $1J(XH)$ and $J(CC)$ is removed from the F_1 dimension.

Related Experiments

2D HSQC experiments

Also see other 2D CT-HSQC to measure coupling constants in labeled proteins.

References

G.W. Vuister & A. Bax, J. Magn. Reson. 98, 428-435 (1992)

2D Constant-time Heteronuclear Correlations

Phase-sensitive Constant-time *ge*-2D HSQC (CT-HSQC)

- Using adiabatic pulses (**hsqcctetgppsp**)
- Using adiabatic pulses without *CO* refocusing (**hsqcctetgppsp.2**)
- Using adiabatic pulses and PEP (**hsqcctetgppisp**)
- 2D H-1/C-13 CT-HSQC for measuring $1J(\text{CH})$ via J-modulation (**hsqcctetgppjc**)


Phase-sensitive Constant-time *ge*-2D HMQC (CT-HMQC)

- Using adiabatic pulses (**hmqcctetgpp**)
- For correlating CH_2 groups (**hmqcctetgpp.2**)


Also see in Coupling Constants in Proteins:

2D spin-echo difference CT-HSQC (**hsqcctetgppjldr** | **HSQCCTETGPPJCLR**) - $3J[\text{CO}-\text{CB}]$ via spin-echo difference

Also see: 2D HMQC and 2D HSQC experiments



Evolution Period in Conventional HSQC




Overall Duration: $2d+2d0$


X Chemical shift: $d0 + d0 + d - d = 2d0$

J(XH) Coupling Constant: $d0 - d0 - d + d = 0$

J(HH) and J(XX) Coupling Constant: $d0 + d0 + d + d = 2d0 + 2d$



Evolution Period in Constant-Time HSQC




Variable Constant-time period
Overall Duration: $2d23$


X Chemical shift: $d0 + d23 - (d23-d0) = 2d0$



J(XH) Coupling Constant: $d0 - d23 + (d23-d0) = 0$

J(HH) and J(XX) Coupling Constant: $d0 + d23 + (d23-d0) = 2d23$



d23: $d23 = T : 13.3$ or 26.6 msec
 $2T$ (constant time period) = $n/J(CC)$






N. Tjandra & A. Bax, *J. Magn. Reson.* 124, 512-515 (1997)

CT-HMQC:


J.P. Marino, J.L. Diener, P.B. Moore & C. Griesinger, J. Am. Chem. Soc. 119, 7361-7366 (1997)

$$d23 = T, 2T \text{ (constant time period)} = n/J(CC) \quad [8.8 \text{ msec}]$$

hmqcctetgp.2



hmqcctetgp



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D INVERSE HETCOR EXPERIMENT

Experiment Description:

The 2D Inverse HETCOR experiment can be understood as the proton-detected version of the HETCOR experiment with the use of gradients for coherence selection. Experimentally, the inverse INEPT experiment is quite similar to an HSQC experiment with the major difference that the pulse scheme starts from carbon magnetization instead of proton magnetization.

Sample and hardware requirements are the same as described for the HSQC experiment.

NMR Experiment Sensitivity

$$S/N \propto n \gamma_e \sqrt{\gamma_d^3 B_o^3 t}$$



	<u>START NUCLEUS</u>	<u>DETECTED NUCLEUS</u>
<u>HETCOR</u>	1H	13C
<u>HMQC/HSQC</u>	1H	1H
<u>INVERSE INEPT</u>	13C	1H



IMPORTANT: The pre-scan delay d1 must optimized as a function of the T1 values for the starting nucleus

Inverse HETCOR Experiments

- 1D Inverse INEPT/DEPT:

1D inverse DEPT with refocusing and no decoupling (**ideptnd**)

1D inverse INEPT without refocusing and without decoupling (**iineptnd**)

1D inverse INEPT with refocusing and decoupling (**iineptrd**)

- 2D Inverse INEPT:

Phase-sensitive *ge*-2D Inverse INEPT using echo-antiecho (**xhcoetgp**)

Also see:

2D Inverse ^1H - ^{31}P INEPT experiments for Nucleic Acids

Using echo-antiecho (**na_xhcoetf3gp**)


Constant-time using echo-antiecho (**na_xhcoctetf3gp**)

2D HETCOR, 2D HMQC and 2D HSQC experiments

Inverse DEPT:

M.R. Bendall, D.T. Pegg, D.M. Doddrell & J. Field, *J. Magn. Reson.* 51, 520 - 526 (1983)

p0:
45 degree - all positive
90 degree - XH only
135 degree - XH, XH3 positive, XH2 negative




Inverse INEPT:


M.R. Bendall, D.T. Pegg, D.M. Doddrell & J. Field, J. Mgn. Reson. 51, 520 - 526 (1983)

d20: $1/(6J(XH))$ for all multiplicities
 d4 : $1/(4J(XH))$


iineptnd



iineptrd



xhcoetgp



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D MULTIPLICITY-EDITED HSQC EXPERIMENT

Experiment Description

The multiplicity edited HSQC experiment is a simple variant of the conventional HSQC pulse scheme designed to obtain simultaneously the heteronuclear correlation between ^1H and X heteronuclei via the scalar coupling constant, $1J(\text{XH})$, and the multiplicity of each cross-peak

Sample Requirements

Multiplicity-edited HSQC experiments can be recorded on any type of sample.

Hardware Requirements

Multiplicity-edited HSQC experiments can be recorded on any probehead but an inverse probe equipped with gradients is strongly recommended.

NMR Spectrum


A 2D Multiplicity-edited HSQC map correlates ^1H and X chemical shifts via $1J(\text{XH})$. CH and CH₃ cross peaks present opposite phase with respect to CH₂ cross peaks.

Related Experiments

1D DEPT and other carbon edited methods
2D HSQC experiments
2D HMQC-DEPT

References:

W. Willker, D. Leibfritz, R. Kerssebaum & W. Bermel, *Magn. Reson. Chem.* 31, 287-292 (1993)



2D Multiplicity-edited HSQC Experiments

- Gradient-enhanced form f2 channel

Phase-sensitive *ge*-2D multiplicity-edited HSQC using z-filter (**hsqcedgpph** | **HSQCEDGPPH**)

Phase-sensitive *ge*-2D multiplicity-edited HSQC using echo-antiecho (**hsqcedetgp** | **HSQCEDETGP**)

Phase-sensitive *ge*-2D multiplicity-edited HSQC using echo-antiecho and adiabatic pulses (**hsqcedetgppsp**)

Phase-sensitive *ge*-2D multiplicity-edited HSQC using echo-antiecho and inversion and matched sweep adiabatic pulses (**hsqcedetgppsp.3**)

Phase-sensitive *ge*-2D multiplicity-edited HSQC using PEP and adiabatic inversion pulses (**hsqcedetgppisp**)

Phase-sensitive *ge*-2D multiplicity-edited HSQC using PEP and adiabatic inversion and refocusing pulses (**hsqcedetgppisp.2**)

Phase-sensitive *ge*-2D multiplicity-edited HSQC using PEP and adiabatic inversion pulses with gradients in back-inept (**hsqcedetgppisp2**)

Phase-sensitive *ge*-2D multiplicity-edited HSQC using PEP and adiabatic inversion and refocusing pulses with gradients in back-inept (**hsqcedetgppisp2.2**)


Phase-sensitive *ge*-2D multiplicity-edited HSQC using PEP and inversion, refocusing and matched sweep adiabatic pulses with gradients in back-inept (**hsqcedetgppisp2.3**)

Phase-sensitive *ge*-2D multiplicity-edited HSQC using PEP and inversion, refocusing and matched sweep adiabatic pulses with gradients in back-inept with HSQC-COSY peak suppression (**hsqcedetgppisp2.4**)


- Gradient-enhanced form f3 channel

Phase-sensitive *ge*-2D ^1H - ^{15}N HSQC-edited using PEP (**hsqcedetf3gpsi**)

Phase-sensitive *ge*-2D ^1H - ^{15}N HSQC-edited using PEP with gradients in back-inept (**hsqcedetf3gpsi2**)



Evolution Period in Conventional HSQC



Overall Duration: $2d+2d_0$


$$\text{X Chemical shift: } d_0 + d_0 + d - d = 2d_0$$

$$\text{J(XH) Coupling Constant: } d_0 - d_0 - d + d = 0$$

$$\text{J(HH) and J(XX) Coupling Constant: } d_0 + d_0 + d + d = 2d_0 + 2d$$



Evolution Period in Multiplicity-edited HSQC




Overall Duration: $2d_{21}+2d_0$

$$\text{X Chemical shift: } (d_{21}+d_0) + d_0 - d_{21} = 2d_0$$



$$\text{J(XH) Coupling Constant: } (d_{21}+d_0) - d_0 + d_{21} = 2d_{21}$$

$$\text{J(HH) and J(XX) Coupling Constant: } (d_{21}+d_0) + d_0 + d_{21} = 2d_{21} + 2d_0$$


d_{21} : set d_{21} according to multiplicity selection
 $1/(2J(\text{XH}))$ XH, XH3 positive, XH2 negative





gp 1 : gp 2
 80:20.1 for C-13
 80:8.1 for N-15




hsqcetdgpsisp




hsqcetdgpsisp2




hsqcedetgpsisp.2




hsqcedetgpsisp2.2





hsqcedetgpsisp2.3




hsqcedetgpsisp2.4





Edited-HSQC in 1D mode:



hsqcedetf3gpsi



hsqcedetf3gpsi2



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D DEPT-HMQC EXPERIMENT

Experiment Description:

The 2D DEPT-HMQC experiment is a variant of the HMQC experiment in which a DEPT pulse train has been incorporated as an initial preparation block to obtain carbon multiplicity information. Also see multiplicity-edited HSQC experiment for similar purposes.

References:

H. Kessler, P. Schmieder & M. Kurz, J. Magn. Reson 85, 400-405 (1989)

HMQC-DEPT Experiments



1D DEPT-HMQC

- 1D DEPT-HMQC with refocusing and decoupling (**indecord1d**)
- 1D DEPT-HMQC using BIRD with refocusing and decoupling (**indecobird1d**)

2D DEPT-HMQC




- Phase-sensitive 2D DEPT-HMQC (**indecoph**)
- Phase-sensitive 2D DEPT-HMQC using BIRD (**indecobiph**)
- Phase-sensitive 2D DEPT-HMQC-TOCSY using BIRD (**indecobimlph**)

Also see:
 2D HMQC Experiments



p0 in 1D DEPT-HMQC:

- 60 degree - XH positive, XH2, XH3 negative
- 90 degree - XH2 only
- 120 degree - XH3 positive, XH, XH2 negative
- 180 degree - XH2 positive, XH, XH3 negative




BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D SPIN-STATE-EDITED HSQC EXPERIMENTS

The HSQC pulse scheme can be modified in a great number of ways in order to afford different versions that can be used for many different purposes.

Of particular interest are the spin-state edited versions that allows to simplify the analysis and interpretation of conventional coupled spectra resulting of recording HSQC experiments without decoupling during acquisition (F2-coupled), during d0 (F1-coupled) or during both periods (F1,F2 coupled).



2D Spin-state Edited HSQC Experiments


From f2 channel

- *ge*-2D ^1H - ^{13}C HSQC-IPAP with adiabatic inversion pulses and IPAP editing in the indirect F1 dimension (**hsqcetgpiasp**)

From f3 channel



- *ge*-2D ^1H - ^{15}N α,β -HSQC with IPAP editing in the F2 detected dimension (**hsqctf3gpss**)
- *ge*-2D ^1H - ^{15}N HSQC-IPAP using watergate with IPAP editing in the indirect F1 dimension (**hsqcf3gpiaphwg**)
- *ge*-2D ^1H - ^{15}N HSQC-IPAP using watergate with similar IP and AP delays and IPAP editing in the indirect F1 dimension (**hsqcf3gpiaphwg.2**)
- *ge*-2D ^1H - ^{15}N HSQC-IPAP using watergate, water flip-back and sensitivity improvement and with IPAP editing in the indirect F1 dimension (**hsqcf3gpiaphsiwg**)

Also see all experiments in "2D TROSY experiments". TROSY is a special case of spin-state selection.





Comparison of ^1H - ^{15}N HSQC spectra of ubiquitin: decoupled (left), F2-coupled (middle) and F1-coupled (right).

Spin-State Selection in Reverse INEPT




Comparison of ^1H - ^{15}N HSQC spectra of ubiquitin: decoupled (left), alpha-spin-edited (middle) and beta-spin-edited (right).




IPAP-HSQC

1. M. Ottiger, F. Delaglio & A. Bax, J. Magn. Reson. 131, 373-378 (1998)
2. F. Cordier, A.J. Dingley & S. Grzesiek, J. Biomol. NMR 13, 175-180 (1999)


hsqcf3gpiaphwg



hsqcf3gpiaphwg.2



hsqcf3gpiaphsiwg



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D TROSY EXPERIMENTS

Experiment Description

The TROSY (Transverse-Relaxation Optimized Spectroscopy) pulse scheme can be understood as a spin-state-edited HSQC experiment in which only one of the four components of a fully-F1,F2-coupled multiplet is observed.

Sample Requirements

TROSY experiments are usually applied in large biomolecules labeled with ^{15}N .

Hardware Requirements

TROSY experiments are recorded in triple-resonance inverse probeheads equipped with gradients.

NMR Spectrum

The TROSY experiment afford a conventional ^1H -Xchemical shift correlation map

Related Experiments

2D HSQC

Also see 3D TROSY-like triple-resonance NMR experiments applied on labeled proteins.

References:

1. M. Czisch & R. Boelens, *J. Magn. Reson.* 134, 158-160 (1998)
2. K. Pervushin, G. Wider & K. Wuethrich, *J. Biomol. NMR* 12, 345-348 (1998)
3. A. Meissner, T. Schulte-Herbrueggen, J. Briand & O.W. Sorensen, *Mol. Phys.* 96, 1137-1142 (1998)
4. J. Weigelt, *J. Am. Chem. Soc.* 120, 10778-10779 (1998)
5. M. Rance, J.P. Loria & A.G. Palmer III, *J. Magn. Reson.* 136, 91-101 (1999)
6. G. Zhu, X.M. Kong & K.H. Sze, *J. Biomol. NMR* 13, 77-81 (1999)

2D TROSY Experiments

From f2 channel:

Phase-sensitive ge -2D TROSY with presaturation (**trosygpphpr**)
 ge -2D TROSY for aromatic residues with WATERGATE (**trosyangpphwg**)


From f3 channel:

Phase-sensitive ge -2D ^1H - ^{15}N TROSY



- Using echo-antiecho (**trosyetf3gpsi** | TROSYETF3GPSI)
- Using echo-antiecho (**trosyetf3gpsi.2**)
- Using echo-antiecho and water flip-back (**trosyetf3gpsi**)
- Using echo-antiecho and different phase cycling (**trosyetf3gpsi2**)
- Using echo-antiecho and different phase cycling to give IPAP TROSY (**trosyetf3gpsi**)
- Using echo-antiecho and different phase cycling to give IPAP TROSY (**trosyetf3gpsi.2**)
- Using WATERGATE (3-9-19) (**trosyf3gpph19** | TROSYF3GPPH19)
- Using WATERGATE and improved sensitivity (**trosyf3gpphsi19** | TROSYF3GPPHSI19)


Phase-sensitive ge -2D ^1H - ^{15}N ZQ-TROSY using WATERGATE (**trosyzqpphwg**)

Also see:
Spin-State Edited HSQC-IPAP experiments
CRINEPT Experiments






PEP vs TROSY and half-TROSY






1. D. Nietlispach, J. Biomol. NMR 31, 161-166 (2005)
2. D.W. Yang & L.E. Kay, J. Biomol. NMR 13, 3-10 (1999)







trocyf3gpqh19




trocyf3gpqhsi19




trocyf3gpqhsi19.2




trosetf3gpiasi




trosetf3gpiasi.2




trosetzqgpphwg




;K.V. Pervushin, G. Wider, R. Riek & K. Wuethrich, PNAS, 96, 9607-9612 (1999)



trosetf3gpsi




trosetf3gpsi.2




T. Schulte-Herbrueggen & O.W. Sorensen, *J. Magn. Reson.* 144, 123 - 128 (2000)

trosetf3gpsi2



trosetf3gpsi



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D CRINEPT EXPERIMENT

Experiment Description:

The CRINEPT (Cross-Correlated Relaxation Enhanced INEPT) experiment has been specifically designed to observe NH correlation peaks in very large ^{15}N -labeled proteins.


References:

1. R. Riek, G. Wider, K. Pervushin & K. Wuethrich, Proc. Natl. Acad. Sci. USA 96, 4918-4923 (1999)

2D CRINEPT Experiments

ge-2D ^1H - ^{15}N CRINEPT using flip-back (`crineptgp`)

Also see: TROSY experiments



use AU-program `splitcrinept` to process data

BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D IDIS-HSQC EXPERIMENT


IDIS-HSQC: Isotopically Discriminated HSQC experiment

References:


1. A.P. Golovanov, R.T. Blankley, J.M. Avis & W. Bermel, *J. Am. Chem. Soc.* 129, 6528-6535 (2007)
2. W. Bermel, E.N. Tkach, A.G. Sobol & A.P. Golovanov, *J. Am. Chem. Soc.* 131, 8564-8570 (2009)

2D IDIS-HSQC Experiments


- Phase-sensitive 2D ^1H - ^{15}N IDIS-HSQC for isotopical discrimination between $^{15}\text{N}(^{13}\text{CO})$ and $^{15}\text{N}(^{12}\text{CO})$ using WATERGATE (**hsqcf3gpiddphwg**)
- Phase-sensitive 2D ^1H - ^{15}N IDIS-TROSY for isotopical discrimination between $^{15}\text{N}(^{13}\text{CO})$ and $^{15}\text{N}(^{12}\text{CO})$ using WATERGATE (**troisyf3gpiddphwg**)
- Phase-sensitive 2D ^1H - ^{15}N IDIS-TROSY-IPAP for isotopical discrimination between $^{15}\text{N}(^{13}\text{CO})$ and $^{15}\text{N}(^{12}\text{CO})$ using WATERGATE (**troisyf3gpiaiddphwg**)



trocyf3gpidphwg



trocyf3gpiaidphwg



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D HMQC-COSY EXPERIMENTS

2D HMQC-COSY Experiments

- **Phase cycled:**

Magnitude-mode 2D HMQC-COSY using BIRD (`hmqcbindqfrl`)
 Phase-sensitive 2D HMQC-COSY using BIRD with decoupling (`hmqcbiphrl`)
 Phase-sensitive 2D HMQC-COSY using BIRD without decoupling (`hmqcbindphrl`)
 Phase-sensitive 2D HMQC-COSY-DQF using BIRD (`hmqcbidfphrl`)



- **Gradient-enhanced:**

Magnitude-mode ge-2D HMQC-COSY (`hmqcgpqfrl`)
 H2BC experiment with a three-low-pass filter (`h2bcetgpl3`)




- **H2BC experiment:**

H2BC experiment with a three-low-pass filter (`h2bcetgpl3`)
 Multiplicity-edited H2BC experiment with a three-low-pass filter (`h2bcdetgpl3`)
 HAT-HMBC (H2BC + edited HMBC) experiment with a three-low-pass filter (`hathmbcetgpl3`)

Also see HMQC and HMQC-TOCSY experiments



1. L. Lerner & A. Bax, J. Magn. Reson. 69, 375-380 (1986)
2. W. Willker, D. Leibfritz, R. Kerssebaum & W. Bermel, Magn. Reson. Chem. 31, 287-292 (1993)



H2BC:



1. N.T. Nyberg, J.O. Duus & O.W. Soerensen, *J. Am. Chem. Soc.* 127, 6154-6155 (2005)

Edited-H2BC:


2. N.T. Nyberg, J.O. Duus & O.W. Soerensen, *Magn. Reson. Chem.* 43, 971-974 (2005)

HAT-HMBC:



3. A.J. Benie & O.W. Soerensen, *J. Magn. Reson.* 184, 315-321 (2007)



H2BC:
 Two-bond correlations
 but no signals from
 quaternary
 carbons




Comparison between HMBC(blue) and H2BC (red) allows differentiation 2-bond vs 3-bond connectivities



Processing:
use AU-program split [ipap 2] to create separate datasets

hathmbcetgp13



Processing:
 use AU-program split [ipap 2] to create separate datasets

BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D HMQC-TOCSY EXPERIMENTS

Experiment Description:

The 2D HMQC-TOCSY experiment is an hybrid experiment consisting of a first refocused HMQC pulse train (optimized to $d_2=1/2J(XH)$) followed by a mixing TOCSY process defined by d9. Other experimental set-up as usually made in HMQC experiments.

Two type of peaks are obtained in a conventional 2D correlation map : i) Direct X-H correlations and ii) Relayed correlations connecting each protonated X nucleus with all 1H belonging to the same spin system.

See 2D TOCSY and 2D HSQC-TOCSY for related experiments.

2D HMQC-TOCSY Experiments

- Phase-cycled:

Phase-sensitive 2D HMQC-TOCSY (**hmqcmlph**)

Phase-sensitive 2D HMQC-TOCSY without decoupling (**hmqcmlndph**)

Magnitude-mode 2D HMQC-TOCSY using BIRD (**hmqcbimlqf**)

Magnitude-mode 2D HMQC-TOCSY using BIRD without decoupling (**hmqcbimlndqf**)

Phase-sensitive 2D HMQC-TOCSY using BIRD (**hmqcbimlph**)

Phase-sensitive 2D HMQC-TOCSY using BIRD without decoupling (**hmqcbimlndph**)

- Phase-cycled and solvent suppression:

Phase-sensitive 2D HMQC-TOCSY with presaturation (**hmqcmlphpr**)

Phase-sensitive 2D HMQC-TOCSY with presaturation and without decoupling (**hmqcmlndphpr**)

- Gradient-enhanced from f2 channel:

Magnitude-mode ge-2D HMQC-TOCSY with MLEV (**hmqcgpmlqf** | **HMQC6PML**)

Phase-sensitive ge-2D HMQC-TOCSY with DIPSI-2 using echo-antiecho (**hmqcdietgp**)

Phase-sensitive ge-2D HMQC-TOCSY with DIPSI-2 using PEP (**hmqcdietgpsi**)

Phase-sensitive ge-2D HMQC-TOCSY with DIPSI-2 using PEP using shorter overall timing (**hmqcdietgpsi.2**)

- Gradient-enhanced from f3 channel:

Phase sensitive ge-2D ^1H - ^{15}N HMQC-TOCSY with DIPSI-2 using echo-antiecho (**hmqcdietf3gp**)

Phase sensitive ge-2D ^1H - ^{15}N HMQC-TOCSY with DIPSI-2 using PEP (**hmqcdietf3gpsi**)

Phase sensitive ge-2D ^1H - ^{15}N HMQC-TOCSY with DIPSI-2 using PEP and shorter overall timing (**hmqcdietf3gpsi.2**)

Also see in 3D HMQC-TOCSY experiments


3D ^1H - ^{13}C HMQC-TOCSY experiment using BIRD (**hmqcmlvbi3d**)

Also see:


2D HMQC and 2D TOCSY experiments

2D HSQC-TOCSY experiments


hmqcmlndph




hmqcmlph



hmqcmlndphpr




hmqcmlphpr




hmqcbimlqf

hmqcbimlph






hmqcbimlndqf

hmqcbimlndph




The delay $d9$ defines the TOCSY mixing period, independently if DIPSI-2 or MLEV are used. Typical values range from 0 to 120 ms and transfer efficiency depends on spin system topologies.




1. L. Lerner & A. Bax, *J. Magn. Reson.* 69, 375-380 (1986)
2. A.G. Palmer III, J. Cavanagh, P.E. Wright & M. Rance, *J. Magn. Reson.* 93, 151-170 (1991)


hmqcdietgp




hmqcdietgps




hmqcdietgps.2




hmqcdietf3gp



hmqcdietf3gpsi



hmqcdietf3gpsi.2



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D HMQC-ROESY EXPERIMENTS

Experiment Description:

The 2D HMQC-ROESY experiment is an hybrid experiment consisting of a first refocused HMQC pulse train (optimized to $d2=1/2J(XH)$) followed by a mixing ROESY process defined by p15. Other experimental set-up as usually made in HMQC experiments.

Two type of peaks are obtained in a conventional 2D correlation map : i) Direct X-H correlations and ii) ROE correlations connecting a protonated X nucleus with a 1H resonance which is close to the directly-attached ^1H-X proton.

The experiment can be recorded without X-decoupling during 1H acquisition (pl12 or pl16 set to 120dB). This is useful to observe NOE between degenerate protons (for instance, symmetrical molecules).

See 2D ROESY and 2D HSQC-ROESY for related experiments.

References:

A. Bax & D.G. Davis, J. Magn. Reson 63, 207-213 (1985).

2D HMQC-ROESY Experiments

- Gradient-enhanced from the f2 channel

Phase-sensitive ge-2D HMQC-ROESY using echo-antiecho (**hmqcetgpro**)

Phase-sensitive ge-2D HMQC-ROESY with T-ROESY using echo-antiecho (**hmqcetgpro.2**)

- Gradient-enhanced from the f3 channel

Phase-sensitive ge-2D $^1H-^{15}N$ HMQC-ROESY using echo-antiecho (**hmqctf3gpro**)


Phase-sensitive ge-2D $^1H-^{15}N$ HMQC-ROESY with T-ROESY using echo-antiecho (**hmqctf3gpro.2**)

Also see:


2D HMQC and 2D ROESY Experiments

2D HMQC-NOESY experiments


The pulse p15 (in microseconds), applied at power level pl11, defines the mixing period in all ROESY experiments.




hmqcetgpro




hmqcetgpro.2






hmqcetf3gpro



hmqcetf3gpro.2



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D HMQC-NOESY EXPERIMENTS

Experiment Description:

The 2D HMQC-NOESY experiment is an hybrid experiment consisting of a first refocused HMQC pulse train (optimized to $d2=1/2J(XH)$) followed by a mixing NOESY building block defined by d8. Other experimental set-up as usually made in HMQC experiments.

Two type of peaks are obtained in a conventional 2D correlation map : i) Direct X-H correlations and ii) ROE correlations connecting a protonated X nucleus with a 1H resonance which is close to the directly-attached 1H-X proton.

The experiment can be recorded without X-decoupling during 1H acquisition (pl12 or pl16 set to 120dB). This is useful to observe NOE between degenerate protons (for instance, symmetrical molecules).

See 2D NOESY, 2D HOESY and 2D HSQC-NOESY for related experiments.

2D HMQC-NOESY Experiments

- Phase cycled:

Phase-sensitive 2D HMQC-NOESY with presaturation (**hmqcnophpr**)

Phase-sensitive 2D HMQC-NOESY using BIRD (**hmqcbinoph**)

- Gradient-enhanced from the f2 channel:

Phase-sensitive ge-2D HMQC-NOESY using echo-antiecho (**hmqcetgpno**)

- Gradient-enhanced from the f3 channel:

Phase-sensitive ge-2D ^1H - ^{15}N HMQC-NOESY using echo-antiecho (**hmqcetf3gpno**)

Also see in "3D HMQC-NOESY experiments":


3D ^1H - ^{13}C HMQC-NOESY experiment using BIRD (**hmqcnoesybi3d**)

Also see:


2D HMQC and 2D NOESY Experiments

2D HMQC-ROESY experiments


The delay d8 defines the mixing period in all NOESY experiments,




hmqcnophpr




hmqcbinoph






hmqcetgpno



hmqcetf3gpno



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D HSQC-TOCSY EXPERIMENTS

Experiment Description:


The 2D HSQC-TOCSY experiment is an hybrid experiment consisting of a first refocused HSQC pulse train (optimized to $d_4=1/4J(XH)$) followed by a mixing TOCSY process defined by d_9 . Other experimental set-up as usually made in regular HSQC experiments.

Two type of peaks are obtained in a conventional 2D correlation map : i) Direct X-H correlations and ii) Relayed correlations connecting each protonated X nucleus with all 1H belonging to the same spin system.

Versions of the basic experiment include:

- i) Use of MLEV or DIPSI as a TOCSY mixing process.
- ii) Use of gradients for coherence selection. *OPTIONAL* combination with sensitivity-improved building block (PEP).
- iii) Use of adiabatic ^{13}C 180 pulses.
- iv) Use of editing of direct responses in order to differentiate direct and relayed correlations .
- v) Use of multiplicity editing to know the multiplicity of the carbon resonance.
- vi) *OPTIONAL* use of the f2 and f3 channel version, for triple resonance probeheads.

Also see 2D TOCSY, 2D HMQC-TOCSY and 3D HSQC-TOCSY for related experiments.



2D 1H - ^{13}C HSQC-TOCSY spectrum of strychnine

2D HSQC-TOCSY Experiments

- **Gradient-enhanced from the f2 channel**

Phase sensitive ge-2D HSQC-TOCSY with MLEV using z-filter (**hsqcgpmlph** | HSQC6PMLPH)

Phase-sensitive ge-2D HSQC-TOCSY with MLEV using echo-antiecho (**hsqcetgpml** | HSQCET6PML)

Phase-sensitive ge-2D HSQC-TOCSY with DIPSI-2 using PEP (**hsqcdietgps**)

Phase-sensitive ge-2D HSQC-TOCSY with DIPSI-2 using PEP and adiabatic inversion pulses (**hsqcdietgpsisp** | HSQCDIET6PSISP)

Phase-sensitive ge-2D HSQC-TOCSY with DIPSI-2 using PEP and adiabatic inversion and refocusing pulses (**hsqcdietgpsisp.2**)

- **Gradient-enhanced with editing from the f2 channel**

Phase sensitive ge-2D HSQC-TOCSY using PEP with editing of multiplicity (**hsqcdiedetgpsisp.1**)

Phase sensitive ge-2D HSQC-TOCSY using PEP with editing of direct responses (**hsqcdiedetgpsisp.2**)

Phase sensitive ge-2D HSQC-TOCSY using PEP with editing of multiplicity and direct responses (**hsqcdiedetgpsisp.3**)

- **Gradient-enhanced from the f3 channel**

Phase sensitive ge-2D ^1H - ^{15}N HSQC-TOCSY with MLEV using echo-antiecho (**hsqcetf3gpml**)

Phase sensitive ge-2D ^1H - ^{15}N HSQC-TOCSY with DIPSI-2 using PEP (**hsqcdietf3gps** | HSQCDIETF3GPSI)

Also see:

HSQC-TOCSY type experiments (HETLOC, HECAD, spin-edited HSQC-TOCSY) for $^1\text{J}_{\text{CH}}$ measurements

3D HSQC-TOCSY experiments:

Phase-sensitive ge-2D HSQC-TOCSY with MLEV using echo-antiecho (**hsqcetgpml3d.2**)


Phase-sensitive ge-2D HSQC-TOCSY with DIPSI-2 using PEP and adiabatic inversion and refocusing pulses (**hsqcdietgpsisp3d.2**)

Also see:


2D HSQC and 2D TOCSY experiments

2D HMQC-TOCSY Experiments


The delay d9 defines the TOCSY mixing period, independently if DIPSI-2 or MLEV are used. Typical values range from 0 to 120 ms and transfer efficiency depends on spin system topologies.





hsqcgpmiph



hsqcetgpmi




hsqc dietgpsi




Expansion region of the 2D ^1H - ^{13}C HSQC-TOCSY spectrum of ubiquitine


hsqcdietgpsisp




hsqcdietgpsisp.2




hsqcdiedetgpsisp.1




hsqcdiedetgpsisp.2





hsqcdiedetgpsisp.3







2D ^1H - ^{13}C HSQC-TOCSY spectrum of strychnine with editing of direct responses




Column expansions of the 2D ^1H - ^{13}C HSQC-TOCSY spectrum of sucrose with several combinations of editing of direct responses and/or editing of carbon multiplicity



hsqcetf3gpml



hsqcdietf3gpsi



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D HSQC-ROESY EXPERIMENT

Experiment Description:

The 2D HSQC-ROESY experiment is a hybrid experiment consisting of a first refocused HSQC pulse train (optimized to $d4=1/4J(XH)$) followed by a mixing ROESY process defined by p15. Other experimental set-up as usually made in HMQC experiments.

Two type of peaks are obtained in a conventional 2D correlation map : i) Direct X-H correlations and ii) ROE correlations connecting a protonated X nucleus with a 1H resonance which is close to the directly-attached ^1H-X proton.

The experiment can be recorded without X-decoupling during 1H acquisition (pl12 or pl16 set to 120dB). This is useful to observe NOE between degenerate protons (for instance, symmetrical molecules).

See 2D ROESY and 2D HMQC-ROESY for related experiments.

References

A. Bax & D.G. Davis, J. Magn. Reson 63, 207-213 (1985)

2D HSQC-ROESY Experiments

- Gradient-enhanced from the f2 channel

Phase-sensitive *ge*-2D HSQC-ROESY using echo-antiecho and adiabatic pulses (**hsqcetgprosp** | **HSQCET#PROSP**)

Phase-sensitive *ge*-2D HSQC-ROESY using echo-antiecho and adiabatic pulses with T-ROESY(**hsqcetgprosp.2**)

- Gradient-enhanced from the f3 channel

Phase-sensitive *ge*-2D 1H - ^{15}N HSQC-ROESY using echo-antiecho (**hsqcetf3gpro**)




Phase-sensitive *ge*-2D 1H - ^{15}N HSQC-ROESY with T-ROESY using echo-antiecho (**hsqcetf3gpro.2**)


Also see:

2D HSQC and 2D ROESY


2D HMQC-ROESY experiments

The pulse p15 (in microseconds), applied at power level pl11, defines the mixing period in all ROESY experiments.






hsqcetf3gpro



hsqcetf3gpro.2



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D HSQC-NOESY EXPERIMENT

Experiment Description:

The 2D HSQC-NOESY experiment is an hybrid experiment consisting of a first refocused HSQC pulse train (optimized to $d4=1/4J(XH)$) followed by a mixing NOESY building block defined by $d8$. Other experimental set-up as usually made in HSQC experiments.

Two type of peaks are obtained in a conventional 2D correlation map : i) Direct X-H correlations and ii) ROE correlations connecting a protonated X nucleus with a $1H$ resonance which is close to the directly-attached $1H-X$ proton.

The experiment can be recorded without X-decoupling during $1H$ acquisition (pl12 or pl16 set to 120dB). This is useful to observe NOE between degenerate protons (for instance, symmetrical molecules).

See 2D NOESY and 2D HMQC-NOESY for related experiments.

2D HSQC-NOESY Experiments

- **Gradient-enhanced from the f2 channel**

Phase-sensitive *ge*-2D HSQC-NOESY using echo-antiecho and adiabatic pulses (**hsqcetgpno**sp | **HSQCETGPNO**SP)

- **Gradient-enhanced from the f3 channel**

Phase-sensitive *ge*-2D 1H - ^{15}N HSQC-NOESY using echo-antiecho (**hsqcetf3gpno** | **HSQCETF3GPNO**)
Phase-sensitive *ge*-2D 1H - ^{15}N HSQC-NOESY using XY16 and WATERGATE (**hsqcf3gpnowgxy** | **HSQCF3GPNOWGXY**)

Also see in "Selective 1D Experiments":

Selective *ge*-1D HSQC-NOESY experiment (**selhsqcgpnosp**)



Selective *ge*-1D HSQC-NOESY experiment without decoupling (**selhsqcgpndnosp**)


Also see:

2D HSQC and 2D NOESY


2D HSQC-ROESY Experiments

The delay $d8$ defines the mixing period in all NOESY experiments,






hsqcetf3gpmo



hsqcf3gpnogxy



1. L. Mueller, P. Legault & A. Pardi, *J. Am. Chem. Soc.* 117, 11043-11048 (1995)
2. F.A.A. Mulder, C.A.E.M. Spronk, M. Slijper, R. Kaptein & R. Boelens, *J. Biomol. NMR* 8, 223-228 (1996)

BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D HMBC EXPERIMENTS

Experiment Description

The HMBC (Heteronuclear Multiple-Bond Correlation) experiment is a proton-detected experiment designed to obtain long-range (usually two-bond and three-bond) heteronuclear correlations between ^1H and X heteronuclei via the scalar coupling constant ($2J(\text{XH})$, $3J(\text{XH} \dots)$).

The experiment can be applied on any heteronucleus, typically ^{13}C , ^{15}N , ^{29}Si , ^{31}P ...

Sample Requirements

HMBC experiments can be recorded on any type of sample.

Hardware Requirements

HMBC experiments can be recorded on any probehead but an inverse probe equipped with gradients is strongly recommended.

NMR Spectrum

A 2D HMBC experiment yields a typical ^1H -X map that correlates ^1H and X chemical shifts via $nJ(\text{XH})$ ($n>1$). However, residual one-bond connectivities can be also present as a large doublets ($1J$ around 130-160 Hz). Versions that use low-pass filters can be used to minimize them.

The most important parameter to optimize is the **d6 delay** that allows the evolution to $nJ(\text{CH})$ in order to achieve antiphase magnetization. Usually, a d6 value of 60-65 ms is a good starting point for routine applications.

Related Experiments

See modified 2D HMBC and long-range HSQC (HSQMBC) pulse sequences to measure quantitatively the value of long-range proton-carbon coupling constants, $nJ(\text{XH})$.

Also see other experiments that provide long-range heteronuclear correlations such as 2D COLOC, 2D ADEQUATE ...

References:

1. A. Bax & M.F. Summers, J. Am. Chem. Soc. 108, 2093 - 2094 (1986)

2D HMBC Experiments

- **Phase cycled:**

Magnitude-mode 2D HMBC using low-pass J-filter (**hmbclpndqf** | HMBCLPND)
Magnitude-mode 2D HMBC with presaturation (**hmbcndprqf**)
Magnitude-mode 2D HMBC with off-resonance presaturation (**hmbcndpsqf**)

- **Gradient-based from f2 channel:**

Magnitude-mode ge-2D HMBC (**hmbcgpndqf** | HMBCGPND)
Magnitude-mode ge-2D HMBC using low-pass J-filter (**hmbcgplpndqf** | HMBCGPLPND)
Magnitude-mode ge-2D HMBC using double low-pass J-filter (**hmbcgpl2ndqf**)

Phase-sensitive ge-2D HMBC using echo-antiecho (**hmbcetgpnd**)
Phase-sensitive ge-2D HMBC using a two-fold low-pass J-filter (**hmbcetgpl2nd**)
Phase-sensitive ge-2D HMBC using a three-fold low-pass J-filter (**hmbcetgpl3nd**)
Phase-sensitive ge-2D Multiplicity-edited HMBC using a three-fold low-pass J-filter (**hmbcetgpl3nd**)

ge-2D Constant-time HMBC (CT-HMBC) using echo-antiecho (**hmbcctetgpnd**)
ge-2D Constant-time HMBC (CT-HMBC) using echo-antiecho and a two-fold low-pass J-filter (**hmbcctetgpl2nd**)

Magnitude-mode band-selective ge-2D HMBC without decoupling (**shmbcgpndqf**)

Magnitude-mode CIGAR-HMBC without decoupling (**hmbcacgplpndqf**)
Magnitude-mode CIGAR-HMBC with decoupling (**hmbcacgplpqf**)

ge-2D 2J,3J HMBC, STAR-HMBC (**hmbcacbigpl2ndqf**)
ge-2D HSMC (**hmscetgpnd**)


- **Gradient-based from f3 channel:**

Magnitude-mode ge-2D HMBC (**hmbcf3gpndqf**)
Magnitude-mode ge-2D HMBC using low-pass J-filter (**hmbcf3gplpndqf**)
Phase-sensitive ge-2D HMBC using echo-antiecho (**hmbcetf3gpnd**)
Phase-sensitive ge-2D HMBC using a three-fold low-pass J-filter (**hmbcetf3gpl3nd**)


- **Gradient-based simultaneous CN-HMBC**

Magnitude-mode ge-2D simultaneous ^1H - ^{15}N / ^1H - ^{13}C HMBC (**hmbcetgpl2ndsm**)


Also see in "Measurement of long-range proton-carbon coupling constants":
ge-2D J-HMBC using a two-fold low-pass J-filter (**hmbcetgpjcl2nd**)
Long-Range HSQC Experiments (HSQMBC)





hmbclpndqf




hmbcndprqf



hmbcndpsqf



2D 1H-13C HMBC spectrum of strychnine




NMR Building Block: A Low-Pass Filter

```

...
1 ze
2 d1
3 p1 ph1
  d2
  p3:f2 ph3
  d6
  p3:f2 ph4
.....

;d2 : 1/(2J)XH
;d6 : delay for evolution
;of long range couplings

```



NMR Building Block: A Two-Fold Low-Pass Filter


```

"DELTA1=1s/(2 * cnst6)-p16-d16"
"DELTA2=1s/(2 * cnst7)-p16-d16"
"DELTA3=d6-p16-d16-4u"


...
3 p1 ph1
  DELTA1 UNBLKGRAD
  p16:gp4
  d16
  p3:f2 ph3
  DELTA2
  p16:gp5
  d16
  p3:f2 ph3
  4u
  p16:gp6
  d16
  DELTA3
  p3:f2 ph4
.....

;cnst6: = 1J(XH)min
;cnst7: = 1J(XH)max
;cnst13: = J(XH) long range
;gpz4: 15%
;gpz5: -10%
;gpz6: -5%

```




hmbcgpndqf
hmbcf3gpndqf




gp 1 : gp 2 : gp 3
50 : 30 : 40.1 for C-13
70 : 30 : 50.1 for N-15

Use the AU program **gradratio** for calculation_of gradient ratios for common inverse gradient pulse programs.


hmbcgp1pndqf
hmbcf3gp1pndqf




hmbcgp12ndqf



hmbcgp1pndprqf







2D 1H-15N HMBC spectrum of strychnine

Phase-Sensitive HMBC:


For Measurement of nJ_{CH} .

D.O. Cicero, G. Barbato & R. Bazzo, *J. Magn. Reson.* 148, 209-213 (2001)


hmbcetgpnd




hmbcetf3gpnd




hmbcetgpl2nd



hmbcetgpl3nd




hmbcetf3gpl3nd



Edited-HMBC:


N.T. Nyberg & O.W. Soerensen, Magn. Reson. Chem. 44, 451-454 (2006)




Processing:
 use AU-program split [ipap 2] to split data

CT-HMBC:

T.D.W. Claridge & I.Perez-Victoria, *Org. Biomol. Chem.* 1, 3632-3534 (2003)




Band-Selective HMBC




CIGAR-HMBC:

1. C.E. Hadden, G.E. Martin & V.V. Krishnamurthy, *Magn. Reson. Chem.* 38, 143-147 (2000)
2. C.E. Hadden, G.E. Martin & V.V. Krishnamurthy, *J. Magn. Reson.* 140, 274-280 (1999)
3. R. Wagner & S. Berger, *Magn. Reson. Chem.* 36, S44-S46 (1998).


hmbcacgplpqf



hmbcacgplpndqf




hmbcacgplpndqf.2



50:30:40.1:5:60:-40:-20:17


2J, 3J-HMBC or STAR-HMBC:

V.V. Krishnamurthy, D.J. Russel, C.E. Hadden & G.E. Martin,
 J. Magn. Reson. 146, 232-239 (2000)



HMSC:


R. Burger, C. Schorn & P. Bigler, J. Magn. Reson. 148, 88-94 (2001)



Processing:
 use AU-program split [ipap 2] to split and recombine data

Simultaneous CN-HMBC:

M. Perez-Trujillo, P. Nolis, W. Bermel & T. Parella, *Magn. Reson. Chem.* 45, 325-9 (2007)



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D EXPERIMENTS TO MEASURE
LONG-RANGE PROTON-CARBON
COUPLING CONSTANTS

2D Experiments to measure $^nJ_{CH}$ measurement

- **J- HMBC experiments**

ge-2D J-HMBC using a two-fold low-pass J-filter (**hmbcetgpjcl2nd**)

- **Long-range optimized ge-2D HSQC**

Phase-sensitive ge-2D long-range optimized HSQC (HSQMBC) (**hsqcetgplrsp**)

Phase-sensitive ge-2D long-range optimized HSQC using G-BIRD (GBIRD-HSQMBC)

(**hsqcetgpjclrmd**)

Phase-sensitive ge-2D long-range optimized HSQC using CPMG-XY16 (GBIRD-HSQMBC)

(**hsqcetgpjclrmdxy**)

ge-2D long-range optimized J-HSQC (EXSIDE) (**hsqcetgplrjrcsp**)

- **ge-2D HSQC-TOCSY type experiments**

ge-2D w1-filtered TOCSY using DIPSI-2 (HETLOC) (**dipsi2etgpjcsix1**)

Phase-sensitive ge-2D HSQC-HECADE (**hsqcdietgpjcnndsisp**)

Phase-sensitive ge-2D F2-spin-edited HSQC-TOCSY using DIPSI-2 and sensitivity improvement

(**hsqcdietgpiasisp**)

- **2D Selective J-Resolved experiments**

Magnitude-mode proton-selective 2D Heteronuclear J-Resolved (**seljresqfsp**)

Also see in 2D HMBC Experiments:

Phase-sensitive ge-2D HMBC using echo-antiecho (**hmbcetgpnd**)

Phase-sensitive ge-2D HMBC using a two-fold low-pass J-filter (**hmbcetgpl2nd**)

Phase-sensitive ge-2D HMBC using a three-fold low-pass J-filter (**hmbcetgpl3nd**)

Phase-sensitive ge-2D HMBC using echo-antiecho from f3 channel (**hmbcetf3gpnd**)

Phase-sensitive ge-2D HMBC using a three-fold low-pass J-filter from f3 channel

(**hmbcetf3gpl3nd**)

Phase-sensitive ge-2D Multiplicity-edited HMBC using a three-fold low-pass J-filter

(**hmbcetgpl3nd**)


ge-2D Constant-time HMBC (CT-HMBC) using echo-antiecho (**hmbcctetgpnd**)

ge-2D Constant-time HMBC (CT-HMBC) using echo-antiecho and a two-fold low-pass J-filter


(**hmbcctetgpl2nd**)

Direct measurement on Multiplets in F2-detected Dimension

HSQC-TOCSY




HMBC



Others: HMQC-TOCSY, HECADE, HETLOC
 Two steps: $^1J_{CH}$ + J_{HH}
 Only for protonated carbons
 In-phase Magnetization

Others: HSQMBC, EXSIDE, J-HMBC
 A single step: $^nJ_{CH}$
 For all carbons
 Anti-phase Magnetization



Measurement from spin-state selective patterns:

alpha vs beta in F2
 TROSY/antiTROSY pattern
 ECOSY pattern


Direct measurement on Multiplets in F1-indirect Dimension

J-Resolved Experiments
 J-HMBC Experiments: Amplification factor




J-HMBC:

A. Meissner & O.W. Soerensen, Magn. Reson. Chem. 39, 49-52 (2001)

J-scaling factor of N
 N depends of sw and td1




hmbcetgpjcl2nd




EXSIDE:

V.V. Krishnamurthy, J. Magn. Reson., Series A 121, 33-41 (1996)

J-scaling factor of N



hsqcetgpljrcsp




cnst16: = J(scale) factor


p12: f1 channel - 180 degree shaped pulse (selective)

HSQMBC:


1. B.L. Marquez, W.H. Gerwick & R.T. Williamson, *Magn. Reson. Chem.* 39, 499-530 (2001)
2. R.T. Williamson, B.L. Marquez, W.H. Gerwick & K.E. Kover, *Magn. Reson. Chem.* 38, 265-273 (2000)



HSQMBC-BIRD




HSQMBC-BIRD-CPMG






80:20.1:17:11

HETLOC:

D. Uhrin, G. Batta, V.J. Hruby, P.N. Barlow & K.E. Kover, *J. Magn. Reson.* 130, 155-161 (1998)




dipsi2etgpjcsix1




HSQC-HECADE:



W. Kozminski & D. Nanz, J. Magn. Reson. 142, 294-299 (2000)



hsqcdietgpjncdisp





cnst16: = J(scale) factor





Spin-state Selective HSQC-TOCSY:

P. Nolis & T. Parella, J. Magn. Reson., 176, 15-26 (2005)




Processing:

use AU-program split [ipap 2] to process data




Selective J-Resolved

R. Freeman and A. Bax, JACS, 1082 (1982)



Selective J-Resolved spectrum after selective inversion of the H1 proton of sucrose



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

2D ADEQUATE EXPERIMENTS

Experiment Description:

The ADEQUATE experiment allows to connect ^1H and ^{13}C separated by two- or more bonds. The basic pulse sequence consists of an INADEQUATE pulse scheme incorporated into an HSQC experiment. Several versions are available as a function how the experiment/delays are optimized. Thus, an 1,1-ADEQUATE experiment is equivalent to an HCC experiment in which the delays are optimized to $1J(\text{CH})$ and $1J(\text{CC})$. On the other hand, an 1,n-ADEQUATE experiment is optimized to $1J(\text{CH})$ and $nJ(\text{CC})$.

The experiment is based on the selection of double-quantum ^{13}C - ^{13}C coherences and therefore it can be considered of very low sensitivity.

Other related experiments: 2D INADEQUATE, 2D HMBC and 2D HSQC experiments

2D ADEQUATE Experiments

- 1,1-ADEQUATE:

Phase-sensitive 1,1 ADEQUATE using adiabatic pulses (`adeq11etgppsp`)

Phase-sensitive w1-refocused 1,1 ADEQUATE using adiabatic pulse (`adeq11etgprdsp`)

Phase-sensitive w1-refocused 1,1 ADEQUATE using adiabatic pulses with evolution of $J(\text{CC})$ (`adeq11etgprdsp.2`)

Phase-sensitive w1-refocused 1,1 J-ADEQUATE using adiabatic pulses for measuring $J(\text{CC})$ (`adeq11etgpjcrdsp.2`)

- 1,n-ADEQUATE:

Phase-sensitive 1,n ADEQUATE (`adeq1netgp`)

- n,1-ADEQUATE:








Phase-sensitive n,1 ADEQUATE (`adeqn1etgp`)

- n,n-ADEQUATE:



Phase-sensitive n,n ADEQUATE (`adeqnnetgp`)

Also see:

2D INADEQUATE





1. B. Reif, M. Koeck, R. Kerssebaum, H. Kang, W. Fenical & C. Griesinger *J. Magn. Reson.* A118, 282-285 (1996).
2. B. Reif, M. Koeck, R. Kerssebaum, J. Schleucher & C. Griesinger *J. Magn. Reson.* B112, 295-301 (1996)
3. M. Koeck, R. Kerssebaum & W. Bermel, *Magn. Reson. Chem.* 41, 65-69 (2003)




2D 1,1-ADEQUATE spectrum of menthol

adeq11etgprdsp



2D Refocused 1,1-ADEQUATE spectrum of menthol

adeq11etgprdsp.2





+78.4:+77.4:-59




J-Adequate:

C.M. Thiele & W. Bermel, *Magn. Reson. Chem.* 45, 889-894 (2007)

K.E. Kover, P. Forgo, *J. Magn. Reson.* 166, 47-52 (2004)




2D 1,1-J-ADEQUATE spectrum of menthol




BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

DIFFUSION/DOSY EXPERIMENTS



DOSY Spectrum:
Chemical Shift vs Diffusion Coefficient



A Description of setting up and running Diffusion/DOSY experiments can be found in the DOSY Manual (see **Help** menu in Topspin)

DOSY/Diffusion

- Conventional 1D:

- 1D Stimulated Echo experiment (STE) (**stegp1s1d**)
- 1D Stimulated Echo experiment (STE) using bipolar gradients (**stebpgp1s1d**)
- 1D LED experiment (**ledgp2s1d**)
- 1D LED experiment using bipolar gradients (**ledbpgp2s1d**)
- 1D LED experiment using bipolar gradients and presaturation (**ledbpgppr2s1d**)
- 1D Double-Stimulated Echo Experiment (DSTE) (**dste gp3s1d**)
- 1D Double-Stimulated Echo Experiment (DSTE) using bipolar gradients (**dste bpgp3s1d**)

- 1D Stimulated Echo experiment using bipolar gradients and WATERGATE (**stebpgp1s191d**)

- 1D STE-INEPT experiment (**stebpgpin1s1d**)

- 2D DOSY maps:

- 2D Stimulated Echo experiment (STE) (**stegp1s**)
- 2D Stimulated Echo experiment using bipolar gradients (**stebpgp1s**)

- 2D Double-Stimulated Echo Experiment (DSTE) (**dste gp3s**)
- 2D Double-Stimulated Echo Experiment (DSTE) using bipolar gradients (**dste bpgp3s**)
- 2D LED experiment (**ledgp2s**)
- 2D LED experiment using bipolar gradients (**ledbpgp2s**)
- 2D LED experiment using bipolar gradients and presaturation (**ledbpgppr2s**)

- 2D Stimulated Echo experiment using bipolar gradients and WATERGATE (**stebpgp1s19**)

- 2D STE-INEPT experiment (**stebpgpin1s**)

- 2D & 3D DOSY related experiments:

- 3D DOSY-COSY using LED with bipolar gradients (**ledbpgpco2s3d**)

- 2D DOSY-TOCSY with LED using bipolar gradients (**ledbpgpml2s2d**)
- 2D DOSY-TOCSY with LED using bipolar gradients and WATERGATE (**ledbpgpml2s192d**)
- 3D DOSY-TOCSY using LED with bipolar gradients (**ledbpgpml2s3d**)

- 3D DOSY-NOESY using LED with bipolar gradients (**ledbpgpno2s3d**)

Automated DOSY Data Acquisition:

Run DOSY experiment with "**dosy**" and answer the questions.

Automated DOSY Data Processing:

Transform the n 1D spectra with **xf2**


Phase the first ser and apply the phase correction with **xf2p**

Baseline correction with **abs2**


Set diffusion parameters with **setdiffparm**

Transform diffusion dimension with **dosy2d**


step1s
step1s1d




stebpgp1s
stebpgp1s1d




ledgp2s
ledgp2s1d



ledbpgp2s
ledbpgp2s1d




ledbpgppr2s
ledbpgppr2s1d



NMR Building Block: A Diffusion Filter.


Any 1D Diffusion Experiment (STE, DSTE, LED or LEDBP) can be used as a diffusion filter in nD Experiments.

In the example, the initial read 90 pulse in a COSY experiment is replaced by a LEDBP building block




ledbpgpml2s3d


ledbpgpml2s2d



ledbpgpml2s192d



ledbpgpno3s3d



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

1D & 2D SATURATION TRANSFER
DIFFERENCE (STD) EXPERIMENTS

Saturation Transfer Difference (STD) Experiments

- **1D STD:**

1D STD (**stddiff**)

1D STD with spoil (**stddiff.2**)

1D STD with spoil and T2 filter (**stddiff.3**)

- **1D STD with solvent suppression:**

1D STD using 3-9-19 WATERGATE (**stddiffgp19**)

1D STD with spoil using 3-9-19 WATERGATE (**stddiffgp19.2**)

1D STD with spoil and T2 filter using 3-9-19 WATERGATE (**stddiffgp19.3**)

1D STD using excitation sculpting (**stddiffesgp**)

1D STD with spoil using excitation sculpting (**stddiffesgp.2**)

1D STD with spoil and T2 filter using excitation sculpting (**stddiffesgp.3**)

- **2D STD-TOCSY:**

2D STD-TOCSY (**stdmlevph**)

2D STD-TOCSY using 3-9-19 WATERGATE (**stdmlevgpph19**)

2D STD-TOCSY using excitation sculpting (**stdmlevesgpph**)

- **2D STD-NOESY:**

2D STD-NOESY with T2 filter in F2 (**stdnoesygpph**)

2D STD-NOESY with T2 filter in F1 and F2 (**stdnoesygpph.2**)

2D STD-NOESY using 3-9-19 WATERGATE with T2 filter in F2 (**stdnoesygpph19**)

2D STD-NOESY using 3-9-19 WATERGATE with T2 filter in F1 and F2 (**stdnoesygpph19.2**)

2D STD-NOESY using excitation sculpting with T2 filter in F2 (**stdnoesyegpph**)

2D STD-NOESY using excitation sculpting with T2 filter in F1 and F2 (**stdnoesyegpph.2**)


- **2D STD-HSQC:**

2D STD-HSQC using echo-antiecho (**stdhsqetgpph**)

2D STD-HSQC with sensitivity-improvement (**stdhsqetgpsisp**)

References:

1. M. Mayer & B. Meyer, *Angew. Chem. Int. Ed.* 38, 1784-1788 (1999)
2. M. Mayer & B. Meyer, *Angew. Chem.* 111, 1902-1906 (1999)



NMR Element: Basic Saturation Loop for STD


```


...
6 (p13:sp13 ph1):f2
  4u
  lo to 6 times 15
....

;sp13: f2 channel - shaped pulse for saturation           [40 - 60 dB]
;p13: f2 channel - shaped pulse for saturation           [50 msec]
;d20: saturation time
;l4: l4 = number of averages = (total number of scans) / NS
;l5: loop for saturation: p13 * 15 = saturation time
;td1: number of experiments
;NBL: NBL = number of irradiation frequencies


;define FQ2LIST (irradiation frequencies)
;          (list has to be stored in "/u/exp/stan/nmr/lists/f1")
    
```

A STD Element can be incorporated just prior the read 90 pulse in any nD Experiments.
In the example, the STD element is incorporated in a TOCSY experiment.







stddiff.2




stddiff.3




stddiffgp19




stddiffgp19.2




stddiffgp19.3




stdiffesgp




stdiffesgp.2




stdiffesgp.3




stdmlevph




stdmlevgpph19




stdmlevsgpph




stdnoesygpph




stdnoesygpph.2





stdnoesygpph19



stdnoesygpph19.2





BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

CLEANEX EXPERIMENTS

CLEANEX: Clean Exchange Spectroscopy.

Water-selective experiments to study exchange processes


References:

1. T.L. Hwang, S. Mori, A.J. Shaka & P.C.M. van Zijl, *J. Am. Chem. Soc.* 119, 6203-6204 (1997)
2. T.-L. Hwang & A.J. Shaka, *J. Magn. Reson., Series A* 112 275-279 (1995)


CLEANEX Experiments

- 1D CLEANEX using 3-9-19 WATERGATE (**zgcxgp19**)
- 1D CLEANEX using excitation sculpting (**zgcxesgp**)
- 2D CLEANEX-Fast HSQC using 3-9-19 WATERGATE (**fhsqccxf3gp**)
- 2D CLEANEX-TROSY using 3-9-19 WATERGATE (**troscopyxf3gp**)


zgcxgp19



zgcxesgp



T.L. Hwang, P.C.M. van Zijl & S. Mori, J. Biomol. NMR 11, 221-226 (1998)



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide


SOLID-STATE NMR EXPERIMENTS



$$H_D^{I,S} = -\frac{\mu_0 \mu_N^2 g_I g_S}{4\pi \hbar^2} \frac{1}{r_{IS}^3} \frac{1}{2} (3\cos^2 \theta - 1) 2I_z S_z$$


⁷⁹Br spectrum
 KBr sample
 Spinning Rate 5 KHz

$\vartheta = 54,7^\circ$



Magic Angle Calibration with Kbr

Enter to Solids Line Shape Analysis Guide
by typing **solaguide**:



A Complete Description of Solid-State NMR Experiments
can be found in the **AVANCE Solids User Manual**
(see **Help** menu in Topspin)

Solid-State NMR Experiments

Basic ssNMR experiments:

- Hahn-Echo experiment to be used in MAS experiments ([hahnecho.av](#))
- 1D acquisition on X with High-power proton decoupling ([hpdec.av](#))
- 2D J-resolved experiment for MAS experiments([jres.av](#))
- CPMG experiment for quadrupolar nuclei([qcpmq.av](#) | [qcpmqall.av](#) | [qcpmqint.av](#))
- Saturation recovery T1 experiment using echo prior to detection ([satrececho1.av](#))
- Saturation recovery T1 experiment ([satrect1.av](#))
- 90-90 solid echo sequence for wide line observation ([solidecho.av](#))
- 1D solid-state APT experiment ([sostapt.av](#))
- 1D solid-state APT experiment using shaped pulse for phase modulated Lee-Goldburg decoupling for homonuclear dipolar ([sostapt_pmlg.av](#))

1D Cross-Polarization Experiments:

- Basic CP experiment ([cp.av](#))
- Basic CP experiment with flipback 90 ([cp90.av](#))
- Basic CP experiment with explicit programming of acquisition ([cpadc.av](#))
- Basic CP experiment with composite 180 pulse, for pulse determination ([cpc180.av](#))
- Basic CP experiment with Hahn echo for stationary CP ([cphahn.av](#))
- T1 inversion recovery with CP detection ([cpht1.av](#))
- Proton T1rho experiment, X detected via cross polarisation ([cpht1rho.av](#))
- Basic CP experiment with NQS, refocussed([cpnqs.av](#))
- CP editing experiment by phase inversion ([cppi.av](#))
- CP editing experiment by phase inversion ([cppircp.av](#))
- CP editing experiment by phase inversion ([cppispi.av](#))
- Basic CP experiment with SELTICS SSB-suppression ([cpseltics.av](#))
- Basic CP experiment with recoupling of CSA ([cpsuper.av](#))
- Basic CP experiment with TOtal Suppresson of Sidebands ([cptoss.av](#))
- Basic CP experiment with TOtal Suppresson of Sidebands and Non Quaternary Suppression ([cptoss_nqs.av](#))
- Basic CP experiment with variable contact time ([cpvc.av](#))
- Basic CP experiment with X T1 inversion recovery ([cpxt1.av](#))
- Basic CP experiment for T1rho measurement of X-nucleus ([cpxt1rho.av](#))
- Frequency shifted Lee-Goldburg decoupling after cp to see X-H J couplings ([fqlg.av](#))
- Basic cp experiment, contains calculations for various decoupling schemes: cw, tppm, xix, and pidedec ([vacp.av](#))

Double-CP experiments:

- Basic double CP experiment ([doubcp2d.av](#))
- Basic double CP experiment ([doubcp.av](#))

Multiple-Quantum Spectroscopy:

Using BABA

- 2D SQ-DQ correlation experiment for 1 rotor period recoupling using BABA. Standard experiment for 1H DQ spectra ([baba1.av](#))
- 2D SQ-DQ correlation experiment for 2 rotor periods recoupling using BABA for weaker dipole dipole interactions, e.g. 19F, compensated for pulse imperfections ([baba2rot.av](#))
- 2D SQ-DQ correlation experiment for 4 rotor periods recoupling using BABA. DQ experiment for weaker interaction hamiltonians, full compensated for pulse imperfections and offset effects ([baba4rot.av](#))
- 2D SQ-TQ correlation experiment for 1 rotor period recoupling using BABA. Standard experiment for 1H DQ spectra ([baba_1rp3qc.av](#))
- 2D SQ-DQ correlation experiment for 2 rotor period recoupling using BABA with cross polarization for weak dipole dipole couplings, compensated for pulse imperfections ([babacp2.av](#))
- 2D SQ-DQ correlation experiment for 4 rotor period recoupling using BABA with cross polarization for weak dipole dipole couplings, full compensated for pulse imperfections and offset effects ([babacp4.av](#))

Using POST-C7

- SQ-TQ experiment using POST_C7 sequence ([pc7_1dtqf.av](#))
- 1D DQ excitation sequence POST_C7 with cross polarization ([pc7cp1d.av](#))
- 2D SQ-DQ correlation experiment with POST_C7 sequence and cross polarization ([pc7cp2d.av](#))
- 2D SQ-DQ correlation experiment with POST_C7 and cross polarization for large sweep width ([pc7cp2dls.w.av](#))
- CS-CS correlation through dipolar couplings NOESY type sequence using POST-C7 for broader DQ excitation range ([pc7cp2dnoe.av](#))
- 1D DQ excitation sequence POST_C7 for setup of the 2D INADEQUATE type experiments ([pc71d.av](#))
- 2D SQ-DQ experiment using POST_C7 sequence ([pc72d.av](#))
- 2D SQ-DQ experiment using POST_C7 sequence ([pc72dls.w.av](#))
- 2D SQ-TQ experiment using POST_C7 sequence ([pc72dls.wtq.av](#))

Using R14_2^6

- 1D DQ excitation sequence R14_2^6 supercycled ([r14_2_1d.av](#))
- 2D SQ-DQ correlation experiment with R14_2^6 supercycled use r14_2_1d.av for setup ([r14_2_2d.av](#))
- 1D DQ excitation sequence R14_2^6 supercycled ([r14_2_cp1d.av](#))
- 2D SQ-DQ correlation through dipole dipole interaction R14_2^6 supercycled ([r14_2_cp2d.av](#))

Using SC14

- SC14 1D sequence to set up SC14 power level pl11 ([sc14cp1d.av](#))
- DQ correlation - inadequate type sequence using SC14 ([sc14cp2d.av](#))

Using SPC5

- 1D SQ-DQ correlation experiment with SPC5 sequence and cross polarization ([spc5cp1d.av](#))
- 2D SQ-DQ correlation experiment with SPC5 sequence and cross polarization ([spc5cp2d.av](#))
- 2D SQ-DQ correlation experiment with SPC5 sequence and cross polarization ([spc5cp2dls.w.av](#))

Homonuclear Correlation / Exchange Experiments:

- 2D exchange NMR in rotating solids ([cpnoesy.av](#))
- 2D INADEQUATE using cross-polarization ([cpinadequate.av](#))
- 2D exchange NMR in rotating solids ([cprfdr.av](#))
- 2D exchange NMR in rotating solids ([rfdrps.av](#))

Heteronuclear Correlation Experiments:

- 2D correlation experiment using heteronuclear multiple quantum coherences ([hmqc.av](#))
- Heteronuclear correlation between protons and X nuclei with or without homonuclear decoupling during t1 possible decoupling schemes are FSLG, PMLG, and DUMBO ([hxhetcor.av](#))
- 2D HETCOR with FLSG during t1([lghe tfq.av](#))
- 2D HETCOR with FLSG during t1([lghe tfqpi.av](#))
- 2D HETCOR with FLSG during t1([lghe tloop.av](#))
- 2D FSLG-HETCOR with shaped pulse during t1 ([lghetshape.av](#))
- 2D MAS-J-HMQC experiment with FSLG decoupling ([masjhmqc.av](#))

CA-CSA Heteronuclear Correlation Experiments:

- PASS CS-CSA experiment ([cppass.av](#))

Recoupling Experiments:

- 2D REDOR experiment ([cpredor.av](#))
- 2D REDOR experiment. Interleaved acquisition of S and S0 signal ([cpredori.av](#))
- 1D Rotational echo double resonance experiment ([cpredxy8.av](#))
- 2D Rotational echo double resonance experiment ([cpredxy82d.av](#))
- REDOR experiment ([redor.av](#))
- REDOR experiment ([redori.av](#))
- 1D rotational echo double resonance experiment ([selredor1d.av](#))
- 2D rotational echo double resonance experiment ([selredor2d.av](#))

MQMAS Experiments for quadrupolar nuclei:

- 3Q MAS pulse program for half integer spin nuclei ([mp3q.av](#))
- 3Q MAS pulse program for odd half integer spin nuclei ([mp3qdfs.av](#))
- 3Q MAS pulse program for odd half integer spin nuclei. 4-pulse experiment with z-filter ([mp3qdfs.z.av](#))
- 3Q MAS pulse program for 3/2 nuclei ([mp3qfam.av](#))
- 3Q MAS pulse program for 3/2 nuclei. 4-pulse experiment with z-filter ([mp3qfam.z.av](#))
- 3Q MAS pulse program for half integer spin nuclei. 4-pulse experiment with z-filter ([mp3qzfil.av](#))
- 3Q MAS pulse program for half integer spin nuclei and with zero quantum filter ([mp3qzqf.av](#))

Heteronuclear Experiments for quadrupolar nuclei:

- 3Q HETCOR experiment for spin-5/2, spin-7/2, and spin-9/2 nuclei. 3Q in the indirect dimension for the quadrupole nucleus and MAS in the direct dimension for the spin-1/2 nucleus ([mqhetcor.av](#))
- MQ-INEPT based 2D heteronuclear correlation experiment between spin 1/2 and half integer nucleus ([mqjhetcor.av](#))

PISEMA Experiments:

- 2D PISEMA experiment to correlate X chemical shifts with 1H-X dipolar couplings, frequency switched Lee-Goldburg during 1H-15N dipolar evolution ([pise.ma.av](#))
- 2D PISEMA experiment 2D pulse sequence to correlate X chemical shifts with 1H-X dipolar couplings and using ramped contact during FSLG period ([pise.maramp.av](#))
- PISEMA experiment based on the Magic-Sandwich Decoupling in the indirect dimension ([sammy.av](#))

STMAS Experiments:

- (satellite transition) STMAS experiment on half integer nuclei ([stmas.av](#))
- (satellite transition) STMAS experiment on half integer nuclei using double quantum filter for CT/CT signal suppression and z-filter ([stmasdqfe.av](#))
- (satellite transition) STMAS experiment on half integer nuclei using double quantum filter for CT/CT signal suppression and z-filter ([stmasdqfz.av](#))


CRAMPS Experiments:

- 1H CRAMPS, BR24 experiment ([br24.av](#))
- 1H CRAMPS, MREV-8 experiment ([mrev8.av](#))


Miscellaneous Experiments:

- Double frequency sweep followed by a selective excitation pulse ([dfs90sel.av](#))
- Pseudo-2d pulse program for acquiring quad echo data with goniometer automation. ([echogon.av](#))
- LGCR ([lgcr1.av](#))
- Scalar Coupling Driven DQF chemical shift correlation ([uc2qfcosyph.av](#))
- Pseudo-2d pulse program for acquiring quad echo data with goniometer automation ([zqgon.av](#))


hpdec.av




zggon.av



ADAMANTANE SAMPLE
Spin rate 3 KHz




jres.av




qcpmgint.av

qcpmgall.av


qcpmg.av




hahnecho.av









solidecho.av




echogon.av







satrect1.av




satrecehot1.av




cpxt1.av




cpxt1rho.av




cpht1.av





cpht1rho.av




cpnqs.av




GLYCINE SAMPLE





cpqi.av




cppircp.av
cppispi.av




sostapt.av




sostapt_pmlg.av




cptoss.av




cptoss_nqs.av




cpseltics.av




cpsuper.av




cpnoesy.av




cprfdr.av




rfdrps.av









cpinadequate.av




cppass.av







baba1.av




baba_1rp3qc.av




baba2rot.av




baba4rot.av










babacp2.av




babacp4.av







r14_2_2d.av




r14_2_1d.av




r14_2_cp2d.av




r14_2_cp1d.av




sc14cp2d.av




sc14cp1d.av







spc5cp2d.av
 spc5cp2dls.w.av




spc5cp1d.av







mp3qzqf.av




mp3q.av




mp3qzfil.av




mp3qdfs.av




mp3qdfs.z.av




mp3qfam.av




mp3qfam.z.av




mqhetcor.av




mjqhetcor.av




stmas.av




stmasdqfe.av




stmasdqfz.av




mrev8.av




br24.av




pisemaramp.av
 pisema.av




sammy.av



uc2qfcosyph.av



dfs90sel.av



BRUKER PULSE PROGRAM CATALOGUE

NMRGuide

APPENDIX

```
;Pulprog.info  
;avance-version (10/01/28)  
;  
; $CLASS=HighRes Info  
; $COMMENT=
```

```
;For a pulseprogram the first characters (usually up to 6, but  
;sometimes more) specify the type of experiment, e.g. DEPT, COSY,  
;NOESY etc.. Further properties of the pulseprogram are  
;indicated by a two-character code, which is added to the name  
;in alphabetical order. For 2D experiments the mode (absolute value,  
;phase sensitive or echo-antischo) is always indicated. H- or X-  
;decoupling is assumed to be default for heteronuclear experiments,  
;but not for homonuclear ones (except inad).  
;In case of redundant information some two-character codes may be  
;omitted.
```

```
;The two-character codes used are the following:
```

```
ac  accordion type experiment  
ad  using adiabatic spinlock  
ar  experiment for aromatic residues  
at  adiabatic TOCSY  
bi  with bird pulse for homonuclear J-decoupling  
bp  using bipolar gradients  
cc  cross correlation experiment  
cn  C13 and N15 dependent information in different indirect dimensions  
co  with COSY transfer  
cp  with composite pulse  
ct  constant time  
cv  convection compensated  
cw  decoupling using cw command  
cx  using CLEANEX_PM  
dc  decoupling using cpd command  
df  double quantum filter  
di  with DIPSI mixing sequence  
dh  homonuclear decoupling in indirect dimension  
dw  decoupling using cpd command only during wet sequence  
dq  double quantum coherence  
ea  phase sensitive using Echo/Antiecho method  
ec  with E.COSY transfer  
ed  with multiplicity editing  
es  excitation sculpting  
et  phase sensitive using Echo/Antiecho-TPPI method  
fb  using f2 - and f3 - channel  
fd  using f1 - and f3 - channel (for presaturation)  
fr  with presaturation using a frequency list  
ft  using f1 -, f2 - and f3 - channel (for presaturation)  
fh  F-19 observe with H-1 decoupling  
fp  using a flip-back pulse  
fl  for F-19 ecoupler  
fw  forward directed type experiment  
f2  using f2 - channel (for presaturation)  
f3  using f3 - instead of f2 - channel
```

f4 using f4 - instead of f2 - channel
gd gated decoupling using cpd command
ge gradient echo experiment
gp using gradients with ":gp" syntax
gr using gradients
gs using shaped gradients
hb hydrogen bond experiment
hc homodecoupling of a region using a cpd-sequence
hd homodecoupling
hf H-1 observe with F-19 decoupling
hs with homospoil pulse
ia InPhase-AntiPhase (IPAP) experiment
id IDIS - isotopically discriminated spectroscopy
ig inverse gated
ii using inverse (invi/HSQC) sequence
im with incremented mixing time
in with INEPT transfer
i4 using inverse (inv4/HMQC) sequence
jc for determination of J coupling constant
jd homonuclear J-decoupled
jr with jump-return pulse
js jump symmetrized (roesy)
lp with low-pass J-filter
lq with Q-switching (low Q)
lr for long-range couplings
l2 with two-fold low-pass J-filter
l3 with three-fold low-pass J-filter
mf multiple quantum filter
ml with MLEV mixing sequence
mq using multiple quantum
nc N15 and C13 dependent information in different indirect dimensions
nd no decoupling
no with NOESY mixing sequence
pc with presaturation and composite pulse
pg power-gated
ph phase sensitive using States-TPPI, TPPI, States or QSEQ
pl preparing a frequency list
pn with presaturation using a 1D NOESY sequence
pp using purge pulses
pr with presaturation
ps with presaturation using a shaped pulse
qf absolute value mode
qn for QNP-operation
qs phase sensitive using qseq-mode
rc for determination of residual dipolar couplings (RDC)/ J couplings
rd refocussed
rl with relay transfer
ro with ROESY mixing sequence
rs with radiation damping suppression using gradients
rt real time
ru using radiation damping compensation unit
rv with random variation
r2 with 2 step relay transfer
r3 with 3 step relay transfer

se spin echo experiment
sh phase sensitive using States et al. method
si sensitivity improved
sm simultaneous evolution of X and Y chemical shift
sp using a shaped pulse
sq using single quantum
ss spin-state selective experiment
st phase sensitive using States-TPPI method
sy symmetric sequence
s3 S3E experiment
tf triple quantum filter
tp phase sensitive using TPPI
tr using TROSY sequence
tz zeroquantum (ZQ) TROSY
ul using a frequency list
us updating shapes
wg watergate using a soft-hard-soft sequence
wt with WET watersuppression
w5 watergate using W5 pulse
xf x-filter experiments
xy with XY CPMG sequence
x1 x-filter in F1
x2 x-filter in F2
x3 x-filter in F3
zf with z-filter
zq zero quantum coherence
zs using a gradient/rf spoil pulse
1d 1D version
1s using 1 spoil gradients
11 using 1-1 pulse
19 using 3-9-19 pulse
19f for F19
2h using 2H lockswitch unit
2s using 2 spoil gradients
3d 3D sequence
3n for E.COSY (3 spins, negative correlation)
3p for E.COSY (3 spins, positive correlation)
3s using 3 spoil gradients
30 using a 30 degree flip angle
45 using a 45 degree flip angle
90 using a 90 degree flip angle
135 using a 135 degree flip angle
180 using a 180 degree pulse

;Typical experiment names would be:

; cosy, dept, dipsi2, hmbc, hmqc, hoesy, hsqc, inad, inept,
; mlev, noesy, roesy or trosy.

;Inverse correlations are denoted as hmbc, hmqc or hsqc.

; Experiments with a BIRD sequence in the beginning
; also contain a bi in the name.

;1D experiments, which are analogues of 2D experiments by virtue of
; a selective pulse, start with sel.
;Semiselective 2D experiments have the same name as the unselective
; version but with an s at the beginning:
;
; scosyph <-> cosyph.

;A phase-sensitive (States-TPPI, TPPI etc.) NOESY experiment with
; presaturation would then be:
;
; noesy + ph + pr = noesyphpr.

;In the other direction the pulseprogram hmbcgpndqf would be
;
; hmbc + gp + lp + nd + qf
;
; and therefor an:
;
; inverse correlation for long-range couplings (HMBC) with
; coherence selection using gradients with ":gp" syntax,
; low-pass J-filter,
; no decoupling
; in absolute value mode.

;The nomenclature of parameters is described in Pulprog.info.

;Comments like:

;
; ;avance-version
; ;begin ____
; ;end ____
;
; with (____ = MLEV17, DIPSI2, ...)
;

;are evaluated by NMRSIM for the pulseprogram display and should
;therefor not be removed. The syntax for begin/end statements allows
;characters, numbers and '_'. Arithmetic operators must not be used.

;
;
;The comments:

;
; ;preprocessor-flags-start
; ;preprocessor-flags-end
;

;are also evaluated to identify flags used in the pulseprogram and
;must also not be removed.

;\$Id: \$

```

;Param.info
;avance-version (10/02/01)
;
;The following convention is used for power levels, pulses, delays
;and loop counters throughout the microprograms:
;
;
;$CLASS=HighRes Info
;$COMMENT=

;pl0 :
;pl1 : f1 channel - power level for pulse (default)           {all, PL90(F1)}
;pl2 : f2 channel - power level for pulse (default)           {all, PL90(F2)}
;pl3 : f3 channel - power level for pulse (default)           {all, PL90(F3)}
;pl4 : f4 channel - power level for pulse (default)           {all, PL90(F4)}
;pl5 : f5 channel - power level for pulse (default)           {}
;pl6 : f6 channel - power level for pulse (default)           {}
;pl7 : f7 channel - power level for pulse (default)           {}
;pl8 : f8 channel - power level for pulse (default)           {}
;pl9 : f1 channel - power level for presaturation              {default+lcnmr+triple+triple2+triple_na,
SQPL[3](F1)}
;pl10: f1 channel - power level for TOCSY-spinlock             {all, SQPL[1](F1)}
;pl11: f1 channel - power level for ROESY-spinlock             {all, SQPL[2](F1)}
;pl12: f2 channel - power level for CPD/BB decoupling          {all, SQPL[0](F2)}
;pl13: f2 channel - power level for second CPD/BB decoupling  {default+lcnmr+triple_c, SQPL[4](F2)}
; or f2 channel - power level for Cbeta/CO decoupling         {triple+triple2, SHPL[22](F2)}
;pl14: f2 channel - power level for cw saturation              {default, SQPL[13](F2)}
; or f2 channel - power level for low power decoupling         {lcnmr+triple+triple2, SQPL[5](F2)}
;pl15: f2 channel - power level for TOCSY-spinlock             {all, SQPL[1](F2)}
;pl16: f3 channel - power level for CPD/BB decoupling          {all, SQPL[0](F3)}
;pl17: f4 channel - power level for CPD/BB decoupling          {all, SQPL[0](F4)}
;pl18: f1 channel - power level for 3-9-19-pulse (watergate)  {default+lcnmr+triple+triple2+triple_na,
PL90(F1)}
;pl19: f1 channel - power level for CPD/BB decoupling          {default+lcnmr+triple+triple2+triple_na,
SQPL[0](F1)}
;pl20: f1 channel - power level for Dante-z pulse              {}
; or f2 channel - power level for TOCSY-spinlock (high sel.)  {triple, SQPL[8](F2)}
; or f2 channel - power level for TOCSY-spinlock (med. sel.)  {triple_na, SQPL[7](F2)}
;pl21: f1 channel - power level for TOCSY-spinlock (med. sel.) {triple, SQPL[7](F1)}
; or f2 channel - power level for presaturation                {default+lcnmr, SQPL[3](F2)}
;pl22: f2 channel - power level for TOCSY-spinlock (med. sel.) {triple, SQPL[7](F2)}
; or f3 channel - power level for presaturation                {lcnmr, SQPL[3](F3)}
; or f3 channel - power level for TOCSY-spinlock (med. sel.)  {triple_na, SQPL[9](F3)}
;pl23: f3 channel - power level for TOCSY-spinlock             {default+lcnmr+triple_c, SQPL[1](F3)}
; or f3 channel - power level for Rexchange/T2 spinlock        {triple+triple2, SQPL[16](F3)}
; or f3 channel - power level for TOCSY-spinlock               {triple_na, SQPL[8](F3)}
;pl24: f2 channel - power level for hd/hc decoupling           {all, SQPL[14](F2)}
;pl25: f1 channel - power level for TOCSY spinlock (higher sel.) {triple_na, SQPL[8](F1)}
; or f2 channel - power level for Eretic                       {default+lcnmr, SHPL[20](F2)}
; or f3 channel - power level for T1rho spinlock               {triple+triple2, SQPL[15](F3)}
;pl26: f2 channel - power level for cw decoupling              {default+lcnmr, SQPL[12](F2)}
; or f2 channel - power level for TOCSY spinlock (higher sel. II) {triple_na, SQPL[8](F2)}
; or f3 channel - power level for low power decoupling         {triple, SQPL[5](F3)}
;pl27: f1 channel - power level for pulsed ROESY-spinlock      {default+lcnmr, SQPL[11](F1)}

```



```

; or f1 channel - power level for cleanex spinlock          {triple2, SQPL[10](F1)}
; or f2 channel - power level for TOCSY spinlock (higher sel. III) {triple_na, SQPL[9](F2)}
; or f3 channel - power level for TOCSY-spinlock          {triple, SQPL[1](F3)}
; pl28: f2 channel - power level for selective Ca or CO decoupling {triple+triple2, SHPL[20](F2)}
; or f2 channel - power level for selective decoupling      {triple_na, SHPL[28](F2)}
; pl29: f1 channel - power level for trim pulse (T1rho filter) in STD {default, SHPL[1](F1)}
; or f2 channel - power level for simultaneous Ca and CO decoupling {triple2, SHPL[21](F2)}
; pl30: f2 channel - power level for bilev decoupling      {default+triple+triple2+triple_na, SQPL[0](F2)}
; pl31: f2 channel - power level for bilev decoupling      {default+triple+triple2+triple_na, SQPL[6](F2)}
; pl32: f1 channel - power level for low power presaturation {default+lcnmr+triple+triple2+triple_na,
SQPL[12](F1)}
;
;
;
; sp0 : f1 channel - shaped pulse 180 degree (adiabatic TOCSY) {}
; or f2 channel - shaped pulse 180 degree (two-fold modulated) {triple_na, SH[29](F2)}
; sp1 : f1 channel - shaped pulse for selective excitation {default, SH[0](F1)}
; or f1 channel - shaped pulse for wet {lcnmr, SH[7](F1)}
; or f1 channel - shaped pulse for water flipback {triple+triple2+triple_na, SH[5](F1)}
; sp2 : f1 channel - shaped pulse 180 degree {default, SH[1](F1)}
; or f1 channel - shaped pulse for wet {lcnmr, SH[7](F1)}
; or f2 channel - shaped pulse 90 degree (on resonance) {triple+triple2, SH[6](F2)}
; or f2 channel - shaped pulse 90 degree (on resonance) {triple_na, SH[23](F2)}
; sp3 : f2 channel - shaped pulse 180 degree (adiabatic) {default+lcnmr, SH[4](F2)}
; or f2 channel - shaped pulse 180 degree (on resonance) {triple+triple2, SH[8](F2)}
; or f2 channel - shaped pulse 180 degree (on resonance) {triple_na, SH[25](F2)}
; sp4 : f2 channel - shaped pulse 90 degree (off resonance) {triple+triple2, SH[6](F2)}
; or f2 channel - shaped pulse 180 degree (short, broadband) {default+lcnmr, SH[14](F2)}
; sp5 : f1 channel - shaped pulse 180 degree (adiabatic) {default, SH[5](F1)}
; or f1 channel - shaped pulse 180 degree (off resonance) {triple+triple2, SH[8](F2)}
; or f2 channel - shaped pulse 180 degree (off resonance) {triple_na, SH[25](F2)}
; sp6 : f1 channel - shaped pulse for presaturation {default+lcnmr+triple2+triple_na, SH[4](F1)}
; or f2 channel - shaped pulse 90 degree (off res., time reversed) {triple, SH[7](F2)}
; sp7 : f1 channel - shaped pulse for wet {lcnmr, SH[7]*0.817(F1)}
; or f1 channel - shaped pulse 180 degree (adiabatic) {triple_c, SH[5](F1)}
; or f2 channel - shaped pulse 180 degree (adiabatic) {default, SH[5](F2)}
; or f2 channel - shaped pulse 180 degree (off resonance2) {triple+triple2, SH[8](F2)}
; or f2 channel - shaped pulse 180 degree (off resonance2) {triple_na, SH[25](F2)}
; sp8 : f1 channel - shaped pulse for wet {lcnmr, SHPL[7]*1.270(F1)}
; or f1 channel - shaped pulse 180 degree (adiabatic) {default, SH[4](F1)}
; or f2 channel - shaped pulse 90 degree (on res., time reversed) {triple+triple2, SH[7](F2)}
; or f2 channel - shaped pulse 90 degree (on res., time reversed) {triple_na, SH[24](F2)}
; sp9 : f1 channel - shaped pulse for wet {lcnmr, SHPL[7]*0.593(F1)}
; or f2 channel - shaped pulse 180 degree (higher selectivity) {triple+triple2, SH[11](F2)}
; or f3 channel - shaped pulse 180 degree (on resonance) {triple_na, SH[9](F3)}
; sp10: f1 channel - shaped pulse for tilted ROESY {}
; or f1 channel - shaped pulse for wet {lcnmr, SHPL[7]*3.198(F1)}
; or f1 channel - shaped pulse 180 degree (excitation sculpting) {default, SH[5](F1)}
; or f2 channel - shaped pulse 90 degree (higher selectivity) {triple+triple2, SH[9](F2)}
; or f2 channel - shaped pulse 90 degree (higher selectivity) {triple_na, SH[26](F2)}
; sp11: f1 channel - shaped pulse for water flipback {default+lcnmr, SH[5](F1)}
; or f1 channel - shaped pulse for water flipback2 {triple+triple2+triple_na, SH[6](F1)}
; or f2 channel - shaped pulse for water flipback {triple_c, SH[5](F2)}

```

```

;sp12: f1 channel - shaped pulse for wet2          {}
; or f2 channel - shaped pulse 90 degree (higher sel., time rev.) {triple+triple2, SH[10](F2)}
; or f2 channel - shaped pulse 90 degree (higher sel., time rev.) {triple_na, SH[26](F2)}
;sp13: f1 channel - shaped pulse for wet2          {}
; or f1 channel - shaped pulse 180 degree (adiabatic)          {triple_c, SH[4](F1)}
; or f2 channel - shaped pulse 180 degree (adiabatic)          {triple+triple2+triple_na, SH[4](F2)}
;sp14: f1 channel - shaped pulse for wet2          {}
; or f2 channel - shaped pulse 180 degree (adiabatic bilev decoupling) {default+triple+triple2, SH[13](F2)}
; or f3 channel - shaped pulse 180 degree (adiabatic)          {triple_na, SH[8](F3)}
;sp15: f2 channel - shaped pulse 180 degree for decoupling (Ca or CO) {triple, SH[20](F2)}
; or f2 channel - shaped pulse 180 degree for decoupling (Cbeta) {triple2, SH[22](F2)}
; or f2 channel - shaped pulse 180 degree for decoupling (C') {triple_na, SH[28](F2)}
;sp16: f2 channel - shaped pulse 180 degree (higher sel., off res.) {triple+triple2, SH[11](F2)}
;sp17: f2 channel - shaped pulse 180 degree (higher sel., off res.) {triple+triple2, SH[11](F2)}
;sp18: f2 channel - shaped pulse 180 degree (adiabatic matched sweep) {default+triple, SH[19](F2)}
;sp19: f1 channel - shaped pulse for wet          {default, SH[7]*0.817(F1)}
; or f2 channel - shaped pulse 180 degree (inversion (sharp)) {triple, SH[31](F2)}
; or f2 channel - shaped pulse 90 degree (NH)          {triple_c, SH[10](F2)}
; or f3 channel - shaped pulse 90 degree (T1rho, adiabatic ramp up) {triple2, SH[6](F3)}
;sp20: f1 channel - shaped pulse for wet          {default, SH[7]*1.270(F1)}
; or f2 channel - shaped pulse 180 degree (off resonance3) {triple, SH[8](F2)}
; or f2 channel - shaped pulse 90 degree (NH, time reversed) {triple_c, SH[11](F2)}
; or f3 channel - shaped pulse 90 degree (T1rho, adiabatic ramp down) {triple2, SH[7](F3)}
;sp21: f1 channel - shaped pulse for wet          {default, SH[7]*0.593(F1)}
; or f1 channel - shaped pulse 180 degree (cleanex, H2O) {triple2, SH[17](F1)}
; or f1 channel - shaped pulse 180 degree (med. selectivity) {triple_c, SH[15](F1)}
; or f2 channel - shaped pulse 180 degree (refocussing (sharp)) {triple, SH[32](F2)}
;sp22: f1 channel - shaped pulse for wet          {default, SH[7]*3.198(F1)}
; or f1 channel - shaped pulse 90 degree (cleanex, H2O) {triple2, SH[16](F1)}
; or f1 channel - shaped pulse 180 degree (off resonance) {triple_c, SH[8](F1)}
;sp23: f1 channel - shaped pulse 90 degree (on resonance) {triple_c, SH[6](F1)}
; or f1 channel - shaped pulse 120 degree (NH, best-) {triple, SH[8](F1)}
; or f1 channel - shaped pulse 180 degree (off resonance) {triple_na, SH[19](F1)}
; or f2 channel - shaped pulse 180 degree (med. selectivity) {triple2, SH[15](F2)}
;sp24: f1 channel - shaped pulse 180 degree (on resonance) {triple_c, SH[8](F1)}
; or f1 channel - shaped pulse 180 degree (NH, best-, I) {triple, SH[9](F1)}
; or f1 channel - shaped pulse 180 degree (off resonance2) {triple_na, SH[19](F1)}
; or f2 channel - shaped pulse 180 degree (high selectivity) {triple2, SH[16](F2)}
;sp25: f1 channel - shaped pulse 90 degree (on res., time reversed) {triple_c, SH[7](F1)}
; or f1 channel - shaped pulse 90 degree (NH, best-, I) {triple, SH[10](F1)}
; or f2 channel - shaped pulse 180 degree (higher selectivity) {triple_na, SH[28](F2)}
; or f2 channel - shaped pulse 90 degree (high selectivity) {triple2, SH[17](F2)}
;sp26: f1 channel - shaped pulse 180 degree (off resonance) {triple_c, SH[8](F1)}
; or f1 channel - shaped pulse 180 degree (NH, best-, II) {triple, SH[12](F1)}
; or f1 channel - shaped pulse 180 degree (C, selective) {default, SH[1](F1)}
; or f2 channel - shaped pulse 90 degree (high selectivity, tr) {triple2, SH[18](F2)}
;sp27: f1 channel - shaped pulse 180 degree (off resonance) {triple_c, SH[8](F1)}
; or f1 channel - shaped pulse 90 degree (NH, best-, I tr) {triple, SH[11](F1)}
; or f2 channel - shaped pulse 90 degree (high selectivity) {triple2, SH[17](F2)}
;sp28: f1 channel - shaped pulse 180 degree (higher selectivity) {triple_c, SH[11](F1)}
; or f1 channel - shaped pulse 90 degree (NH, best-, II) {triple, SH[13](F1)}
; or f2 channel - shaped pulse 180 degree (higher selectivity) {triple2, SH[11](F2)}
;sp29: f1 channel - shaped pulse 180 degree (off resonance) {triple_c, SH[8](F1)}
; or f1 channel - shaped pulse 180 degree (adiabatic sweep: z-spoil) {default, SH[18](F1)}

```

```

; or f1 channel - shaped pulse 90 degree (NH, best-, II tr)      {triple, SH[14](F1)}
; or f2 channel - shaped pulse 180 degree (high selectivity)    {triple2, SH[16](F2)}
;sp30: f1 channel - shaped pulse 180 degree (sim. Ca + CO)      {triple_c, SH[21](F1)}
; or f1 channel - shaped pulse 180 degree (broadband, best-)    {triple, SH[15](F1)}
; or f2 channel - shaped pulse 180 degree for decoupling (sim. Ca + CO) {triple2, SH[21](F2)}
;sp31: f2 channel - shaped pulse 180 degree (adiabatic bilev decoupling) {default+triple+triple2+triple_na,
SH[12](F2)}
;
;
;
;
;p0 :                                                              {all, PW90(F1)}
;p1 : f1 channel - 90 degree high power pulse                    {all, PW90(F1)}
;p2 : f1 channel - 180 degree high power pulse                   {all, PW90*2(F1)}
;p3 : f2 channel - 90 degree high power pulse                    {all, PW90(F2)}
;p4 : f2 channel - 180 degree high power pulse                   {all, PW90*2(F2)}
;p5 : f1 channel - 60 degree low power pulse                     {all, SQPW[1]*0.66(F1)}
;p6 : f1 channel - 90 degree low power pulse                     {all, SQPW[1](F1)}
;p7 : f1 channel - 180 degree low power pulse                     {default+lcnmr+triple+triple_c, SQPW[1]*2(F1)}
; or f1 channel - 180 degree shaped pulse (cleanex sel. H2O)    {triple2, SHPW[17](F1)}
; or f2 channel - 90 degree pulse at pl20 (TOCSY, higher sel.) {triple_na, SQPW[7](F2)}
;p8 : f2 channel - 60 degree low power pulse                      {}
; or f1 channel - 90 degree shaped pulse (wet)                   {default, SHPW[7](F1)}
; or f1 channel - 180 degree shaped pulse (adiabatic)           {triple_c, SHPW[4](F1)}
; or f2 channel - 180 degree shaped pulse (adiabatic)           {triple+triple2+triple_na, SHPW[4](F2)}
;p9 : f2 channel - 90 degree low power pulse (TOCSY)            {all, SQPW[1](F2)}
;p10: f1 channel - 90 degree low power pulse (cleanex spinlock) {triple2, SQPW[10](F1)}
; or f2 channel - 180 degree low power pulse                     {default+lcnmr+triple+triple_c, SQPW[1]*2(F2)}
; or f2 channel - 180 degree shaped pulse (higher selectivity)  {triple_na, SHPW[28](F2)}
;p11: f1 channel - 90 degree shaped pulse (selective excitation) {default, SHPW[0](F1)}
; or f1 channel - 90 degree shaped pulse (selective excitation) {triple_c, SHPW[6](F1)}
; or f1 channel - 90 degree shaped pulse (wet)                   {lcnmr, SHPW[7](F1)}
; or f1 channel - 90 degree shaped pulse (water flipback/watergate) {triple+triple2+triple_na, SHPW[5](F1)}
;p12: f1 channel - 180 degree shaped pulse (H, selective)        {default+lcnmr, SHPW[1](F1)}
; or f1 channel - 180 degree shaped pulse (C, selective)        {triple_c, SHPW[8](F1)}
; or f1 channel - 180 degree shaped pulse (excitation sculpting) {triple+triple2, SHPW[5]*2(F1)}
; or f1 channel - 180 degree shaped pulse (H, selective)        {triple_na, SHPW[19](F1)}
;p13: f1 channel - 180 degree shaped pulse (C, adiabatic)        {default+lcnmr, SHPW[5](F1)}
; or f2 channel - 90 degree shaped pulse                         {triple+triple2, SHPW[6](F2)}
; or f2 channel - 90 degree shaped pulse                         {triple_na, SHPW[23](F2)}
; or f2 channel - 90 degree shaped pulse (H, selective)         {triple_c, SHPW[10](F2)}
;p14: f2 channel - 180 degree shaped pulse (adiabatic)           {default+lcnmr, SHPW[4](F2)}
; or f2 channel - 180 degree shaped pulse (selective)           {triple+triple2, SHPW[8](F2)}
; or f2 channel - 180 degree shaped pulse (selective)           {triple_na, SHPW[25](F2)}
;p15: f1 channel - pulse for ROESY spinlock                      {default+lcnmr, SQPW[2](F1)}
; or f1 channel - 90 degree shaped pulse (cleanex sel. H2O)    {triple2, SHPW[16](F1)}
; or f2 channel - 180 degree shaped pulse (adiabatic matched sweep) {triple, SHPW[19](F2)}
; or f2 channel - 90 degree shaped pulse (higher selectivity)   {triple_na, SHPW[26](F2)}
;p16: homospoil/gradient pulse                                   {all, P_grad1}
;p17: f1 channel - trim pulse at pl10 or pl15                    {all, P_mlev(F1)}
;p18: f1 channel - shaped pulse (off resonance presaturation)    {default+lcnmr+triple+triple2+triple_na,
SHPW[4](F1)}
;p19: homospoil/gradient pulse 2                                 {all, P_grad2}
;p20: f2 channel - trim pulse                                    {all, P_mlev(F2)}

```

```

;p21: f3 channel - 90 degree high power pulse      {all, PW90(F3)}
;p22: f3 channel - 180 degree high power pulse     {all, PW90*2(F3)}
;p23: f1 channel - 180 degree shaped pulse (med. selectivity) {triple_c, SHPW[15](F1)}
; or f2 channel - 90 degree shaped pulse (higher selectivity) {triple+triple2, SHPW[9](F2)}
; or f2 channel - 90 degree shaped pulse (twofold modulated) {triple_na, SHPW[29](F2)}
; or f4 channel - 90 degree high power pulse      {default, PW90(F4)}
;p24: f1 channel - 180 degree shaped pulse (adiabatic) {triple_c, SHPW[5](F1)}
; or f2 channel - 180 degree shaped pulse (adiabatic) {default, SHPW[5](F2)}
; or f2 channel - 180 degree shaped pulse (higher selectivity) {triple+triple2, SHPW[11](F2)}
; or f3 channel - 90 degree pulse at pl22 (TOCSY, higher sel.) {triple_na, SQPW[9](F3)}
; or f4 channel - 180 degree high power pulse      {}
;p25: f1 channel - 90 degree pulse at pl27 (pulsed ROESY) {default+lcnmr, SQPW[11]*2(F1)}
; or f1 channel - 90 degree shaped pulse (higher selectivity) {triple_c, SHPW[11](F1)}
; or f3 channel - pulse for t1rho experiment        {pp}
; or f3 channel - pulse for TOCSY-spinlock experiment {triple, SQPW[1]}
; or f3 channel - 90 degree pulse at pl23 (TOCSY) {triple_na, SQPW[8](F3)}
; or f3 channel - 180 degree low power pulse (Rexchange) {triple2, SQPW[16]*2(F3)}
;p26: f1 channel - 90 degree pulse at pl19          {triple*, SQPW[0](F1)}
; or f1 channel - 180 degree shaped pulse (adiabatic) {default, SHPW[4](F1)}
;p27: f1 channel - 90 degree pulse at pl18 (3-9-19 watergate) {default+lcnmr+triple+triple2+triple_na, PW90(F1)}
;p28: f1 channel - trim pulse at pl1                {all, P_hsqc(F1)}
;p29: f1 channel - 90 degree shaped pulse (water flipback) {default, SHPW[5](F1)}
; or f1 channel - 90 degree shaped pulse (water flipback2) {triple+triple_na, SHPW[6](F1)}
; or f2 channel - 90 degree shaped pulse (water flipback) {triple_c, SHPW[5](F2)}
; or f3 channel - 90 degree shaped pulse (T1rho adiabatic ramp) {triple2, SHPW[6](F3)}
; or homospoil/gradient pulse 3                    {pp}
;p30: f1 channel - 180 degree shaped pulse (sim. Ca + CO) {triple_c, SHPW[21](F1)}
; or f2 channel - 180 degree shaped pulse (sim. Ca + CO decoupling) {triple2, SHPW[21](F2)}
; or f3 channel - 180 degree pulse at pl23 (T2)      {triple, SQPW[16]*2(F3)}
; or f3 channel - 180 degree shaped pulse          {triple_na, SHPW[5](F3)}
; or homospoil/gradient pulse 4                    {pp}
; or gradient pulse for diffusion (dosy)            {}
;p31: f2 channel - 180 degree shaped pulse (adiabatic matched sweep) {default, SHPW[19](F2)}
; or f2 channel - 180 degree shaped pulse (sel. Ca or CO decoupling) {triple, SHPW[20](F2)}
; or f2 channel - 180 degree shaped pulse (Cbeta decoupling) {triple2, SHPW[22](F2)}
; or f2 channel - 180 degree shaped pulse (sel. C decoupling) {triple_na, SHPW[28](F2)}
; or f2 channel - 90 degree pulse (low power decoupling) {lcnmr, SQPW[5](F2)}
; or homospoil/gradient pulse 5                    {pp}
;p32: f1 channel - 180 degree shaped pulse (adiabatic sweep: z-spoil) {default, SHPW[18](F1)}
; or f3 channel - 180 degree shaped pulse (adiabatic) {triple_na, SHPW[8](F3)}
;p33: f2 channel - 180 degree shaped pulse (med. selectivity) {triple2, SHPW[15](F2)}
; or f3 channel - trim pulse                        {triple_na, P_mlev(F3)}
;p34: f2 channel - 180 degree shaped pulse (high selectivity) {triple2, SHPW[16](F2)}
;p35: f2 channel - 90 degree shaped pulse (high selectivity) {triple2, SHPW[17](F2)}
;p36: f2 channel - 180 degree shaped pulse (C, selective) {default, SHPW[1](F2)}
;
;p39: f1 channel - 120 degree shaped pulse for excitation (best-) {triple, SHPW[8](F1)}
; or f2 channel - 180 degree shaped pulse (short, broadband) {default, SHPW[14](F2)}
;p40: f1 channel - 180 degree shaped pulse for refocussing (best-) {triple, SHPW[9](F1)}
; or f1 channel - 180 degree shaped pulse (excitation sculpting) {default, SHPW[5]*2(F1)}
;p41: f1 channel - 90 degree shaped pulse for refocussing (best-) {triple, SHPW[10](F1)}
;p42: f1 channel - 180 degree shaped pulse for refocussing (best-) {triple, SHPW[12](F1)}
;p43: f1 channel - 90 degree shaped pulse for refocussing (best-) {triple, SHPW[13](F1)}

```

```

;p44: f1 channel - 180 degree shaped pulse for refocussing (best-) {triple, SHPW[15](F1)}
;p45: f2 channel - 90 degree pulse at pl20 ((hetero) TOCSY high sel.) {triple, SQPW[8](F2)}
;p46: f2 channel - 90 degree pulse at pl20 ((hetero) TOCSY med. sel.) {triple, SQPW[7](F2)}
;p47: f2 channel - 180 degree shaped pulse for inversion (sharp) {triple, SHPW[31](F2)}
;p48: f2 channel - 180 degree shaped pulse for refocussing (sharp) {triple, SHPW[32](F2)}
;p49: f1 channel - 180 degree shaped pulse (H2O) {triple, SHPW[21](F1)}
;
;p61: f2 channel - 90 degree pulse (low power decoupling) {triple+triple2, SQPW[5](F2)}
;p62: f3 channel - 90 degree pulse (low power decoupling) {triple+triple2, SQPW[5](F3)}
;p63: f2 channel - 180 degree shaped pulse (adiabatic biley sweep) {default+triple+triple2+triple_na,
SHPW[12](F2)}
;
;
;
;d0 : incremented delay (2D or 3D) [3 usec]
;d1 : relaxation delay; 1-5 * T1
;d2 : 1/(2J)
;d3 : 1/(3J) or 1/(6J)
;d4 : 1/(4J)
;d5 : DE/2
;d6 : delay for evolution of long range couplings
;d7 : delay for inversion recovery
;d8 : NOESY mixing time
;d9 : TOCSY mixing time {all, TTOC(F1)}
;d10: incremented delay (3D)
;d11: delay for disk I/O [30 msec]
;d12: delay for power switching [20 usec]
;d13: short delay [4 usec]
;d14: delay for evolution after shaped pulse
;d15: TOCSY mixing time (CC) {triple*, TTOC(F2)}
;d16: delay for homospoil/gradient recovery {all, D_grad}
;d17: delay for DANTE pulse-train
;d18: delay for evolution of long range couplings
;d19: delay for binomial water suppression
;d20: for different applications
;d21: for different applications
;d22: 1/(2J(XY))
;d23: 1/(4J(XY)) or 1/(2J(XY))
;d24: for different applications
;d25: 1/(6J(YH)) or 1/(8J(XY))
;d26: 1/(4J(YH))
;d27: for different applications
;d28: for different applications
;d29: for different applications
;d30: for different applications
;d31: incremented delay (> 3D, t1)
;d32: incremented delay (> 3D, t2)
;d33: incremented delay (> 3D, t3)
;d34: incremented delay (> 3D, t4)
;d35: incremented delay (> 3D, t5)
;d41: decremented delay (> 3D, t1)
;d42: decremented delay (> 3D, t2)
;d43: decremented delay (> 3D, t3)

```

```

;d44: decremented delay (> 3D, t4)
;d45: decremented delay (> 3D, t5)
;d51: incremented delay (> 3D, t1)
;d52: incremented delay (> 3D, t2)
;d53: incremented delay (> 3D, t3)
;d54: incremented delay (> 3D, t4)
;d55: incremented delay (> 3D, t5)
;
;
;
;
;cnst0 : for protein experiments - N chemical shift (offset, in ppm)
; or   for na experiments - calculated chemical shift (offset, in ppm)
; or   for na experiments - N(aromatic) chemical shift (offset, in ppm)    [195 ppm]
;cnst1 : J (HH)
;cnst2 : J (XH)
;cnst3 : J (XX)
;cnst4 : J (YH)
;cnst5 : J (XY)
;cnst6 : J (XH)min
;cnst7 : J (XH)max
;cnst8 : bandwidth of excitation for Dante-z pulse
;cnst9 : for different applications as J
;cnst10: for different applications as J
;cnst11: for multiplicity selection
;cnst12: for multiplicity selection
;cnst13: J (XH) long range
;cnst14: J (XH) long range (min)
;cnst15: J (XH) long range (max)
;cnst16: J-scale factor
; or   for na experiments - H6/8 and/or H1' chemical shift (offset, in ppm)
;cnst17: factor to compensate for coupling evolution during a pulse
; or   for na experiments - H1' chemical shift (offset, in ppm)
;cnst18: for protein experiments - H2O chemical shift (offset, in ppm)
; or   for na experiments - H2O chemical shift (offset, in ppm)
;cnst19: for protein experiments - H(N) chemical shift (offset, in ppm)
; or   : for na experiments - H(N) chemical shift (offset, in ppm)
;cnst20: for protein experiments - Haliphatic chemical shift (offset, in ppm)
;cnst21: for na experiments - C1' chemical shift (offset, in ppm)    [90 ppm]
; or   for protein experiments - CO chemical shift (offset, in ppm)
;cnst22: for protein experiments - Calpha chemical shift (offset, in ppm)
; or   for na experiments - C6/8 chemical shift (offset, in ppm)    [137 ppm]
;cnst23: for protein experiments - Caliphatic chemical shift (offset, in ppm)
; or   for na experiments - C2' chemical shift (offset, in ppm)    [72 ppm]
;cnst24: for protein experiments - Caromatic chemical shift (offset, in ppm)
; or   for na experiments - C4 (C/U) chemical shift (offset, in ppm)  [169 ppm]
;cnst25: for protein experiments - flag for cross peak / reference experiments
; or   for na experiments - C6 (A) chemical shift (offset, in ppm)  [160 ppm]
;cnst26: for protein experiments - Call chemical shift (offset, in ppm)
; or   for na experiments - C5 (G) chemical shift (offset, in ppm)  [119 ppm]
;cnst27: for protein experiments - ( Cgamma chemical shift (offset, in ppm) )
; or   for na experiments - C2/4 chemical shift (offset, in ppm)  [152 ppm]
;cnst28: for protein experiments - Haromatic chemical shift (offset, in ppm)
; or   for na experiments - C5 (C/U) chemical shift (offset, in ppm)  [105 ppm]

```

```

;cnst29: for protein experiments - N(H) chemical shift (offset, in ppm)
; or   for na experiments - C(aro) chemical shift (offset, in ppm) [145 ppm]
;cnst30: for protein experiments - Cbeta chemical shift (offset, in ppm)
; or   for na experiments - N(H) chemical shift (offset, in ppm) [151 ppm]
;cnst31: scaling factor
; or   for na experiments - N(H2) chemical shift (offset, in ppm) [81 ppm]
;cnst34: for music experiments - Cab(Leu) chemical shift (offset, in ppm) [48.0 ppm]
;cnst35: for music experiments - Cbgd(Leu) chemical shift (offset, in ppm) [35.0 ppm]
;cnst36: for music experiments - Cgd(Leu) chemical shift (offset, in ppm) [19.0 ppm]
;cnst37: for music experiments - Cb(Val)/Cg(Ile) chem. shift (offset,in ppm) [27.0 ppm]
;cnst38: for music experiments - Cab(Asp)Cd(Arg) chem. shift (offset,in ppm) [42.0 ppm]
;cnst39: for music experiments - Cg(Arg)/Cd(Lys) chem. shift (offset,in ppm) [29 ppm]
;cnst40: compensation of chemical shift evolution during p42 (best-)
; or   for music experiments - Ce(Lys) chemical shift (offset, in ppm) [44.5 ppm]
;cnst41: compensation of chemical shift evolution during p41 (best-)
; or   for music experiments - Cg(Phe,Tyr,His) chem. shift (offset, in ppm)[136.0 ppm]
;cnst42: compensation of chemical shift evolution during p42 (best-)
; or   for music experiments - Cg(Trp) chemical shift (offset, in ppm) [104.0 ppm]
;cnst43: compensation of chemical shift evolution during p43 (best-)
; or   for music experiments - Cb(Thr) chemical shift (offset, in ppm) [68.5 ppm]
;cnst44: for music experiments - Cab(Ser) chemical shift (offset, in ppm) [60.5 ppm]
;cnst45: for music experiments - Cd(Ile) chemical shift (offset, in ppm) [11.0 ppm]
;cnst46: for music experiments - Cbg(Ile) chemical shift (offset, in ppm) [33.0 ppm]
;cnst47: for music experiments - N(H) chemical shift (offset, in ppm) [117.0 ppm]
;cnst48: for music experiments - N(Pro) chemical shift (offset, in ppm) [134.0 ppm]
;cnst49: for protein experiments - H(methyl) chemical shift (offset, in ppm) [0.5 ppm]
;
;
;
;
;l1 : loop for spinlock cycle
;l2 : loop for GARP cycle: l2 * 31.75 * 4 * p9 => AQ
;l3 : loop for phase sensitive 2D or 3D using
;   States et al. or States-TPPI method: l3 = td1/2
;l4 : for different applications
;l5 : for different applications
;l6 : loop for shaped pulse presaturation during relaxation delay
;l7 : loop for shaped pulse presaturation during mixing time
;l8 : number of frequencies for multiple presaturation
;l11: loop for spinlock cycle 2
;
;
;
;
;vc : variable loop counter, taken from vc-list
;vd : variable delay, taken from vd-list

;$Id: $

```

```

;Relations.info
;avance-version (10/02/01)
;
;$CLASS=HighRes Info
;$COMMENT=
;The following convention is used for power levels, pulses, delays
;and loop counters in the different relation files for prosol:
;
;all = default + lcnmr + triple + triple2 + triple_c + triple_na
;triple* = triple + triple2 + triple_c + triple_na
;!__ = except
;
;prosol par.    rel. file        pulseprogram parameter
;
;DE            all            de
;D_grad        all            d16: delay for homospoil/gradient recovery
;
;PW90(F1)     all            p0 :
;PW90(F1)     all            p1 : f1 channel - 90 degree high power pulse
;PW90(F1)     all(!triple_c) p27: f1 channel - 90 degree pulse at pl18 (3-9-19 watergate)
;PW90*2(F1)   all            p2 : f1 channel - 180 degree high power pulse
;PW90(F2)     all            p3 : f2 channel - 90 degree high power pulse
;PW90*2(F2)   all            p4 : f2 channel - 180 degree high power pulse
;PW90(F3)     all            p21: f3 channel - 90 degree high power pulse
;PW90*2(F3)   all            p22: f3 channel - 180 degree high power pulse
;PW90(F4)     default        p23: f4 channel - 90 degree high power pulse
;
;PL90(F1)     all            pl1 : f1 channel - power level for pulse (default)
;PL90(F1)     all(!triple_c) pl18: f1 channel - power level for 3-9-19-pulse (watergate)
;PL90(F2)     all            pl2 : f2 channel - power level for pulse (default)
;PL90(F3)     all            pl3 : f3 channel - power level for pulse (default)
;PL90(F4)     all            pl4 : f4 channel - power level for pulse (default)
;
;SH[0](F1)   default        sp1 : f1 channel - shaped pulse for selective excitation
;SH[1](F1)   default        sp2 : f1 channel - shaped pulse 180 degree
;SH[1](F2)   default        sp26: f2 channel - shaped pulse 180 degree
;SH[4](F1)   all(!triple+triple_c) sp6 : f1 channel - shaped pulse for presaturation
;SH[4](F1)   default        sp8 : f1 channel - shaped pulse 180 degree (adiabatic)
;SH[4](F1)   triple_c       sp13: f1 channel - shaped pulse 180 degree (adiabatic)
;SH[4](F2)   default+lcnmr  sp3 : f2 channel - shaped pulse 180 degree (adiabatic)
;SH[4](F2)   triple+triple2+triple_na sp13: f2 channel - shaped pulse 180 degree (adiabatic)
;SH[5](F1)   triple+triple2+triple_na sp1 : f1 channel - shaped pulse for water flipback
;SH[5](F1)   default        sp5 : f1 channel - shaped pulse 180 degree (adiabatic)
;SH[5](F1)   triple_c       sp7 : f1 channel - shaped pulse 180 degree (adiabatic)
;SH[5](F1)   default        sp10: f1 channel - shaped pulse for water flipback
;SH[5](F1)   default+lcnmr  sp11: f1 channel - shaped pulse for water flipback
;SH[5](F2)   default        sp7 : f2 channel - shaped pulse 180 degree (adiabatic)
;SH[5](F2)   triple_c       sp11: f2 channel - shaped pulse 90 degree (water flipback)
;SH[6](F1)   triple+triple2+triple_na sp11: f1 channel - shaped pulse for water flipback2
;SH[6](F1)   triple_c       sp23: f1 channel - shaped pulse 90 degree (on resonance)
;SH[6](F2)   triple+triple2  sp2 : f2 channel - shaped pulse 90 degree (on resonance)
;SH[6](F2)   triple+triple2  sp4 : f2 channel - shaped pulse 90 degree (off resonance)
;SH[6](F3)   triple2        sp19: f3 channel - shaped pulse 90 degree (T1rho, adiab. ramp up)

```


:SH[7](F1) lcnmr	sp1 : f1 channel - shaped pulse for wet
:SH[7](F1) lcnmr	sp2 : f1 channel - shaped pulse for wet
:SH[7]*0.817(F1) lcnmr	sp7 : f1 channel - shaped pulse for wet
:SH[7]*0.817(F1) default	sp19 : f1 channel - shaped pulse for wet
:SH[7]*1.270(F1) default	sp20: f1 channel - shaped pulse for wet
:SH[7]*0.593(F1) default	sp21: f1 channel - shaped pulse for wet
:SH[7]*3.198(F1) default	sp22: f1 channel - shaped pulse for wet
:SH[7](F1) triple_c	sp25: f1 channel - shaped pulse 90 degree (on resonance)
:SH[7](F2) triple	sp6 : f2 channel - shaped pulse 90 degree (off res., time reversed)
:SH[7](F2) triple+triple2	sp8 : f2 channel - shaped pulse 90 degree (on res., time reversed)
:SH[7](F3) triple2	sp20: f3 channel - shaped pulse 90 degree (T1rho, diab. ramp down)
:SH[8](F1) triple_c	sp22: f1 channel - shaped pulse 180 degree (off resonance)
:SH[8](F1) triple	sp23: f1 channel - shaped pulse 120 degree (NH, best-)
:SH[8](F1) triple_c	sp24: f1 channel - shaped pulse 180 degree (on resonance)
:SH[8](F1) triple_c	sp26: f1 channel - shaped pulse 180 degree (off resonance)
:SH[8](F1) triple_c	sp27: f1 channel - shaped pulse 180 degree (off resonance)
:SH[8](F1) triple_c	sp29: f1 channel - shaped pulse 180 degree (off resonance)
:SH[8](F2) triple+triple2	sp3 : f2 channel - shaped pulse 180 degree (on resonance)
:SH[8](F2) triple+triple2	sp5 : f2 channel - shaped pulse 180 degree (off resonance)
:SH[8](F2) triple+triple2	sp7 : f2 channel - shaped pulse 180 degree (off resonance2)
:SH[8](F2) triple	sp20: f2 channel - shaped pulse 180 degree (off resonance3)
:SH[8](F3) triple_na	sp14: f3 channel - shaped pulse 180 degree (_NA: N, adiabatic)
:SH[9](F1) triple	sp24: f1 channel - shaped pulse 180 degree (NH, best-, I)
:SH[9](F2) triple+triple2	sp10: f2 channel - shaped pulse 90 degree (higher selectivity)
:SH[9](F3) triple_na	sp9 : f3 channel - shaped pulse 180 degree (_NA: N)
:SH[10](F2) triple+triple2	sp12: f2 channel - shaped pulse 90 degree (higher sel., time rev.)
:SH[10](F2) triple_c	sp19: f2 channel - shaped pulse 90 degree (HN)
:SH[10](F1) triple	sp25: f1 channel - shaped pulse 90 degree (NH, best-, I)
:SH[11](F1) triple	sp27: f1 channel - shaped pulse 90 degree (NH, best-, I tr)
:SH[11](F1) triple_c	sp28: f1 channel - shaped pulse 180 degree (higher selectivity)
:SH[11](F2) triple+triple2	sp9 : f2 channel - shaped pulse 180 degree (higher selectivity)
:SH[11](F2) triple+triple2	sp16: f2 channel - shaped pulse 180 degree (higher selectivity)
:SH[11](F2) triple+triple2	sp17: f2 channel - shaped pulse 180 degree (higher selectivity)
:SH[11](F2) triple_c	sp20: f2 channel - shaped pulse 90 degree (HN tr)
:SH[11](F2) triple2	sp28: f2 channel - shaped pulse 180 degree (higher selectivity)
:SH[12](F1) triple	sp26: f1 channel - shaped pulse 180 degree (NH, best-, II)
:SH[12](F2) default+triple+triple2	sp31: f2 channel - shaped pulse 180 degree (adiabatic decoupling)
: +triple_na	
:SH[13](F1) triple	sp28: f1 channel - shaped pulse 90 degree (NH, best-, II)
:SH[13](F2) default+triple+triple2	sp14: f2 channel - shaped pulse 180 degree (adiabatic bilev decoupling)
:SH[14](F1) triple	sp29: f1 channel - shaped pulse 90 degree (NH, best-, II tr)
:SH[14](F2) default+lcnmr	sp4 : f2 channel - shaped pulse 180 degree (short, broadband)
:SH[15](F1) triple	sp30: f1 channel - shaped pulse 180 degree (broadband, best-)
:SH[15](F1) triple_c	sp21: f1 channel - shaped pulse 180 degree (med. selectivity)
:SH[15](F2) triple2	sp23: f2 channel - shaped pulse 180 degree (med. selectivity)
:SH[16](F1) triple2	sp22: f1 channel - shaped pulse 90 degree (H2O on resonance)
:SH[16](F2) triple2	sp24: f2 channel - shaped pulse 180 degree (high selectivity)
:SH[16](F2) triple2	sp29: f2 channel - shaped pulse 180 degree (high selectivity)
:SH[17](F1) triple2	sp21: f1 channel - shaped pulse 180 degree (H2O on resonance)
:SH[17](F2) triple2	sp25: f2 channel - shaped pulse 90 degree (high selectivity)
:SH[17](F2) triple2	sp27: f2 channel - shaped pulse 90 degree (high selectivity)
:SH[18](F1) default	sp29: f1 channel - shaped pulse 180 degree (adiabatic: z-spoil)
:SH[18](F2) triple2	sp26: f2 channel - shaped pulse 90 degree (high selectivity, tr)
:SH[19](F1) triple_na	sp23: f1 channel - shaped pulse 180 degree (_NA: H)
:SH[19](F1) triple_na	sp24: f1 channel - shaped pulse 180 degree (_NA: H)

;SH[19](F2)	default+triple	sp18: f2 channel - shaped pulse 180 degree (adiabatic matched sweep)
;SH[20](F2)	triple	sp15: f2 channel - shaped pulse 180 degree for decoupling (Ca or CO)
;SH[21](F1)	triple_c	sp30: f1 channel - shaped pulse 180 degree (sim. Ca + CO decoupling)
;SH[21](F1)	triple	sp32: f1 channel - shaped pulse 180 degree (H2O)
;SH[21](F2)	triple2	sp30: f2 channel - power level for simultaneous Ca and CO decoupling
;SH[22](F2)	triple2	sp15: f2 channel - shaped pulse 180 degree for decoupling (Cbeta)
;SH[23](F2)	triple_na	sp2 : f2 channel - shaped pulse 90 degree (_NA: C)
;SH[24](F2)	triple_na	sp8 : f2 channel - shaped pulse 90 degree (_NA: C, tr)
;SH[25](F2)	triple_na	sp3 : f2 channel - shaped pulse 180 degree (_NA: C)
;SH[25](F2)	triple_na	sp5 : f2 channel - shaped pulse 180 degree (_NA: C)
;SH[25](F2)	triple_na	sp7 : f2 channel - shaped pulse 180 degree (_NA: C)
;SH[26](F2)	triple_na	sp10: f2 channel - shaped pulse 90 degree (_NA: C, higher sel.)
;SH[27](F2)	triple_na	sp12: f2 channel - shaped pulse 90 degree (_NA: C, higher sel., tr)
;SH[28](F2)	triple_na	sp15: f2 channel - shaped pulse 180 degree (_NA: C, decoupling)
;SH[28](F2)	triple_na	sp25: f2 channel - shaped pulse 180 degree (_NA: C, higher sel.)
;SH[29](F2)	triple_na	sp0 : f2 channel - shaped pulse 180 degree (_NA: C, twofold mod)
;SH[31](F2)	triple	sp19: f2 channel - shaped pulse 180 degree (inversion (sharp))
;SH[32](F2)	triple	sp21: f2 channel - shaped pulse 180 degree (refocussing (sharp))
:		
;SHPL[7]*1.270(F1) lcnmr		sp8 : f1 channel - shaped pulse for wet
;SHPL[7]*0.593(F1) lcnmr		sp9 : f1 channel - shaped pulse for wet
;SHPL[7]*3.198(F1) lcnmr		sp10: f1 channel - shaped pulse for wet
;SHPL[20](F2) default+lcnmr		pl25: f2 channel - shaped pulse (Eretic)
;SHPL[20](F2) triple+triple2		pl28: f2 channel - power level for selective Ca or CO decoupling
;SHPL[21](F2) triple2		pl29: f2 channel - power level for simultaneous Ca and CO decoupling
;SHPL[22](F2) triple+triple2		pl13: f2 channel - power level for Cbeta/CO decoupling
;SHPL[28](F2) triple_na		pl28: f2 channel - shaped pulse 180 degree (_NA: C, decoupling)
:		
;SHPW[0](F1) default		p11: f1 channel - 90 degree shaped pulse (selective excitation)
;SHPW[1](F1) default+lcnmr		p12: f1 channel - 180 degree shaped pulse (H, selective)
;SHPW[1](F2) default		p36: f2 channel - 180 degree shaped pulse (C, selective)
;SHPW[4](F1) triple_c		p8 : f1 channel - 180 degree shaped pulse (adiabatic)
;SHPW[4](F1) all(!triple_c)		p18: f1 channel - shaped pulse (off resonance presaturation)
;SHPW[4](F2) triple+triple2+triple_na		p8 : f2 channel - 180 degree shaped pulse (adiabatic)
;SHPW[4](F2) default+lcnmr		p14: f2 channel - 180 degree shaped pulse (adiabatic)
;SHPW[4](F1) default		p26: f1 channel - 180 degree shaped pulse (adiabatic)
;SHPW[5](F1) triple+triple2+triple_na		p11: f1 channel - 90 degree shaped pulse (water flipback/watergate)
;SHPW[5]*2(F1) triple+triple2		p12: f1 channel - 180 degree shaped pulse (excitation sculpting)
;SHPW[5](F1) default+lcnmr		p13: f1 channel - 180 degree shaped pulse (C, adiabatic)
;SHPW[5](F1) triple_c		p24: f1 channel - 180 degree shaped pulse (adiabatic)
;SHPW[5](F1) default		p29: f1 channel - 90 degree shaped pulse (water flipback)
;SHPW[5]*2(F1) default		p40: f1 channel - 180 degree shaped pulse (excitation sculpting)
;SHPW[5](F2) default+lcnmr		p24: f2 channel - 180 degree shaped pulse (adiabatic)
;SHPW[5](F2) triple_c		p29: f2 channel - 90 degree shaped pulse for inversion (water flipback)
;SHPW[5](F3) triple_na		p30: f3 channel - 180 degree shaped pulse (_NA: N)
;SHPW[6](F1) triple_c		p11: f1 channel - 90 degree shaped pulse
;SHPW[6](F1) triple+triple_na		p29: f1 channel - 90 degree shaped pulse (water flipback2)
;SHPW[6](F2) triple+triple2		p13: f2 channel - 90 degree shaped pulse
;SHPW[6](F3) triple2		p29: f3 channel - shaped pulse for adiabatic ramping
;SHPW[7](F1) default		p8 : f1 channel - 90 degree shaped pulse (wet)
;SHPW[7](F1) lcnmr		p11: f1 channel - 90 degree shaped pulse (wet)
;SHPW[8](F1) triple_c		p12: f1 channel - 180 degree shaped pulse (selective)
;SHPW[8](F1) triple		p39: f1 channel - 120 degree shaped pulse (NH, best-)
;SHPW[8](F2) triple+triple2		p14: f2 channel - 180 degree shaped pulse (selective)
;SHPW[8](F3) triple_na		p32: f3 channel - 180 degree shaped pulse (_NA: N, adiabatic)

:SHPW[9](F1)	triple	p40: f1 channel - 180 degree shaped pulse (NH, best-, I)
:SHPW[9](F2)	triple+triple2	p23: f2 channel - 90 degree shaped pulse (higher selectivity)
:SHPW[10](F1)	triple	p41: f1 channel - 90 degree shaped pulse (NH, best-, I)
:SHPW[10](F2)	triple_c	p13: f2 channel - 90 degree shaped pulse (H, selective)
:SHPW[11](F1)	triple_c	p25: f1 channel - 180 degree shaped pulse (higher selectivity)
:SHPW[11](F2)	triple+triple2	p24: f2 channel - 180 degree shaped pulse (higher selectivity)
:SHPW[12](F1)	triple	p42: f1 channel - 180 degree shaped pulse (NH, best-, II)
:SHPW[12](F2)	default+triple+triple2	p63 : f2 channel - 180 degree shaped pulse (adiabatic decoupling)
:SHPW[13](F1)	triple	p43: f1 channel - 90 degree shaped pulse (NH, best-, II)
:SHPW[14](F2)	default	p39: f2 channel - 180 degree shaped pulse (short, broadband)
:SHPW[15](F1)	triple_c	p23: f1 channel - 180 degree shaped pulse (med. selectivity)
:SHPW[15](F1)	triple	p44: f1 channel - 180 degree shaped pulse (broadband, best-)
:SHPW[15](F2)	triple2	p33: f2 channel - 180 degree shaped pulse (med. selectivity)
:SHPW[16](F1)	triple2	p15: f1 channel - 90 degree shaped pulse (H2O on resonance)
:SHPW[16](F2)	triple2	p34: f2 channel - 180 degree shaped pulse (high selectivity)
:SHPW[17](F1)	triple2	p7 : f1 channel - 180 degree shaped pulse (H2O on resonance, cleanex))
:SHPW[17](F2)	triple2	p35: f2 channel - 90 degree shaped pulse (high selectivity)
:SHPW[18](F1)	default	p32: f1 channel - 180 degree shaped pulse for inversion (adiabatic: z-spoil)
:SHPW[19](F1)	triple_na	p12: f1 channel - 180 degree shaped pulse (_NA: H)
:SHPW[19](F1)	default	p31: f1 channel - 180 degree shaped pulse for inversion (adiabatic matched sweep)
:SHPW[19](F2)	triple	p15: f2 channel - 180 degree shaped pulse for inversion (adiabatic matched sweep)
:SHPW[20](F2)	triple	p31: f2 channel - 180 degree shaped pulse (sel. Ca or CO decoupling)
:	+triple_na	
:SHPW[21](F1)	triple_c	p30: f1 channel - 180 degree shaped pulse (sim. Ca + CO decoupling)
:SHPW[21](F1)	triple	p49: f1 channel - 180 degree shaped pulse (H2O)
:SHPW[21](F2)	triple2	p30: f2 channel - 180 degree shaped pulse (sim. Ca + CO decoupling)
:SHPW[22](F2)	triple2	p31: f2 channel - 180 degree shaped pulse (Cbeta decoupling)
:SHPW[22](F2)	triple*	pcpd8: f2 channel - 180 degree shaped pulse (Cbeta decoupling)
:SHPW[23](F2)	triple_na	p13: f2 channel - 90 degree shaped pulse (_NA: C)
:SHPW[25](F2)	triple_na	p14: f2 channel - 180 degree shaped pulse (_NA: C)
:SHPW[26](F2)	triple_na	p15: f2 channel - 90 degree shaped pulse (_NA: C, higher sel.)
:SHPW[28](F2)	triple_na	p10: f2 channel - 180 degree shaped pulse (_NA: C, higher sel.)
:SHPW[28](F2)	triple_na	p31: f2 channel - 180 degree shaped pulse (_NA: C, decoupling)
:SHPW[29](F2)	triple_na	p23: f2 channel - 180 degree shaped pulse (_NA: C, twofold mod)
:SHPW[31](F2)	triple	p47: f2 channel - shaped pulse 180 degree (inversion (sharp))
:SHPW[32](F2)	triple	p48: f2 channel - shaped pulse 180 degree (refocussing (sharp))
:		
:SQPL[0](F1)	all(!triple_c)	pl19: f1 channel - power level for CPD/BB decoupling
:SQPL[0](F2)	all	pl12: f2 channel - power level for CPD/BB decoupling
:SQPL[0](F2)	default+triple+triple2	pl30: f2 channel - power level for CPD/BB decoupling
:	+triple_na	
:SQPL[0](F3)	all	pl16: f3 channel - power level for CPD/BB decoupling
:SQPL[0](F4)	all	pl17: f4 channel - power level for CPD/BB decoupling
:SQPL[1](F1)	all	pl10: f1 channel - power level for TOCSY-spinlock
:SQPL[1](F1)	default	pl29: f1 channel - power level for TOCSY-spinlock (STD SL filter)
:SQPL[1](F2)	all	pl15: f2 channel - power level for TOCSY-spinlock
:SQPL[1](F3)	default+lcnmr+triple_c	pl23: f3 channel - power level for TOCSY-spinlock
:SQPL[1](F3)	triple	pl27: f3 channel - power level for TOCSY-spinlock
:SQPL[2](F1)	all	pl11: f1 channel - power level for ROESY-spinlock
:SQPL[3](F1)	all(!triple_c)	pl9 : f1 channel - power level for presaturation
:SQPL[3](F2)	default+lcnmr	pl21: f2 channel - power level for presaturation
:SQPL[3](F3)	lcnmr	pl22: f3 channel - power level for presaturation
:SQPL[4](F2)	default+lcnmr+triple_c	pl13: f2 channel - power level for second CPD/BB decoupling

```

;SQPL[5](F2) lcnmr+triple+triple2 pl14: f2 channel - power level for low power decoupling
;SQPL[5](F3) triple pl26: f3 channel - power level for low power decoupling
;SQPL[6](F2) default+triple+triple2 pl31: f2 channel - power level for bilev dec. (cw part)
;
+triple_na
;SQPL[7](F1) triple pl21: f1 channel - power level for TOCSY med. sel.
;SQPL[7](F2) triple_na pl20: f2 channel - power level for TOCSY med. sel.
;SQPL[7](F2) triple pl22: f2 channel - power level for TOCSY med. sel.
;SQPL[8](F1) triple_na pl25: f1 channel - power level for hetero TOCSY
;SQPL[8](F2) triple pl20: f2 channel - power level for hetero TOCSY
;SQPL[8](F2) triple_na pl26: f2 channel - power level for hetero TOCSY
;SQPL[8](F3) triple_na pl23: f3 channel - power level for hetero TOCSY
;SQPL[9](F2) triple_na pl27: f2 channel - power level for hetero TOCSY higher sel.
;SQPL[9](F3) triple_na pl22: f3 channel - power level for hetero TOCSY higher sel.
;SQPL[10](F1) triple2 pl27: f1 channel - power level for CLEANEX spinlock
;SQPL[11](F1) default pl27: f1 channel - power level for pulsed ROESY-spinlock
;SQPL[12](F1) all(!triple_c) pl32: f1 channel - power level for low power presaturation
;SQPL[12](F2) lcnmr pl26: f2 channel - power level for cw decoupling
;SQPL[13](F2) default pl14: f2 channel - power level for cw saturation
;SQPL[14](F2) all pl24: f2 channel - power level for hd/hc decoupling
;SQPL[15](F3) triple* pl25: f3 channel - power level for T1rho spinlock
;SQPL[16](F3) triple+triple2 pl23: f3 channel - power level for Rexchange/T2
;
;SQPW[0](F1) triple+triple2+triple_na p26: f1 channel - 90 degree pulse at pl19
;SQPW[0](F1) all pcpd1: f1 channel - 90 degree pulse for CPD decoupling
;SQPW[0](F2) all pcpd2: f2 channel - 90 degree pulse for CPD decoupling
;SQPW[0](F3) all pcpd3: f3 channel - 90 degree pulse for CPD decoupling
;SQPW[0](F4) triple* pcpd4: f4 channel - 90 degree pulse for CPD decoupling
;SQPW[1](F1) all p6 : f1 channel - 90 degree low power pulse
;SQPW[1]*0.66(F1) all(!triple_na) p5 : f1 channel - 60 degree low power pulse
;SQPW[1]*2(F1) all(!triple2+triple_na) p7 : f1 channel - 180 degree low power pulse
;SQPW[1](F2) all p9 : f2 channel - 90 degree low power pulse
;SQPW[1]*2(F2) all(!triple2+triple_na) p10: f2 channel - 180 degree low power pulse
;SQPW[1](F3) triple p25: f3 channel - 90 degree pulse at pl23
;SQPW[2](F1) default+lcnmr p15: f1 channel - pulse for ROESY spinlock
;SQPW[5](F2) lcnmr p31: f2 channel - 90 degree low power pulse (decoupling)
;SQPW[5](F2) triple+triple2 p61: f2 channel - 90 degree low power pulse (decoupling)
;SQPW[5](F3) triple+triple2 p62: f3 channel - 90 degree low power pulse (decoupling)
;SQPW[7](F2) triple_na p7 : f2 channel - 90 degree low power pulse (TOCSY med. sel.)
;SQPW[7](F2) triple p46: f2 channel - 90 degree low power pulse (TOCSY med. sel.)
;SQPW[8](F2) triple p45: f2 channel - 90 degree low power pulse ((hetero) TOCSY high sel.)
;SQPW[8](F3) triple_na p25: f3 channel - 90 degree low power pulse ((hetero) TOCSY high sel.)
;SQPW[9](F3) triple_na p24: f3 channel - 90 degree low power pulse ((hetero) TOCSY very high sel.)
;SQPW[10](F1) triple2 p10: f1 channel - 180 degree low power pulse (CLEANEX spinlock)
;SQPW[11](F1)*2 default+lcnmr p25: f1 channel - 90 degree pulse at pl27 (pulsed ROESY)
;SQPW[16]*2(F3) triple p30: f3 channel - 180 degree pulse at pl23 (T2)
;SQPW[16](F3)*2 triple2 p25: f3 channel - 180 degree low power pulse (Rexchange)
;
;P_grad1 all p16: homospoil/gradient pulse
;P_grad2 all p19: homospoil/gradient pulse 2
;P_hsqc all p28: f1 channel - trim pulse at pl1
;P_mlev all p17: f1 channel - trim pulse at pl10
;P_mlev all p20: f2 channel - trim pulse at pl15
;P_mlev triple_na p33: f3 channel - trim pulse at pl23

```

Complete list of pulse programs/files included into the \$HOME/exp/stan/nmr/lists/pp directory (Topspin v3.0)

Avance.incl	c_caco_ia	c_hncaco_s33d
Daz.incl	c_caco_s3	c_hnco_ia3d
De.incl	c_caco	c_hncoca2_ia3d
Delay.incl	c_can_iasq	c_hncoca_ia3d
Grad.incl	c_can_mq.2	calibgp
Param.info	c_can_mq	cbcaconhgp3d
Pulprog.info	c_canco_ia3d.2	cbcaconhgpwg3d.2
README	c_canco_ia3d	cbcaconhgpwg3d
Relations.info	c_cancoi_ia3d	cbcaconhgpwg4d
Sysconf.incl	c_cbcaco_ia3d	cbcanhgp3d
Update.info	c_cbcaco_s33d	cbcanhgpwg3d
adeq11etgpjcrdsp	c_cbcaco_ia3d	ccaconhgp2h3d
adeq11etgprdsp.2	c_cbcanco_ia3d	ccaconhgp3d.2
adeq11etgprdsp	c_ccco_ia3d	ccaconhgp3d
adeq11etgpsp	c_ccco_s33d	ccanhgp2h3d
adeq1netgp	c_cccon_ia3d	ccanhgp3d.2
adeqn1etgp	c_ccflopsy16_ctia	ccanhgp3d
adeqnnetgp	c_ccflopsy16_ct	ccconhgp2h3d
aptjc	c_ccflopsy16_ia	ccconhgp3d
apt	c_ccflopsy16	clmlevphpr
aring2	c_ccnoesy2	colocqf
aring	c_ccnoesy_ct	cosycwgppsqf
atocsygpph19	c_ccnoesy	cosycwphps
b_hncacbgp3d	c_coca_ia	cosydclrqf
b_hncacbigp3d	c_coca_mq.2	cosydcphwt
b_hncacogp3d	c_coca_mq	cosydcph
b_hncacogp4d	c_coca	cosydcqf
b_hncagp3d	c_con_iasq	cosydfesgpphpp
b_hncaigp3d	c_con_mqia	cosydfesgpph
b_hncocacbgp3d	c_con_mq	cosydfetgp.1
b_hncocacbgp4d	c_con_sq	cosydfetgp.2
b_hncocagp3d	c_cosy2_ct	cosydfetgppp.2
b_hncocagp4d	c_cosy_ct	cosydfgpph19
b_hncogp3d	c_cosy	cosydfphpp
b_hncoigp3d	c_hacaco_3d	cosydfphpr
b_hsqcetf3gpsi	c_hcaco_ia3d	cosydfph
b_trhncacbgp3d	c_hcaco_s33d	cosydfqf
b_trhncacbigp3d	c_hcan_ia3d	cosyetgp
b_trhncacogp3d	c_hcanco_ia3d	cosygpmfphpp
b_trhncagp3d	c_hcancoi_ia3d	cosygpmfph
b_trhncaignp3d	c_hcbca_ia3d	cosygpmfppqf
b_trhncocacbgp3d	c_hcbcaco_ia3d	cosygpmfqf
b_trhncocagp3d	c_hcbcaco_s33d	cosygpphzzfs
b_trhncogp3d	c_hcbcan_ia3d	cosygpppqf
b_trhncoigp3d	c_hccflopsy16_3d	cosygpprqf
b_trosyetf3gpsi	c_hnca_ia3d	cosygpqf
b_trosyf3gpph	c_hncaco_ia3d	cosyjdqf

cosylrqf	dipsi2gpphpr	hcaconhgp3d
cosyphpr	dipsi2gpphwg1	hcaconhgpwg3d
cosyphtf	dipsi2gpphwg2	hcangp2h
cosyph	dipsi2gpphwgxf	hcangp3d
cosyppqf	dipsi2gpphzs	hcannhgp3d
cosyprqf	dipsi2phpr	hcbcacong3d
cosyqf45	dipsi2ph	hccccosygp3d
cosyqf90	dipsihsqcf3gpsi3d	hcccgcbcaconhgp3d
cosyqfr2	dipsitretf3gp3d	hcccgcbcaconhgpwg3d
cosyqfr3	dqs1d	hccconhgp3d1
cosyqfrl	dqseagp135	hccconhgp3d2
cosyqftf	dqseagp90	hccconhgp3d3
cosyqf	dqsphpr	hccconhgpwg3d2
cpmg1d	dqsph	hccconhgpwg3d3
cpmgpr1d	dqsqf	hccconhgpwg4d
cpmg	dstebpgp3s1d	hcchatgp3d2
crineptgpph	dstebpgp3s	hcchatgp3d
croesyphpr	dstegp3s1d	hcchcogp3d2
croesyph	dstegp3s	hcchcogp3d
dazzg	ecos3cph	hcchcogp4d
dec180f3pr	ecos3nph	hcchcosygp3d
dec180pr	ecosygpgh	hcchdigp3d2
dec180sp	eretic30	hcchdigp3d
dec180	ereticgppn	hcchdigp4d
decp902hf4	ereticpr30	hcchecosgp3d
decp902h	ereticpr	hcchetgplr
decp90f3	eretic	hetsfhmqcf3gpph
decp90sp	fhsqccxf3gpph	hfcoqfn
decp90	fhsqcf3gpph	hhxcoqf.2
dept135	h2bcedetgpl3	hhxcoqf
dept45	h2bcetgpl3	hjrseqf
dept90	hacahbcosygp3d	hmbcacbigpl2ndqf
deptnd	haconhgpwg3d	hmbcacgplpndqf.2
deptppnd	hanhgpwg3d	hmbcacgplpndqf
deptqgppsp.2	hathmbcetgpl3	hmbcacgplpqf
deptqgppsp	hcbcbgcchargp	hmbcctetgpl2nd
deptqsp	hcbcbgcdcehegp	hmbcctetgpnd
deptsp135	hcbcbgcdhdgp	hmbcedetgpl3nd
deptsp45	hbhacbcconhgpwg5d	hmbcctf3gpl3nd
deptsp90	hbhacbcconhgpwg4d	hmbcctf3gpnd
deptsp	hbhaconhgp3d	hmbcctgpjcl2nd
dept	hbhaconhgpwg3d	hmbcctgpl2ndsm
dipsi2esfbgpph	hbhaconhgpwg4d	hmbcctgpl2nd
dipsi2esgpph	hbhanhgp3d	hmbcctgpl3nd
dipsi2etgpjcsix1	hbhanhgpwg3d	hmbcctgpnd
dipsi2etgpsi19	hcacogp3d	hmbcf3gplpndqf
dipsi2etgpsi	hcacogpjc3d	hmbcf3gpndqf
dipsi2etgp	hcacong2h	hmbcgppl2ndqf
dipsi2gpph19	hcacong3d	hmbcgpplpndprqf

hmbcgplpndqf	hmqcetgp	hncacosygp3d
hmbcgppndqf	hmqcf3gpph19.2	hncadhgp3d
hmbclpndqf	hmqcf3gpph19	hncadqzqgp3d
hmbcndprqf	hmqcf3gpphwg	hncacosgp3d2
hmbcndpsqf	hmqcf3phpr	hncacosgp3d
hmcmbcbacogpwg3d	hmqcf3ph	hncagp2h3d
hmcmbcbaconhgpwg3d	hmqcfbph	hncagp3d
hmcmbcbagpwg3d	hmqcgpmlqf	hncagprc3d1
hmcmbcbcanhgpwg3d2	hmqcgpnd1d	hncagprcwg3d1
hmcmbcbcanhgpwg3d	hmqcgpqfrl	hncagpwg2h3d
hmqcbidfphrl	hmqcgpqf	hncagpwg3d
hmqcbif3ph	hmqcjgp	hncah3d
hmqcbimlndph	hmqcmlevbi3d	hncahagp3d
hmqcbimlndqf	hmqcmIndphpr	hncaign2h3d
hmqcbimlph	hmqcmIndph	hncaign3d
hmqcbimlqf	hmqcmIphpr	hncaignpwg2h3d
hmqcbind1d	hmqcmIph	hncaignpwg3d
hmqcbindphrl	hmqcnd1d	hncajcgp3d
hmqcbindph	hmqcndph	hncannhgp2h3d.2
hmqcbindqfrl	hmqcndqf	hncannhgp2h3d
hmqcbindqf	hmqcndrd1d	hncannhgp3d.2
hmqcbindr1d	hmqcnoesybi3d	hncannhgp3d
hmqcbinoph	hmqcnoesyhmqcnngp4d	hncannhgpwg2h3d.2
hmqcbiphpr2	hmqcnoesyhsqcncgp4d	hncannhgpwg2h3d
hmqcbiphpr	hmqcnophpr	hncannhgpwg3d.2
hmqcbiphrl	hmqcph11	hncannhgpwg3d
hmqcbiph	hmqcphpr	hncggp3d.1
hmqcbiqf	hmqcph	hncggp3d.2
hmqcbird1d	hmqcqf	hncocacagp3d
hmqcctetgp.2	hmqcrd1d	hncocacbgp2h3d
hmqcctetgp	hmscetgpnd	hncocacbgp3d
hmqcdietf3gpsi.2	hncacbgp2h3d	hncocacbgpwg2h3d
hmqcdietf3gpsi	hncacbgp3d	hncocacbgpwg3d
hmqcdietf3gp	hncacbgpj3d	hncocacgp3d
hmqcdietgpsi.2	hncacbgpwg2h3d	hncocagp2h3d
hmqcdietgpsi	hncacbgpwg3d	hncocagp2h4d
hmqcdietgp	hncacbigp2h3d	hncocagp3d
hmqcetf3gpno	hncacbigp3d	hncocagp4d
hmqcetf3gpro.2	hncacbigpwg2h3d	hncocagpwg2h3d
hmqcetf3gpro	hncacbigpwg3d	hncocagpwg2h4d
hmqcetf3gpsi.2	hncacogp2h3d.2	hncocagpwg3d
hmqcetf3gpsi	hncacogp2h3d	hncocagpwg4d
hmqcetf3gp	hncacogp2h4d	hncocannhgp2h3d.2
hmqcetgp.2	hncacogp3d	hncocannhgp2h3d
hmqcetgpno	hncacogp4d	hncocannhgp3d.2
hmqcetgpro.2	hncacogpwg2h3d	hncocannhgp3d
hmqcetgpro	hncacogpwg2h4d	hncocannhgpwg2h3d.2
hmqcetgpsi.2	hncacogpwg3d	hncocannhgpwg2h3d
hmqcetgpsi	hncacogpwg4d	hncocannhgpwg3d.2

hncocannhgpwg3d	hsqcdietgpsisp.2	hsqcetgpsisp
hncocgp3d	hsqcdietgpsisp3d.2	hsqcetgpsiw
hncocogp3d	hsqcdietgpsisp	hsqcetgps
hncocogsp3d	hsqcdietgps	hsqcetgpsp.2
hncogp3d	hsqcedetf3gps	hsqcetgpsp.3
hncogphb3d.2	hsqcedetf3gpsi	hsqcetgpsp
hncogphb3d2	hsqcedetgpsisp.2	hsqcetgp
hncogphb3d	hsqcedetgpsisp2.2	hsqcexf3gpwgp
hncogpia	hsqcedetgpsisp2.3	hsqcf3gpaphsiwg
hncogprc3d1	hsqcedetgpsisp2.4	hsqcf3gpaphwg.2
hncogprc3d2	hsqcedetgpsisp2	hsqcf3gpaphwg
hncogprc3d3	hsqcedetgpsisp	hsqcf3gpdpwg
hncogprc3d4	hsqcedetgpsp.3	hsqcf3gpnowgx33d
hncogprc3d5	hsqcedetgpsp	hsqcf3gpnowgxy
hncogprc3d6	hsqcedetgp	hsqcf3gp19
hncogprc3d7	hsqcedgpph	hsqcf3phpr
hncogprcw3d1	hsqcetexf3gp	hsqcf3ph
hncogprcw3d2	hsqcetf3gpml	hsqcfpf3gp19
hncogprcw3d3	hsqcetf3gpno	hsqcgpmlph
hncogprcw3d4	hsqcetf3gpro.2	hsqcgpnd1d
hncogprcw3d5	hsqcetf3gpro	hsqcgpnowgx33d
hncogprcw3d6	hsqcetf3gps	hsqcgpph2
hncogprcw3d7	hsqcetf3gpsi	hsqcgpph
hncogpwg3d	hsqcetf3gpss	hsqcnoef3gps
hnconoef3gps	hsqcetf3gpxy.2	hsqcnoegpsi.2
hncorexf3gps	hsqcetf3gpxy	hsqcnoegpsi
hncot1f3gps	hsqcetf3gp	hsqcnoesyhmccngp4d
hncot2f3gps	hsqcetf3gpjcsi	hsqcnoesyhsqcccgp4d
hncotr3f3gps	hsqcetf3gps	hsqcnoesyhsqccngp4d
hnhagp3d	hsqcetf3gp	hsqcnoesyhsqcncgp4d
hnhbecogsp3d	hsqcetf3gp	hsqcnoesyhsqcngp4d
hnhbgp3d	hsqcetf3gppsi2	hsqcphpr
hnnhdigp3d	hsqcetgpasp	hsqcphps
hoesyfhqfqrnv	hsqcetgpjclrndxy	hsqcph
hoesyph	hsqcetgpjclrnd	hsqcrexetf3gps3d
hoesyqfrv	hsqcetgplrjcs	hsqct1etf3gps.2
hsqcccf3gp19	hsqcetgplrsp	hsqct1etf3gps3d.2
hsqccetgpjclr	hsqcetgpml3d.2	hsqct1etf3gps3d
hsqccetgpjc	hsqcetgpml	hsqct1etf3gps
hsqccetgpsisp	hsqcetgpno	hsqct1etgps3d.2
hsqccetgpsp.2	hsqcetgpprsisp2.2	hsqct1etgps3d
hsqccetgpsp	hsqcetgprosp.2	hsqct2etf3gps3d
hsqcdhetgpsp	hsqcetgprosp	hsqct2etf3gps
hsqcdiedetgpsisp.1	hsqcetgps	hsqct2etgps3d.2
hsqcdiedetgpsisp.2	hsqcetgpsismsp	hsqctretf3gps3d.2
hsqcdiedetgpsisp.3	hsqcetgpsisp.2	hsqctretf3gps3d
hsqcdietf3gps	hsqcetgpsisp2.2	hsqctretf3gps
hsqcdietgpiasp	hsqcetgpsisp2.3	hxcobicpqf
hsqcdietgpjcnisp	hsqcetgpsisp2	hxcobiqf

hxcocpqf	lc1pncwfd	mlevhmqcpr3d
hxcoqf2h	lc1pncwps	mlevhsqcetf3gp3d
hxcoqf	lc1pncw	mlevhsqcetgp3d
hxdeptbiph	lc1pnf2	mlevnoesy3d
hxdeptbiqf	lc1pnfr	mlevphpp
hxdeptmlph	lc1pnft	mlevphpr.2
hxdeptph	lc1pnps	mlevphpr
hxdeptqf	lc1prf2	mlevph
hxinepph	lc1prft	mlevroesy3d
hxinepqf	lc2grdonflow	mqsqgp1d2
hxxcoqf	lc2pg	mqsqgp1d
ideptnd	lc2pnf2	music_cm_3d_2
iineptnd	lc2pnps	music_cm_3d
iineptrd	lc2pn	music_de_3d_2
imgegp1d19f	lc2prf2	music_de_3d
imgegp1d2h	lc2pr	music_fhyw_3d_2
imgegp1d	lc2ps	music_fhyw_3d
imgegp2d	lc2wetdcus	music_gly_3d_2
imgegp3d	lc2wetdc	music_gly_3d
imgegpcv1d	lc2	music_ile_3d_2
imgegpcvsp1d2h	lcjrescwfdprqf	music_ile_3d
imgegpsp1d	lcjrescwprqf	music_kr_3d_2
imrcbgegp3d	lcjrescwpsqf	music_kr_3d
inad1d	lcjresf2prqf	music_lavia_3d_2
inadcp1d	lcjresprqf	music_lavia_3d
inadgpqfsp	lcjrespsqf	music_pro_1_3d.2
inadphsp	lcmlevcwfdpcph	music_pro_1_3d
inadph	lcmlevcwpcphps	music_pro_2_3d.2
inadqf.2	lcmlevf2pcph	music_pro_2_3d
inadqfsy	lcmlevf2phpr	music_qn_3d_2
inadqf	lcmlevpcphps	music_qn_3d
inadrd1d	lcmlevpcph	music_ser_3d_2
indecobimlph	ledbpgp2s1d	music_ser_3d
indecobiph	ledbpgp2s	music_tavi_3d_2
indecobird1d	ledbpgpco2s3d	music_tavi_3d
indecoph	ledbpgpml2s192d	music_trpe_2d
indecord1d	ledbpgpml2s2d	na_c6noesyhsqcgp3d
inepin	ledbpgpml2s3d	na_ghcchfwdigp3d
ineptnd	ledbpgpno3s3d	na_h56c56c4n3h_3d2
ineptpnd	ledbpgppr2s1d	na_h56c56c4n3h_3d
ineptrdsp	ledbpgppr2s	na_h5c5c4n3h_3d
ineptrd	ledgp2s1d	na_h68c68n19c42
ineptrl1	ledgp2s	na_hcchecgp3d
ineptrl2	mlevdcphwt	na_hcchfwdiecgp3d
jmod	mlevesgpph	na_hcchfwdiecgpjr3d
jresdcqf	mlevetgp	na_hcchfwdigp3d
jresgpprqf	mlevgpph19	na_hccnhdigpww3d
jresqf	mlevgpphprzf	na_hcncchcomq3d
lc1cpcwps	mlevgpphw5	na_hcncchdimq3d

selcogp	stddiffesgp.2	trhncacbgp2h3d
selcogp	stddiffesgp.3	trhncacbgp3d
selcorl	stddiffesgp	trhncacbietgp2h3d
selco	stddiffgp19.2	trhncacbietgp3d
selcssf.2	stddiffgp19.3	trhncacbigp2h3d
selcssfdizs.2	stddiffgp19	trhncacbigp3d
selcssfdizs	stddiff	trhncacoetgp2h3d
selcssfnozs.2	stdhsqcetgpsisp	trhncacoetgp2h4d
selcssfnozs	stdhsqcetgpsp	trhncacoetgp3d
selcssfzs.2	stdmlevesgpph	trhncacoetgp4d
selcssfzs	stdmlevgpph19	trhncacogp2h3d
selcssf	stdmlevph	trhncacogp2h4d
seldigpzs	stdnoesyegpph.2	trhncacogp3d
seldigp	stdnoesyegpph	trhncacogp4d
selgpse	stdnoesygpph.2	trhncacetgp2h3d
selhsqcgplrndsp	stdnoesygpph19.2	trhncacetgp3d
selhsqcgpndnosp	stdnoesygpph19	trhncagp2h3d2
selhsqcgpnosp	stdnoesygpph	trhncagp2h3d
selhsqcgpsisp	stebpgp1s191d	trhncagp3d.2
selhsqcgpsp	stebpgp1s19	trhncagp3d2.2
selina	stebpgp1s1d	trhncagp3d2
selineptlrrdsp	stebpgp1s	trhncagp3d
seljresqfsp	stebpgpin1s1d	trhncagprc3d1
selmlgp.2	stebpgpin1s	trhncaietgp2h3d
selmlgp	stegp1s1d	trhncaietgp3d
selmlzf	stegp1s	trhncaignp2h3d
selnogpzs	t1ir1d	trhncaignp3d
selnogp	t1irpg	trhncannhetgp2h3d.2
selnozf	t1ir	trhncannhetgp2h3d
selno	trcbcaconhetgp3d	trhncannhetgp3d.2
selrogp.2	trcbcaconhgp3d	trhncannhetgp3d
selrogp	trcbcanhetgp3d	trhncannhgp2h3d.2
selro	trcbcanhgp3d	trhncannhgp2h3d
selzggp	tretextf3gp	trhncannhgp3d.2
selzg	trhbhaconhetgp3d	trhncannhgp3d
seqtrcbcanhgp3d	trhbhaconhgp3d	trhncocacbetgp2h3d
seqtrhncacbetgp2h3d	trhbhanhetgp3d	trhncocacbetgp3d
seqtrhncacbetgp3d	trhbhanhgp3d	trhncocacbgp2h3d
seqtrhncacbgp2h3d	trhcccgcbaconhetgp3d	trhncocacbgp3d
seqtrhncacbgp3d	trhcccgcbaconhgp3d	trhncocaetgp2h3d
seqtrhnaetgp2h3d	trhccconhetgp3d2	trhncocaetgp2h4d
seqtrhnaetgp3d	trhccconhetgp3d3.2	trhncocaetgp3d
seqtrhncagp2h3d	trhccconhetgp3d3	trhncocaetgp4d
seqtrhncagp3d	trhccconhgp3d2	trhncocagp2h3d
sfhmqcf3gpphiasi	trhccconhgp3d3.2	trhncocagp2h4d
sfhmqcf3gpph	trhccconhgp3d3	trhncocagp3d
shmbcgpndqf	trhcchcogp3d	trhncocagp4d
stddiff.2	trhncacbetgp2h3d	trhncocannhetgp2h3d.2
stddiff.3	trhncacbetgp3d	trhncocannhetgp2h3d

trhncocannhetgp3d.2	troseyzqgpphwg	zgesfpgp
trhncocannhetgp3d	trrexetf3gpsi3d	zgesgp
trhncocannhgp2h3d.2	trt1etf3gpsi3d.2	zgf2hfpr
trhncocannhgp2h3d	trt1etf3gpsi3d	zgf2pr
trhncocannhgp3d.2	trt1etf3gpsi	zgf3ig
trhncocannhgp3d	trt2etf3gpsi3d	zgfbig
trhncogp3d	trt2etf3gpsi	zgfhighqn.2
trhncogp2h3d	trtretf3gpsi3d	zgfhighqn
trhncogp3d	undeft	zgf1hdqn
trhncogphb3d	wetdc	zgf1qn
trhncogprc3d1	wetdw	zggd30
trhncogprc3d2	wet	zggd
trhncogprc3d3	xhcoetgp	zggegp
trhncogprc3d4	zg0cw	zggppr
trhncogprc3d5	zg0dc	zggpse
trhncogprc3d6	zg0f2pr	zggpw5
trhncogprc3d7	zg0gd	zggpwg
trhncorexf3gp	zg0hd	zghc.2
trhncot1f3gp	zg0ig	zghc.3
trhncot2f3gp	zg0pg	zghc
trhncotr3gp	zg0pr	zghd.2
trnoef3gpsi	zg0	zghd
troesyphpr	zg19f	zghfigqn30
troesyph	zg2d	zghfigqn
troseyargpphwg	zg2h.2	zgif2hf4
troseyxf3gpphsi19	zg2h	zgif2h
troseyetf3gpiasi.2	zg30	zgif30
troseyetf3gpiasi	zgcd	zgif2igf3
troseyetf3gpsi.2	zgs	zgif
troseyetf3gpsi2	zgcpgppr	zgif30
troseyetf3gpsi	zgcppr	zgif
troseyetf3gpsi	zgcw30	zgif
troseyf3gpiaidphwg	zgcw	zgif
troseyf3gpiddphwg	zgcxesgp	zgif
troseyf3gpph19	zgcxgp19	zgif
troseyf3gpphsi19.2	zgcd30	zgif
troseyf3gpphsi19	zgcdf2igf3	zgif
troseygpphpr	zgcd	zgif
		zgs
		zgps
		zgpurge
		zgs
		zs_setup